



A VCL Model of a Cardiac Pacemaker

Jérôme Leemans and Nuno Amálio
Laboratory for Advanced Software Systems
University of Luxembourg
6, rue R. Coudenhove-Kalergi
L-1359 Luxembourg

TR-LASSY-12-04

Contents

Contents	2
1 Introduction	5
1.1 Heart and Pacemakers	5
1.2 The Visual Contract Language	6
1.3 Outline	6
2 The CommonPM package	8
3 The Battery package	10
3.1 State	10
3.2 Overall behaviour	10
4 The RateSmoothing package	12
4.1 State	12
4.2 Overall Behaviour	12
4.3 Global Operations	13
5 The AVDelay package	17
5.1 State	17
5.2 Overall Behaviour	18
5.3 Global Operations	19
6 The RateModulation Package	28
6.1 State	28
6.2 Overall Behaviour	30
6.3 Global Operations	31
7 The Pacing package	43
7.1 State	43
7.2 Overall Behaviour	44
7.3 Local operations of blob ChamberPacing	47
7.4 Global Operations	53
8 The Sensing Package	63
8.1 State	63
8.2 Overall Behaviour	64
8.3 Operations	65

9 The ResponseCommon Package	75
10 The AtrialResponse Package	76
10.1 State	76
10.2 Overall Behaviour	78
10.3 Global Operations	79
11 The VentricularResponse Package	94
11.1 State	94
11.2 Overall Behaviour	95
11.3 Global Operations	95
12 The Response package	101
12.1 State	101
12.2 Overall Behaviour	101
12.3 Global Operations	102
13 The HistoryCommon Package	106
14 The History Package	107
14.1 State	107
14.2 Overall Behaviour	108
14.3 Local operations of blob HistoryEvent	109
14.4 Local operations of blob AtrialEvent	109
14.5 Local operations of blob VentricularEvent	110
14.6 Global Operations	111
15 The PacingWithHistoryI package	112
15.1 State	112
15.2 Overall Behaviour	112
15.3 Global Operations	113
16 The ResponseWithHistoryI package	115
16.1 State	115
16.2 Overall Behaviour	115
16.3 Global Operations	117
17 The PacemakerWithoutHistory package	122
17.1 State	122
17.2 Overall Behaviour	122
18 The Pacemaker package	124
18.1 State	124
18.2 Overall Behaviour	124
18.3 Global Operations	125
19 Snapshot Analysis	127
References	131
A Requirements Issues	133

B Requirements Assumptions	135
C Generated Z Specification	137
C.1 Preamble	137
C.2 Package CommonPM	137
C.3 Package Battery	138
C.4 Package RateSmoothing	139
C.5 Package AVDelay	143
C.6 Package RateModulation	151
C.7 Package Pacing	164
C.8 Package Sensing	190
C.9 Package ResponseCommon	202
C.10 Package AtrialResponse	202
C.11 Package VentricularResponse	220
C.12 Package Response	227
C.13 Package HistoryCommon	236
C.14 Package History	236
C.15 Package PacingWithHistoryI	241
C.16 Package ResponseWithHistoryI	247
C.17 Package PacemakerWithoutHistory	256
C.18 Package Pacemaker	266
D ZOO Toolkit	277

Chapter 1

Introduction

Embedded software in medical devices is becoming ubiquitous and increasing in content and complexity. Because of this complexity and the fact that many failures in existing devices have been attributed to coding errors, regulators demand rigorous assurance that software running on medical devices is safe, reliable and functions correctly [JIJ06]. It is therefore not a surprise that there has been a recent interest in applying formal methods to this domain [JIJ06, JJA10].

This document presents a visual model of a cardiac pacemaker device [Bos07b] expressed in the Visual Contract Language (VCL)[AK10, AKMG10, AGK11]. The model is built according to an informal specification [Bos07b], released into the public domain by Boston Scientific, and which constitutes a pilot case study of the software verification challenge [WLBF09].

The model presented here has been developed using VCL's tool, the Visual Contract Builder (VCB) [AGK11]¹. It is a revised version of the model presented in [Lee11]; at the time [Lee11] was developed, VCB did not support all VCL diagram types (in particular, contract diagrams).

During the development of the VCL model, we encountered many issues with the requirements of the pacemaker device [Bos07b]. These issues are documented in appendix A. To proceed with the modelling effort, assumptions regarding these issues were taken; these are documented in appendix B.

1.1 Heart and Pacemakers

The heart is a beating muscle that continuously pumps blood to the rest of the body. It comprises four chambers (Fig. 1.1): two *atria* (top) and two *ventricles* (bottom). The heartbeat is the rhythmic contraction of the heart's four chambers, which is stimulated by electrical signals that travel through a specific nerve pathway in the heart. These electrical signals begin at the *sinoatrial node* (SA Node), which is located in the right atrium (see Fig. 1.1). The signal then travels through the atria, causing their contraction, to reach the *atrioventricular node* (AV node), which amplifies the signal sending it through the ventricles, causing their contraction and the consequent pumping of blood into the rest of the body. The heart's conduction system functions as the body's




Figure 1.1: Human heart and its conduction system

¹<http://vcl.gforge.uni.lu>

own pacemaker. If it does not work properly, due to heart damage or other medical conditions, it can cause an irregular heartbeat, compromising the flow of blood to the brain and other parts of the body.

A cardiac pacemaker (PM) is an electronic device implanted in the body to regulate the heartbeat [BSS10]. It should be used when the heart's natural conduction system is damaged: either due to a sick SA node (bradycardia), or a block that prevents the signal from reaching the AV node (heart block), or other heart malfunctions [MM07]. This document presents a VCL model of a cardiac pacemaker [Bos07b]. This is a pilot case study of the software verification challenge.

1.2 The Visual Contract Language

The Visual Contract Language (VCL) [AK10, AKMG10, AGK11] is a formal language designed for the abstract description of software systems. Its modelling paradigms are set theory, object-orientation and design-by-contract (pre- and post-conditions). VCL distinguishes itself from other visual modelling languages, such as UML, by its capacity to describe predicates visually.

VCL takes a generative (or translational) approach to semantics. Currently, VCL has a Z semantics; the Z specification generated from the diagrams constituting the formal semantics of a VCL model. The Z enables formal model analysis through theorem proving. The Z specifications generated from VCL follow ZOO [APS05, Amá07], an approach to build object-oriented models in Z (a Z style).

VCL's tool, the Visual Contract Builder (VCB) [AGK11]², an Eclipse plug-in, enables the construction of VCL diagrams, checks their well-formedness through type-checking, and generates Z specifications from them.

VCL's diagram suite comprises: *package*, *structural*, *behaviour*, *assertion* and *contract* diagrams. Package diagrams (PDs) define VCL packages, coarse-grained modules, and their dependencies with other packages (e.g. Figs. 2.1 and 4.1(a)). Structural diagrams (SDs) define structures and their relations that together make the state space of a package (e.g. Figs. 2.2 and 4.1(b)). Behaviour diagrams provide a map over the behaviour units of a package (e.g. Fig. 4.2). Assertion (or constraint) diagrams (ADs) define predicates over a single state, which can be used to define invariants and observe operations (or queries). Finally, contract diagrams (CDs) describe operations that change state through a contract (a pre- and a post-condition).

1.3 Outline

This document develops a VCL model of a cardiac pacemaker device [Bos07b], which is formally validated using snapshot analysis. This chapter introduced the case study and gave some background on VCL. The subsequent chapters are as follows:

- Chapters 2 to 18 present the packages that make the VCL model of the cardiac pacemaker.
- Chapter 19 presents the snapshot analysis performed on the Pacemaker VCL model.
- Appendix A documents the issues that we found with the requirements of the cardiac Pacemaker device documented in [Bos07b].
- Appendix B documents the assumptions that have been taken concerning the requirements issues that we found. The VCL model presented here is based on these assumptions.

²<http://vcl.gforge.uni.lu>

- Appendix C presents the Z generated from the VCL model developed in this document.
- Appendix D presents the ZOO toolkit (VCL's semantic domain), on which the Z specification rests.

Chapter 2

The CommonPM package

This package contains sets that are common across the different packages of the Pacemaker model. The package's PD (Fig. 2.1) defines `CommonPM` as a container package.



Figure 2.1: VCL Package diagram of package `CommonPM`

In the SD (Fig. 2.2), `BradycardiaMode` is a primary blob defining an enumerated set containing all possible operating modes of a DDD pacemaker [Bos07b, BSS10]; they are as follows: *O* (none), *A* (atrial), *V* (ventriculum), *D* (dual), *T* (triggered), *I* (inhibited), *D* (tracked), *R* (rate modulation). Other primary blobs defined this way are `Bool`, `Switch`, `BatteryStatus` and `PMEVENT`. Blob `BatteryStatus` holds all possible status of the device's battery as defined in [Bos07b, BSS10]; namely: `BOL` (beginning of life), `ERN` (elective replacement near), `ERT` (elective replacement time) and `ERP` (elective replacement past). Blob `PMEVENT` represents internal events that need to be propagated and result in some event action to be performed by some device component.

The derived blobs `ChambersPaced`, `ChambersSensed` and `ResponseToSensing` are defined extensionally from the `BradycardiaMode` primary blob. The derived blobs `MicroSec` (microseconds, μs), `MilliSec` (milliseconds, ms), `MicroVolt` (microvolts, μv), `MilliVolt` (millivolts, mv) and `PPM` (pulses per minute) are defined from the primitive blob `Nat` (represents set of natural numbers).




Figure 2.2: VCL Structural diagram of package `CommonPM`

Chapter 3

The Battery package

3.1 State

Package **Battery** is responsible for managing the status of the device's battery. The package's PD (Fig. 3.1(a)) imports package **CommonPM** in order to access the set **BatteryStatus**.

Package **Battery**'s SD (Fig. 3.1(b)) adds the property edge **bs** of imported set **BatteryStatus** to the package blob in order to keep the current status of the battery. At any one time, the battery status may have the values: **BOL** (beginning of life), **ERN** (elective replacement near), **ERT** (elective replacement time) and **ERP** (elective replacement past).




Figure 3.1: VCL Package, structural and behaviour diagram of package **Battery**

3.2 Overall behaviour

BD of package **Battery** (Fig. 3.1(c)) defines the following behavioural units:

- **Init** define the package's initialisation. Initialisation's defining AD (Fig. 3.2(a)) says that initially the battery is set to **BOL** (according to requirements assumption A7, appendix B).
- **UpdateBatteryStatus** is the operation that is called when the battery status level changes. Its definition (CD of Fig. 3.2(b)) receives the current battery status level as an input and changes the current status of the battery to the one received as input (according to requirements assumption A9, appendix B).




Figure 3.2: Assertion diagrams of Initialisation and contract diagram of `UpdateBatteryStatus` in package `Battery`




Figure 3.3: Assertion diagram of operation `SetBradycardiaMode` in package `Battery`

- `SetBradycardiaMode` describes the effect of this operation on this package. Its definition (CD of Fig. 3.3) describes a precondition that will affect the overall operation when composed; the precondition is described in AD `IsBatteryLevelHigh`, which says that `SetBradycardiaMode` may operate (changing the current bradycardia mode of the device) provided the battery status is `BOL` or `ERN`.

Chapter 4

The RateSmoothing package

Rate Smoothing is a programmable feature of Pacemakers designed to prevent sudden changes in the rhythm of pacing [BSS10]. A rate smoothing algorithm ensures gradual slowing or speeding of the pacemaker rate based on a percentage of a preceding cardiac interval [MM07, TLMB08]. Package **RateSmoothing** localises this functionality of the pacemaker device.

4.1 State




Figure 4.1: VCL Package and structural diagram of package **RateSmoothing**

PD of package **RateSmoothing** is given Fig. 4.1(a). It defines ensemble package **RateSmoothing**, saying that it imports package **CommonPM**. This import enables use of common sets defined in **CommonPM**.

Package **RateSmoothing**'s SD is given in Fig. 4.1(b). It introduces the derived blob **RateSmoothingValue**, which defines a percentage between 0 and 25 according to the “programmable parameters table” of [Bos07b, p. 34].

In Fig. 4.1(b), the property edges of package blob **RateSmoothing** define the package's state:

- **rsUp** and **rsDown** define the programmable parameters, representing, respectively, the limit for the increase and decrease in the current cardiac cycle interval.

4.2 Overall Behaviour

BD of package **RateSmoothing** (Fig. 4.2) introduces the following behavioural units:




Figure 4.2: Behavioural Diagram of package **RateSmoothing**

- **Init** (Fig. 4.3) define the package's initialisation.
- **IncRSDown** (Fig. 4.4) and **DecRSDown** (Fig. 4.5) increment and decrement the programmed value of the rate smoothing down percentage.
- **IncRSUp** (Fig. 4.6) and **DecRSUp** (Fig. 4.7) increment and decrement the programmed value of the rate smoothing up percentage.
- **CalcSmoothedCCI** (Fig. 4.8) calculates a new cardiac cycle interval based on the rate smoothing settings.

4.3 Global Operations




Figure 4.3: Initialisation of package **RateSmoothing**




Figure 4.4: Operation `IncRSDown` of package `RateSmoothing`




Figure 4.5: Operation `DecRSDown` of package `RateSmoothing`




Figure 4.6: Operation `IncRSUp` of package `RateSmoothing`




Figure 4.7: Operation `DecRSUp` of package `RateSmoothing`




Figure 4.8: Assertion Diagram of observe operation **CalcSmoothedCCI**




Figure 4.9: Assertion Diagrams of operations **CalcLBSmoothedCCI** and **CalcUBSmoothedCCI** of package **RateSmoothing**

Chapter 5

The AVDelay package

The atrioventricular (AV) Delay or atrioventricular interval (AVI) is the interval between an atrial event (either sensed or paced) and the scheduled delivery of a ventricular stimulus [MM07]. The package `AVDelay` manages this functionality of the pacemaker device.

5.1 State




Figure 5.1: VCL Package diagram of package `AVDelay`

PD of package `AVDelay` is given Fig. 5.1. It defines ensemble package `AVDelay`, saying that it imports package `CommonPM`. This import enables use of common sets defined in `CommonPM`.

Package `AVDelay`'s SD (Fig. 5.2) introduces the following derived blobs:

- `MinDynAVDelay` represents the allowed values for the minimum AV delay as defined in the “programmable parameters table” of [Bos07b, p. 34].
- `MaxDynAVDelay` represents the allowed values of both fixed AV delay and maximum dynamic AV delay. This is defined according to “Fixed AV Delay” in the “programmable parameters table” of [Bos07b, p. 34] and results from requirements assumption A1 (appendix B), which resolves issue RI1 in the requirements (appendix A).
- `SensedAVDelayOffset` represents the allowed values of the parameter sensed AV delay offset, expressed in milliseconds and which is used to shorten the AV Delay following a tracked atrial sense. This blob is defined according to “Sensed AV Delay Offset” interval defined in the “programmable parameters table” of [Bos07b, p. 34]

The SD's package blob (Fig. 5.2) has property edges `fixedAVI` and `maxDynDelay`, `minDynDelay`, `minDynDelay`, `dynMode` and `savDelayOffset`; all programmable features of the device to be set




Figure 5.2: VCL structural diagram of package AVDelay




Figure 5.3: Behavioural Diagram of package AVDelay

by physicians. Property **dynMode** (of imported set **Switch**, which has values **On** and **Off**) indicates the dynamic mode is on or off; when off the fixed AVI is used. The constant **dynDelayFactor**, a natural number whose value is undefined, represents the “factor” (mentioned in [Bos07b, p. 29]) that is used to calculate the dynamic AV delay; this constant emerges from requirements assumption A3 (appendix B).

5.2 Overall Behaviour

BD of package AVDelay (Fig. 5.3) introduces the following behavioural units:

- **Init** (Fig. 5.4(a)) define the package’s initialisation.
- **GetAVDelay** (Fig. 5.4(b)) defines an observe operation to yield the current AV delay. We consider that there is just one AV Delay, and not a sensed and paced AV Delay. This follows assumption A3 (appendix B).
- **DynamicDelayOn** and **DynamicDelayOff** turn the dynamic delay on and off.
- **IncFixedAVI** and **DecFixedAVI** increment and decrement the fixed AVI parameter (property edge **fixedAVI** above).

- `IncMinDynAvDely` and `DecMinDynAvDely` increment and decrement the parameter minimum AV Delay (property edge `minDynDelay` above).
- `IncMaxDynAVDelay` and `DecMaxDynDelay` increment and decrement the parameter maximum AV Delay (property edge `maxDynDelay` above).
- `IncSAVDelayOffset` and `DecSAVDelayOffset` increment and decrement the parameter sensed AV Delay offset (property edge `savDelayOffset` above).

5.3 Global Operations




Figure 5.4: Assertion Diagrams of Initialisation and observe operation `GetAVDelay` in package `AVDelay`




Figure 5.5: Contract Diagrams of operations `DynamicDelayOn` and `DynamicDelayOff`




Figure 5.6: Contract Diagrams of operation `IncFixedAVI`




Figure 5.7: Contract Diagrams of operation `DecFixedAVI`




Figure 5.8: Contract Diagrams of operation `IncMinDynAVDelay`




Figure 5.9: Contract Diagrams of operation `DecMinDynAVDelay`




Figure 5.10: Contract Diagrams of operation `IncMaxDynAVDelay`




Figure 5.11: Contract Diagrams of operation `DecMaxDynAVDelay`




Figure 5.12: Contract Diagrams of operation `IncSAVDelayOffset`




Figure 5.13: Contract Diagrams of operation `DecSAVDelayOffset`

Chapter 6

The RateModulation Package

Rate modulation or *rate adaptive* pacing is the ability of pacemakers to increase the pacing rate in response to physical activity or metabolic demand [MM07]. Package `RateModulation` deals with this functionality in the pacemaker device. In [Bos07b, sec. 5.7, p.32], this functionality is identified as rate adaptive pacing, which concerns the “ability to adjust the cardiac cycle in response to metabolic need as measured from body motion using an accelerometer.

6.1 State




Figure 6.1: VCL Package diagram of package `RateModulation`

`RateModulation`'s PD (Fig. 6.1) defines `RateModulation` as an ensemble package, which imports container package `CommonPM`.

SD of `RateModulation` is given in Fig. 6.2. It introduces several derived blobs representing the allowed values of programmable parameters; they are as follows:

- `MaximumSensorRate` (MSR) represents the possible maximum pacing rates allowed under sensor control [Bos07b, sec. 5.7.1, p. 32]. Blob `MaximumSensorRate` defines values that represent units as paces per minute (ppm); it is defined according to the “programmable parameters table” of [Bos07b, p. 34]. This is used when the patient increases activity and the sensor responds, the heart rate then rises with a programmed response up to a maximal rate defined by the MSR [MM07].
- `ResponseFactor` is used to calculate the pacing rate at various levels of steady state patient activity; depending on the programmed value of response factor, the incremental change of rate may be higher or lower [Bos07b, sec. 5.7.3, p. 32]. This is defined according to the “programmable parameters table” of [Bos07b, p. 34].




Figure 6.2: VCL structural diagram of package **RateModulation**

- **ReactionTime** is the time required to drive the rate from the lower rate limit (LRL) to MSR [Bos07b, sec. 5.7.4, p. 33]. Blob **ReactionTime** defines defines values that represent units as miliseconds (ms) and is defined according to the “programmable parameters table” of [Bos07b, p. 34].
- **RecoveryTime** is the time required for the rate to fall from MSR to LRL [Bos07b, sec. 5.7.5, p. 33]. Blob **RecoveryTime** defines defines values that represent units as miliseconds (ms) and is defined according to the “programmable parameters table” of [Bos07b, p. 34].
- **Activity** represents the various levels of patient activity as output by the accelerometer [Bos07b, sec. 5.7.2, p. 32]. This is defined according to the parameter *activity threshold* of the “programmable parameters table” of [Bos07b, p. 34] and following requirements assumption A5 (appendix B). The different values specified in [Bos07b, p. 34] are here in the model interpreted as natural numbers (1: V-Low; 2: Low; 3: Med-Low; 4: Med; 5:Med-High; 6: High; 7:V-High).
- **LowerRateLimit** represents the allowed values of *lower rate limit* (LRL); this is expressed in ppm and defined according to the “programmable parameters table” of [Bos07b, p. 34]. LRL represents the number of generator pace pulse delivered per minute (atrium or ventricle) in the absence of sensed intrinsic activity or sensor-controlled pacing at a higher rate [Bos07b, p. 28]. **LowerRateLimit**.

In Fig. 6.2, the property edges of the package blob **RateModulation** represent the state of the package. They are as follows:

- **msr**, **responseF**, **reactionTm** and **recoveryTm** represent the programmed values for parameters MSR, response factor, and reaction and recovery times.

- rm and prm represent the current and previous bradycardia mode for rate modulation.
- **threshold** and **lastActivity** are related to the accelerometer output. **threshold** is the value that enables adaptive pacing; if the accelerometer output exceeds this value then the pacemaker's rate is affected by activity data [Bos07b, sec. 5.7.2, p. 32]. **lastActivity** captures the last accelerometer output, following requirements assumption A5 (appendix B).
- **lastSwitchTm** is used by the rate-adaptation algorithm. It marks the time where there is a switch of the target cardiac cycle interval from the MSR to LRL or the other way round.

6.2 Overall Behaviour




Figure 6.3: VCL behavioural diagram of package RateModulation

BD of package RateModulation is given in Fig. 6.3. It introduces the following behavioural units:

- **Init** (Fig. 6.4(a)) is the package's initialisation.
- **SetBradycardiaMode** (Fig. 6.4(b)) is the operation that sets the bradycardia mode for rate modulation.
- **SetTargetCCI** (Fig. 6.5) is the operation that sets the target cardiac cycle interval taking into account rate modulation. This involves a complex processing according to what is described in [Bos07b, sec. 5.7, p. 32]. The diagrams from figures 6.5 to 6.9 are involved in this.
- **SetTargetCCIToLRL** (Fig. 6.10) is the operation that sets the target cardiac cycle interval to the LRL.
- **IncLRL** (Fig. 6.11) and **DecLRL** (Fig. 6.12) increment and decrement the programmed value of the lower rate limit.
- **IncMSR** (Fig. 6.13) and **DecMSR** (Fig. 6.14) increment and decrement the programmed value of MSR.
- **IncRF** (Fig. 6.15) and **DecRF** (Fig. 6.16) increment and decrement the programmed value of response factor.
- **IncReactionTm** (Fig. 6.17) and **DecReactionTm** (Fig. 6.18) increment and decrement the programmed value of reaction time.

- `IncRecoveryTm` (Fig. 6.19) and `DecRecoveryTm` (Fig. 6.20) increment and decrement the programmed value of recovery time.
- `IncActThreshold` (Fig. 6.21) and `DecActThreshold` (Fig. 6.22) increment and decrement the activity threshold.

6.3 Global Operations




Figure 6.4: Package `RateModulation`: Assertion diagram of Initialisation and contract diagram of operation `SetBradycardiaMode`




Figure 6.5: VCL contract diagram of operation `SetTargetCCI` in package `RateModulation`




Figure 6.6: VCL contract and assertion diagrams of operation `SetTargetCCISwitching` in package `RateModulation`




Figure 6.7: VCL assertion diagrams of predicate `IsAboveThreshold` in package `RateModulation`




Figure 6.8: VCL contract and assertion diagrams of operation `SetTargetCCIIIncreasing` in package `RateModulation`




Figure 6.9: VCL contract and assertion diagrams of operation `SetTargetCCIDecreasing` in package `RateModulation`




Figure 6.10: VCL contract diagram of operation `SetTargetCCIToLRL` in package `RateModulation`




Figure 6.11: Contract and assertion diagrams of global operation `IncLRL`




Figure 6.12: Contract and assertion diagrams of global operation `DecRL1`




Figure 6.13: Contract diagrams of operation `IncMSR` of package `RateModulation`




Figure 6.14: Contract diagrams of operation `DecMSR` of package `RateModulation`




Figure 6.15: Contract diagrams of operation `IncRF` of package `RateModulation`




Figure 6.16: Contract diagrams of operation `DecRF` of package `RateModulation`




Figure 6.17: Contract diagrams of operation `IncReactionTm` of package `RateModulation`




Figure 6.18: Contract diagrams of operation `DecReactionTm` of package `RateModulation`




Figure 6.19: Contract diagrams of operation `IncRecoveryTm` of package `RateModulation`




Figure 6.20: Contract diagrams of operation `DecRecoveryTm` of package `RateModulation`




Figure 6.21: Contract diagrams of operation `IncActThreshold` of package `RateModulation`




Figure 6.22: Contract diagrams of operation `DecActThreshold` of package `RateModulation`

Chapter 7

The Pacing package

7.1 State

Pacing is the functionality of the Pacemaker device that stimulates the heart through electric pulses. A pulse has an amplitude or peak value (a voltage) and a width or duration [BSS10]. The **Pacing** package presented here manages the pacing feature of the Pacemaker device; it is responsible for emitting electrical pulses to both the ventricular and atrial chambers of the heart.




Figure 7.1: VCL Package diagram of package **Pacing**

PD of **Pacing** (Fig. 7.1) package says that it imports container package **CommonPM** and ensemble packages **RateSmoothing** and **AVDelay**, which are sub-features of pacing.

Figure 7.2 presents SD of **Pacing**, which introduces the following derived blobs:

- **PulseWidth** represents the allowed values of pulse width, which is defined from set **MicroSec** defined in **CommonPM** as described in the “programmable parameters table” of [Bos07b, p. 34]. **PulseWidth** is expressed in microseconds (μV).
- **PulseAmplitudeRegulated** represents the allowed values of the pulse amplitude during the regulated mode. It is defined according to the “programmable parameters table” of [Bos07b, p. 34]. Pulse amplitudes are expressed in millivolts (mv).
- **PulseAmplitudeUnregulated** represents the allowed values of the pulse amplitude during the unregulated mode. It is defined according to the “programmable parameters table” of [Bos07b, p. 34].

The remaining elements of the SD of Fig. 7.2 are as follows:




Figure 7.2: VCL structural diagram of package Pacing

- Domain blob **ChamberPacing** factors the pacing properties of a chamber (atrial or ventricular); this includes the programming parameters pulse width and pulse amplitude (regulated and unregulated), the last time a pacing event occurred (**lastEvTm**) and the current amplitude mode (**amplMode**, which can either be regulated or unregulated).
- Package blob **Pacing** has the following property edges: **aP** and **vP** of blob **ChamberPacing** (atrial and ventricular pacing), **lrl** (lower rate limit), **hrl** (hysteresis rate limit), **rm** and **prm** of blob **RateModulation** (current and previous rate modulation mode), **hystMode** of blob **Switch** (hysteresis mode, which can be on or off), **cp** and **pcp** of blob **ChambersPaced** (current and previous pacing mode, which can either be atrial, ventricular or dual), and **currentCCI** and **targetCCI** (current and target cardiac cycle intervals).
- Assertion **TargetCCInv** represents an invariant describing constraint of property edge **targetCCI**, which is described in AD of Fig. 7.3. The AD describes three conditions: when rate modulation is on, the target CCI must be less than the interval calculated from the lower rate limit; when the rate modulation is off, then the target CCI is either calculated from the lower rate limit (if the hysteresis mode is off) or the hysteresis rate limit (if the hysteresis mode is on).

7.2 Overall Behaviour

BD of **Pacing** is given in Fig. 7.4. It includes two integral extensions to bring in operations from imported packages: one brings into the new context all operations from package **RateSmoothing**, except **Init** and **CalcSmoothedCCI**, and the other all operations from package **AVDelay**, except **GetAVDelay** and **Init**.

The BD (Fig. 7.4) defines several local operations for the blob **ChamberPacing**:

- **New** (Fig. 7.5) to create a new **ChamberPacing** object.




Figure 7.3: VCL assertion diagram of package TargetCCIIInv

- `SetLastEventTm` (Fig. 7.6(a)) sets the time of the last heart event (either sensed or paced).
- `IncPulseAmplReg` (Fig. 7.7) and `DecPulseAmplReg` (Fig. 7.8) increment and decrement the regulated pulse amplitude.
- `IncPulseAmplUReg` (Fig. 7.9) and `DecPulseAmplUReg` (Fig. 7.10) increment and decrement the unregulated pulse amplitude.
- `IncPulseWidth` (Fig. 7.11) and `DecPulseWidth` (Fig. 7.12) increment and decrement the pulse width.
- Observe operations `GetLastEventTm` (Fig. 7.6(b)) gets time of the heart pacing event, `GetPulseWidth` (Fig. 7.13(a)) gets current pulse's width, and `GetPulseAmpl` (Fig. 7.13(b)) gets the pulse amplitude.

The following global operations are defined in the BD of Fig. 7.4:

- `Init` (Fig. 7.14) is the package's initialisation.
- `SetBradycardiaMode` (Fig. 7.15) sets the bradycardia mode for pacing.
- `UpdateBatteryStatus` (Fig. 7.16) does the required action for pacing when the battery status is updated.
- `IncAPulseAmplReg` (Fig. 7.17(a)) and `DecAPulseAmplReg` (Fig. 7.17(b)) to increment and decrement the regulated atrial pulse amplitude.
- `IncVPulseAmplReg` (Fig. 7.18(a)) and `DecVPulseAmplReg` (Fig. 7.18(b)) to increment and decrement the regulated ventricular pulse amplitude.




Figure 7.4: Behavioural Diagram of package Pacing

- **IncAPulseAmplUReg** (Fig. 7.19(a)) and **DecAPulseAmplUReg** (Fig. 7.19(b)) to increment and decrement the unregulated atrial pulse amplitude.
- **IncVPulseAmplUReg** (Fig. 7.20(a)) and **DecVPulseAmplUReg** (Fig. 7.20(a)) to increment and decrement the unregulated ventricular pulse amplitude.
- **IncAPulseWidth** (Fig. 7.21(a)) and **DecAPulseWidth** (Fig. 7.21(b)) to increment and decrement the atrial pulse width.
- **IncVPulseWidth** (Fig. 7.22(a)) and **DecVPulseWidth** (Fig. 7.22(b)) to increment and decrement the ventricular pulse width.
- **IncHRL** (Fig. 7.23) and **DecHRL** (Fig. 7.24) increment and decrement the hysteresis rate limit.
- **SetTargetCCI** (Fig. 7.28) sets the value of the state property `targetCCI`.
- **ExecuteEvent** (Fig. 7.29) executes some internal pacemaker event in the context of package **Pacing**.
- **TriggerAPace** (Fig. 7.30(a)) and **TriggerVPace** (Fig. 7.30(b)) to trigger a pace on atrial (A) and ventricular (V) chambers.
- **TimedAPace** (Fig. 7.32(a)) and **TimedVPace** (Fig. 7.32(b)) to issue a timed pace (or stimulus) in the atrial (A) and ventricular (V) chambers.

- **InhibitAPace** (Fig. 7.31(a)) and **InhibitVPace** (Fig. 7.31(b)) to inhibit or not issue a pace in the atrial (A) and ventricular (V) chambers; the event is recorded as having happened, but no stimulus is issued.
- **TimedAPace** (Fig. 7.32(a)) and **TimedVPace** (Fig. 7.32(b)) delivers a scheduled stimulus in the atrial (A) or ventricular (V) chambers.
- **TrackASense** (Fig. 7.34) is the pacing action of a tracked sensing in the atrial chamber.

7.3 Local operations of blob ChamberPacing




Figure 7.5: Contract diagrams of operation `New` of blob `ChamberPacing`




Figure 7.6: Contract and assertion diagrams of operations `SetLastEventTm` and `GetLastEventTm` of blob `ChamberPacing`




Figure 7.7: Contract diagrams of operation `IncPulseAmplReg` of blob `ChamberPacing`




Figure 7.8: Contract diagrams of operation `DecPulseAmpReg` of blob `ChamberPacing`




Figure 7.9: Contract diagrams of operation `IncPulseAmplUReg` of blob `ChamberPacing`




Figure 7.10: Contract diagrams of operation `DecPulseAmplUReg` of blob `ChamberPacing`




Figure 7.11: Contract diagrams of operation `IncPulseWidth` of blob `ChamberPacing`




Figure 7.12: Contract diagrams of operation `DecPulseWidth` of blob `ChamberPacing`




Figure 7.13: Assertion diagrams of observe operations `GetPulseWidth` and `GetPulseAmp` of blob `ChamberPacing`

7.4 Global Operations




Figure 7.14: Initialisation of package `Pacing`




Figure 7.15: Contract Diagrams of global operation `SetBradycardiaMode`




Figure 7.16: Contract Diagrams of global operation `UpdateBatteryStatus`




Figure 7.17: Contract Diagrams of global operations `IncAPulseAmplReg` and `DecAPulseAmplReg`

The global operation `UpdateBatteryStatus` (Fig. 7.16) considers two cases: (a) when the bradycardia mode needs to be changed because the battery becomes close to depletion (contract `UpdateBatteryChangeMode`), (b) and when there is no change (contract `UpdateBatterySame`). Contract `UpdateBatteryChangeMode` says that if the battery is close to depletion (assertion `IsBatteryDegraded`) then special mode concerning ERT is applied: operation `SetBradycardiaMode`, which paces on ventricular chamber only without rate modulation. Assertion `IsBatteryDegraded` defines a predicate indicating whether the battery is ERT or ERP. The operation `UpdateBatteryStatus` takes into account assumption A8 (appendix B); we assume that the battery always degrades, and that the only way to take it out of the ventricular mode with rate modulation turned off is by reprogramming the device, which will also require that the battery is changed.




Figure 7.18: Contract diagrams of global operations `IncVPulseAmplReg` and `DecVPulseAmplReg`




Figure 7.19: Contract diagrams of global operations `IncAPulseAmplUReg` and `DecAPulseAmplUReg`




Figure 7.20: Contract diagrams of global operations `IncVPulseAmplUReg` and `DecVPulseAmplUReg`




Figure 7.21: Contract diagrams of global operations `IncAPulseWidth` and `DecAPulseWidth`




Figure 7.22: Contract diagrams of global operations `IncVPulseWidth` and `DecVPulseWidth`




Figure 7.23: Contract and assertion diagrams of global operation `IncHRL`




Figure 7.24: Contract and assertion diagrams of global operation **DecHRL**




Figure 7.25: Contract diagrams of global operations **APaceCommon** and **VPaceCommon**




Figure 7.26: Contract diagrams of global operation `UpdateCCIIfDoesRateSmoothing`




Figure 7.27: Contract diagrams of global operations `UpdateCCI` and `CurrentCCISame`




Figure 7.28: Contract diagram of global operation `SetTargetCCI`




Figure 7.29: Contract diagrams of global operation `ExecuteEvent`




Figure 7.30: Contract diagrams of global operations **TriggerAPace** and **TriggerVPace**




Figure 7.31: Contract diagrams of global operations **InhibitAPace** and **InhibitVPace**




Figure 7.32: Contract diagrams of global operations **TimedAPace** and **TimedVPace**




Figure 7.33: Assertion diagrams of global observe operations `TimeToAPace`, `TimeToVPace` and `TimeToAVDelay`




Figure 7.34: Contract diagram of global operations `TrackASense`

Chapter 8

The Sensing Package

8.1 State

The pacemaker device detects what is going on inside the heart by measuring the potential difference (voltage) between the two electrodes used for pacing [BSS10]. The package **Sensing** is responsible for this sensing functionality; it senses heart signals through the leads, receiving a sensed amplitude from the heart's sensors. In addition, it manages the *blanking period*: an interval of time during which the pacemaker is unable to sense (or does not respond to sensed events) [MM07].




Figure 8.1: VCL Package and structural diagram of package **Sensing**

PD of **Sensing** (Fig. 8.1(a)) package defines **Sensing** as an ensemble package; **Sensing** imports container package **CommonPM**.

Figure 8.1(b) presents the SD of package **Sensing**. It is as follows:

- Derived blob **Sensitivity** captures the allowed values of the sensitivity threshold parameter of the device, expressed in micro-volts (μV); the bounds of this set are taken from the “programmable parameters table” in [Bos07b, p. 34]. This parameter indicates the sensitivity level of sensing.




Figure 8.2: Behavioural Diagram of package **Sensing**

- Derived blob **Blanking** captures the allowed values for both atrial and ventricular blanking, as defined by the “Ventricular Blanking” interval the “programmable parameters table” in [Bos07b, p. 34]. This blob is related to requirements assumption A4 (appendix B).
- The Package blob comprises property edges to represent the current and previous heart chamber being sensed (**cs** and **pcs** of set **ChambersSensed**, which can have values atrial, ventricular or dual), the atrial and ventricular blanking periods (**aBlanking** and **vBlanking**), and the atrial and ventricular sensitivity (**aSensitivity** and **vSensitivity**).

8.2 Overall Behaviour

BD of package **Sensing** is given in Fig. 8.2. It introduces the following operations:

- **Init** (Fig. 8.3(a)) is the package’s initialisation .
- **SetBradycardiaMode** (Fig. 8.3(b)) is the operation that sets the bradycardia mode for sensing.
- **UpdateBatteryStatus** (Fig. 8.4) does the required action for sensing when the battery status is updated.
- **ExecuteEvent** (Fig. 8.5) executes some internal pacemaker event in the context of package **Sensing**.
- **IncASensitivity** (Fig. 8.6) and **DecASensitivity** (Fig. 8.7) increment and decrement the value of atrial sensitivity.
- **IncVSensitivity** (Fig. 8.8) and **DecVSensitivity** (Fig. 8.9) increment and decrement the value of atrial sensitivity.
- **IncABlanking** (Fig. 8.10) and **DecABlanking** (Fig. 8.11) increment and decrement the value of atrial blanking.
- **IncVBlanking** (Fig. 8.12) and **DecVBlanking** (Fig. 8.13) increment and decrement the value of ventricular blanking.
- **TriggerAPace** (Fig. 8.15(a)) and **TriggerVPace** (Fig. 8.15(b)) correspond to the triggered atrial and ventricular pacing. In the triggered mode of response to sensing, a sense in a

chamber shall trigger an immediate pace in that chamber [Bos07b, p. 17]. These operations are called when this occurs. Here, package **Sensing** describes the effect of this package on these operations (it adds further pre-conditions).

- **InhibitAPace** (Fig. 8.16(a)) and **InhibitVPace** (Fig. 8.16(b)) correspond to the inhibited atrial and ventricular pacing. In the inhibited mode of response to sensing, a sense in a chamber shall inhibit a pending pace in that chamber [Bos07b, p. 17]. These operations are called when this occurs. Here, package **Sensing** describes the effect of this package on these operations (it adds further pre-conditions).
- **TrackASense** (Fig. 8.16(c)) corresponds to the tracked atrial sensing. During tracked pacing, an atrial sense shall cause a tracked ventricular pace after a programmed AV delay unless a ventricular sense was detected beforehand [Bos07b, p. 17].
- **ASenseRefractory** (Fig. 8.17(a)) and **VSenseRefractory** (Fig. 8.17(b)) correspond to an atrial and ventricular sense during the atrial or ventricular refractory period.
- **ASenseOnly** (Fig. 8.14(a)) and **VSenseOnly** (Fig. 8.14(b)) correspond to a normal atrial and ventricular sense.

8.3 Operations




Figure 8.3: Package **Sensing**: Assertion diagram of Initialisation and contract diagram of operation **SetBradycardiaMode**




Figure 8.4: Contract Diagrams of global operation `UpdateBatteryStatus`

AD of initialisation (Fig. 8.3(a)) sets the properties of the package blob `Sensing` to the nominal values taken from the “programmable parameters table” [Bos07b, p. 34]. The only exception is the initial value of atrial blanking (not defined in [Bos07b, p. 34]); this is assumed to be as the ventricular blanking according to requirements assumption A4 (appendix B).

CD of operation `SetBradycardiaMode` (Fig. 8.3(b)) sets the property `cs` to the programmed value given as input (`cs?`).

The global operation `UpdateBatteryStatus` (Fig. 8.4) receives the new battery status as an input and considers two cases: (a) when the bradycardia mode needs to be changed because the battery becomes close to depletion (contract `UpdateBatteryChangeMode`), (b) and when there is no change (contract `UpdateBatterySame`). Contract `UpdateBatteryChangeMode` says that if the battery is close to depletion (assertion `IsBatteryDegraded`) then special mode concerning ERT is applied: only the ventricular chamber is paced and the rate modulation is turned off (postcondition). Assertion `IsBatteryDegraded` defines a predicate indicating whether the battery is ERT or ERP. The operation `UpdateBatteryStatus` takes into account assumption A8 (appendix B); we assume that the battery always degrades, and that the only way to take it out of the ventricular mode with rate modulation turned off is by reprogramming the device, which will also require that the battery is changed.

Figures 8.14(a) and 8.14(b) present the ADs of operations `ASenseOnly` and `VSenseOnly`, which are responsible for sensing in the atrial and ventricular chambers. Both operations receive an amplitude in millivolts (`aSAmp1?` and `aSVamp1?`), the current time in milliseconds (`now?`) and the time of the last atrial or ventricular event (`lastAEvTm?` and `lastVEvTm?`). There is a sense if the received amplitude is higher than the atrial or ventricular sensitivity (a sense signal has been detected), if the blanking period has elapsed, and if the mode allows sensing in the atrial or ventricular chambers. These operations make use of assumption A10 (appendix B), which says that sense signals are amplitudes.




Figure 8.5: Contract Diagrams of global operation **ExecuteEvent**




Figure 8.6: Contract Diagrams of global operation `IncASensitivity`




Figure 8.7: Contract Diagrams of global operation `DecASensitivity`




Figure 8.8: Contract Diagrams of global operation `IncVSensitivity`




Figure 8.9: Contract Diagrams of global operation `DecVSensitivity`




Figure 8.10: Contract Diagrams of global operation `IncABlanking`




Figure 8.11: Contract Diagrams of global operation **DecABlanking**




Figure 8.12: Contract Diagrams of global operation **IncVBlanking**




Figure 8.13: Contract Diagrams of global operation **DecVBlanking**




Figure 8.14: Assertion Diagram of predicates **ASenseOnly** and **VSenseOnly**




Figure 8.15: Assertion Diagrams of operations **TriggerAPace** and **TriggerVPace**




Figure 8.16: Assertion Diagrams of operations `InhibitAPace`, `InhibitVPace` and `TrackASense`




Figure 8.17: Assertion Diagrams of operations `ASenseRefractory` and `VSenseRefractory`

Chapter 9

The ResponseCommon Package

The `ResponseCommon` package introduces sets common across the packages that deal with the response-to-sensing functionality of the pacemaker. This is to be used by packages `AtrialResponse` and `VentricularResponse`.




Figure 9.1: VCL Package and structural diagram of package `ResponseCommon`

`ResponseCommon`'s PD (Fig. 9.1(a)) defines `ResponseCommon` as a container package. `ResponseCommon` imports package `CommonPM` in order to access the set `MilliSec`.

The package's SD (Fig. 9.1(b)) defines derived blob `RefractoryPeriod` from the blob `MilliSec` of package `CommonPm` according to the “programmable parameters table” of [Bos07b, p. 34]: an interval between 150 and 500 milliseconds. The *refractory period* refers to the period after either a sensed or paced beat in which sensing is disabled [MM07].

Chapter 10

The AtrialResponse Package

10.1 State

Package `AtrialResponse` is responsible for the response-to-sensing functionality of the pacemaker in what concerns the atrial chamber.




Figure 10.1: VCL Package diagram of package `AtrialResponse`

`AtrialResponse`'s PD (Fig. 10.1) defines `AtrialResponse` as an ensemble package. It imports packages `CommonPM` and `ResponseCommon`.

SD of package `AtrialResponse` (Fig. 10.2) introduces the following derived blobs:

- `UpperRateLimit` represents the allowed values of the *upper rate limit* (URL), expressed in ppm, and defined according to the “programmable parameters table” of [Bos07b, p. 34]. URL represents the maximum rate at which the paced ventricular rate will track atrial events [Bos07b, p. 29]. It represents the speed limit to control the response of the ventricular channel to sensed atrial activity [BSS10].
- `ATRFallbackTime` represents allowed values of *atrial tachycardia response (ATR) fallback time*, expressed in milliseconds (ms) and defined according to the “programmable parameters table” of [Bos07b, p. 34] (in the table the unit is expressed in minutes, which has been converted here to ms). A tachycardia is a rapid heart rate; it typically refers to a heart rate that exceeds the normal range for a resting heart rate (heart rate in an inactive or sleeping individual). ATR prevents long term pacing of a patient at unacceptably high rates during atrial tachycardia [Bos07b, p. 31]. When this happens, if the pacemaker is in the mode ATR, the Pacemaker should drop the current rate to the lower rate limit (LRL); the fallback time is the total time required to drop the rate to the LRL.




Figure 10.2: VCL structural diagram of package **AtrialResponse**

- **ATRDURATION** represents the allowed values of ATR duration; it is a set derived from the natural numbers, representing the number of cardiac cycles, and defined according to the “programmable parameters table” of [Bos07b, p. 34]. ATR duration is the state the Pacemaker device gets into after an atrial tachycardia is detected; it describes a delay (in terms of number of cardiac cycles) before the device enters ATR fallback [Bos07b, p. 31].
- **ATRDetectionCount** represents the allowed values for the ATR entry and exit detection counts, which are used to detect the beginning and the end of atrial tachycardia. This is defined according to requirements assumptions A11 and A12 (appendix B).

The primary blob **ATRStatus** contains values representing the possible status of ATR. The values are: **NoATR** (not responding to atrial tachycardia), **Duration** (ATR is in the *duration* state [Bos07b, p. 31]), **Fallback** (ATR is in the *fallback* state [Bos07b, p. 32]), **FallbackTmElapsed** (the fallback time has elapsed [Bos07b, p. 32]).

SD of Fig. 10.2 introduces package blob **AtrialResponse**, which has the property edges:

- **atrStatus** of blob **ATRStatus** indicates current status of ATR.
- **url**, **arp**, **atrMode**, **atrDuration**, **atrFbTm**, **atrEntryC** and **atrExitC** represent the current value for the programmable features URL [Bos07b, p. 29], atrial refractory period [Bos07b, p. 30], ATR mode (can be **On** or **Off**) [Bos07b, p. 31], ATR duration [Bos07b, p. 31], ATR fallback time [Bos07b, p. 32], and number of cardiac cycles required to detect the beginning and end of atrial tachycardia (assumptions A11 and A12, appendix B). The atrial refractory period is the programmed time interval following an atrial event during which time atrial events shall neither inhibit nor trigger pacing [Bos07b, p. 30]; available in single chamber atrial models only.

- `atrDetC`, `atrDurC`, `atrStartDetTm`, `atrStartFbTm` are used to ATR. `atrDetC` counts the current number of cycles detected for beginning or ending ATR. `atrDurC` counts the current number of cycles since duration started. `atrStartDetTm` indicates the starting time of ATR detection (either for beginning or ending ATR). `atrStartFbTm` indicates time when the device entered the fallback state.

10.2 Overall Behaviour




Figure 10.3: Behavioural Diagram of package `AtrialResponseResponse`

BD of package `AtrialResponse` is given in Fig. 10.3. It introduces the following behavioural units:

- `Init` is the package's initialisation (Fig. 10.4).
- `InATRFallBackMode` (Fig. 10.5) indicates whether the device is in ATR fallback mode.
- `ASenseOnly` (Fig. 10.13) and `ASenseRefractory` (Fig. 10.14) do, respectively, the normal atrial sensing and the the atrial sensing during the refractory period. These operations use the operation `PerformATR`, which is responsible for managing ATR (Fig. 10.6); this involves a complex processing according to what is described in [Bos07b, ps. 31–32]. The diagrams from figures 10.6 to 10.12 are involved in processing ATR.
- `IncURL` (Fig. 10.15) and `DecURL` (Fig. 10.16) increment and decrement the programmed value of URL.
- `IncARP` (Fig. 10.17) and `DecARP` (Fig. 10.18) increment and decrement the programmed value of atrial refractory period.
- `ATROn` (Fig. 10.19(a)) and `ATROff` (Fig. 10.19(b)) turn on and off the programmed value ATR mode.
- `IncATRDuration` (Fig. 10.20) and `DecATRDuration` (Fig. 10.21) increment and decrement the programmed value of ATR duration.
- `IncATRFbTm` (Fig. 10.22) and `DecATRFbTm` (Fig. 10.23) increment and decrement the ATR fallback time.
- `IncATREntryC` (Fig. 10.24) and `DecATREntryC` (Fig. 10.25) increment and decrement the ATR detection counter for beginning ATR.
- `IncATRExitC` (Fig. 10.26) and `DecATREntryC` (Fig. 10.27) increment and decrement the ATR detection counter for ceasing ATR.

10.3 Global Operations




Figure 10.4: Initialisation of package `AtrialResponse`




Figure 10.5: Assertion diagrams of operation `InATRFallBackMode` of package `AtrialResponse`




Figure 10.6: Contract diagrams of operation `PerformATR` of package `AtrialResponse`




Figure 10.7: Contract diagram of operation `ATRDetect` of package `AtrialResponse`




Figure 10.8: Contract diagrams of operation `ATRDetectEntry` of package `AtrialResponse`




Figure 10.9: Contract diagrams of operation `ATRDetectExit` of package `AtrialResponse`




Figure 10.10: Contract diagrams of operation `ATRDoDuration` of package `AtrialResponse`




Figure 10.11: Contract diagrams of operation `ATRFallback` of package `AtrialResponse`




Figure 10.12: Contract diagrams of operation `ATRDoNothing` of package `AtrialResponse`




Figure 10.13: Contract and assertion diagrams of operation `ASenseOnly` of package `AtrialResponse`




Figure 10.14: Contract diagram of operation `ASenseRefractory` of package `AtrialResponse`




Figure 10.15: Contract diagrams of operation `IncURL` of package `AtrialResponse`




Figure 10.16: Contract diagrams of operation `DecURL` of package `AtrialResponse`




Figure 10.17: Contract diagrams of operation **IncARP** of package **AtrialResponse**




Figure 10.18: Contract diagrams of operation **DecARP** of package **AtrialResponse**




Figure 10.19: Contract diagrams of operations `ATROn` and `ATROff` of package `AtrialResponse`




Figure 10.20: Contract diagrams of operations `IncATRDur` of package `AtrialResponse`




Figure 10.21: Contract diagrams of operations DecATRDur of package AtrialResponse




Figure 10.22: Contract diagrams of operations IncATRFbTm of package AtrialResponse




Figure 10.23: Contract diagrams of operations `DecATRFbTm` of package `AtrialResponse`




Figure 10.24: Contract diagrams of operations `IncATREntryC` of package `AtrialResponse`




Figure 10.25: Contract diagrams of operations `DecATREntryC` of package `AtrialResponse`




Figure 10.26: Contract diagrams of operations `IncATRExitC` of package `AtrialResponse`




Figure 10.27: Contract diagrams of operations **DecATRExitC** of package **AtrialResponse**

Chapter 11

The VentricularResponse Package

11.1 State

Package `VentricularResponse` is responsible for the response-to-sensing functionality of the pacemaker in what concerns the ventricular chamber.




Figure 11.1: VCL Package and structural diagram of package `VentricularResponse`

`VentricularResponse`'s PD (Fig. 11.1(a)) defines `VentricularResponse` as an ensemble package. It imports packages `CommonPM` and `ResponseCommon`.

SD of package `VentricularResponse` is given in Fig. 11.1(b). It is as follows:

- The package blob `VentricularResponse` introduces the property edges `vrp` and `pvarp` (of blob `RefractoryPeriod` defined in `ResponseCommon`), representing the *ventricular refractory period* (VRP) and the *post ventricular refractory period* (PVARP). VRP is the programmed time interval following a ventricular event (paced or sensed) during which ventricular senses shall not inhibit nor trigger pacing [Bos07b, p. 30]; its main goal is to eliminate sensing of the T wave and identifying that, incorrectly, as a ventricular event [MM07]. PVARP is the programmable time interval following a ventricular event, when an atrial cardiac cycle event shall not inhibit nor trigger a ventricular pace [Bos07b, p. 30]; this is the period when there is no sensing in the atrium [MM07], it avoids inappropriate atrial sensing of ventricular events and retrograde P waves [BSS10].

- Property edge `pvarpext` of package blob `VentricularResponse` represents the PVARP extension. Its target is the derived blob `PVARPExtension`. Blob `PVARPExtension` represents the allowed values for PVARP (post ventricular atrial refractory period) extension expressed in milliseconds (ms) and defined according to the “programmable parameters table” of [Bos07b, p. 34]; the value 0 represents the fact that PVARP extension is disabled. PVARP extension defines a value to be used in the place of the PVARP when premature ventricular contraction occurs [Bos07b, p. 30].
- Property edge `pvc` of package blob `VentricularResponse` indicates whether a premature ventricular contraction (PVC) has been detected and the PVARP extension needs to be applied.

11.2 Overall Behaviour




Figure 11.2: VCL behavioural diagram of package `VentricularResponse`

BD of package `AtrialResponse` is given in Fig. 11.2. It introduces the following behavioural units:

- `Init` (Fig. 11.3(a)) is the package’s initialisation.
- `GetPVARP` (Fig. 11.3(b)) gets the current value of the post-ventricular atrial period. This can be the normal PVARP or the extended PVARP if a PVC has been detected.
- `IsPVC` (Fig. 11.4) indicates whether a PVC has been detected or not.
- `IncVRP` (Fig. 11.5) and `DecVRP` (Fig. 11.6) increment and decrement the programmed value of the ventricular refractory period (VRP).
- `IncPVARP` (Fig. 11.7) and `DecPVARP` (Fig. 11.8) increment and decrement the programmed value of the post ventricular atrial refractory period (PVARP).
- `IncPVARPExt` (Fig. 11.9) and `DecPVARPExt` (Fig. 11.10) increment and decrement the programmed value of the extended post ventricular atrial refractory period (PVARP).
- `VSenseOnly` (Fig. 11.11(a)) is the operation doing the normal ventricular sensing.
- `VSenseRefractory` (Fig. 11.11(b)) is the operation doing the ventricular sensing during the refractory period.

11.3 Global Operations




Figure 11.3: Assertion diagrams of initialisation and `GetPVARP` of package `VentricularResponse`




Figure 11.4: Assertion diagrams of operation `IsPVC` of package `VentricularResponse`




Figure 11.5: Contract diagrams of operation `IncVRP` of package `VentricularResponse`




Figure 11.6: Contract diagrams of operation `DecVRP` of package `VentricularResponse`




Figure 11.7: Contract diagrams of operation `IncPVARP` of package `VentricularResponse`




Figure 11.8: Contract diagrams of operation `DecPVARP` of package `VentricularResponse`




Figure 11.9: Contract diagrams of operation `IncPVARPExt` of package `VentricularResponse`




Figure 11.10: Contract diagrams of operation `DecPVARPExt` of package `VentricularResponse`




Figure 11.11: Contract diagrams of operations `VSenseOnly` and `VSenseRefractory` of package `VentricularResponse`




Figure 11.12: Contract diagrams of operation `PVCUpdate` of package `VentricularResponse`

Chapter 12

The Response package

12.1 State

Package `AtrialResponse` is responsible for the response-to-sensing functionality of the pacemaker. It puts together packages `AtrialResponse` (chapter 10) and `VentricularResponse` (chapter 11).




Figure 12.1: VCL Package and structural diagram of package `VentricularResponse`

Response's PD (Fig. 12.1(a)) defines `Response` as an ensemble package. It imports packages `CommonPM`, `AtrialResponse` and `VentricularResponse`.

12.2 Overall Behaviour

BD of `Response` is given in Fig. 12.2. It includes two integral extensions to bring in operations from imported packages: one brings into the new context all operations from package `AtrialResponse`, except `Init`, `ASenseRefractory` and `ASenseOnly`; the other all operations from package `VentricularResponse`, except `Init`, `GetVARP` and `VSenseOnly`.

BD of package `Response` is given in Fig. 12.2. It introduces the following global operations:

- `Init` (Fig. 12.3(a)) is the package's initialisation.
- `SetBradycardiaMode` (Fig. 12.3(b)) is the operation that sets the bradycardia mode for response-to-sensing.




Figure 12.2: VCL behavioural diagram of package **Response**

- **UpdateBatteryStatus** (Fig. 12.4) does the required action for response-to-sensing when the battery status is updated.
- **ExecuteEvent** (Fig. 12.5) executes some internal pacemaker event in the context of response-to-sensing.
- **ASenseRefractory** (Fig. 12.6(a)) is the operation doing the atrial sensing during the refractory period.
- **ASenseOnly** (Fig. 12.6(b)) is the operation doing the normal atrial sensing.
- **VSenseOnly** (Fig. 12.7(a)) is the operation doing the normal ventricular sensing.
- **TrackASense** (Fig. 12.7(b)) corresponds to the tracked atrial sensing in response-to-sensing. During tracked pacing, an atrial sense shall cause a tracked ventricular pace after a programmed AV delay unless a ventricular sense was detected beforehand [Bos07b, p. 17].
- **TriggerAPace** (Fig. 12.8(a)) and **TriggerVPace** (Fig. 12.8(b)) correspond to the triggered atrial and ventricular pacing. In the triggered mode of response to sensing, a sense in a chamber shall trigger an immediate pace in that chamber [Bos07b, p. 17]. These operations are called when this occurs. Here, package **Response** describes the effect of this package on these operations (it adds further pre-conditions).
- **InhibitAPace** (Fig. 12.9(a)) and **InhibitVPace** (Fig. 12.9(b)) correspond to the inhibited atrial and ventricular pacing. In the inhibited mode of response to sensing, a sense in a chamber shall inhibit a pending pace in that chamber [Bos07b, p. 17]. These operations are called when this occurs. Here, package **response** describes the effect of this package on these operations (it adds further pre-conditions).

12.3 Global Operations




Figure 12.3: Package **response**: Assertion diagram of Initialisation and contract diagram of operation **SetBradycardiaMode**




Figure 12.4: Contract Diagrams of global operation **UpdateBatteryStatus**




Figure 12.5: Contract Diagrams of global operation **ExecuteEvent**




Figure 12.6: Contract diagrams of operations **ASenseRefractory** and **ASenseOnly** of package **Response**




Figure 12.7: Contract diagrams of operations **VSenseOnly** and **TrackASense** of package **Response**




Figure 12.8: Contract diagrams of operations **TriggerAPace** and **TriggerVPace** of package **Response**




Figure 12.9: Contract and assertion diagrams of operations **InhibitAPace** and **InhibitVPace** of package **Response**

Chapter 13

The HistoryCommon Package

The `HistoryCommon` package introduces common set for recording event markers (`History`).




Figure 13.1: VCL Package and structural diagram of package `HistoryCommon`

`HistoryCommon`'s PD (Fig. 13.1(a)) defines `HistoryCommon` as a container package.

The package's SD (Fig. 13.1(b)) defines enumerated blobs `AugmentationMarker`, `VentricularMarker` and `AtrialMarker`, which are defined according to sections 4.8.1, 4.8.2 and 4.8.3 of [Bos07b, p. 26].

Chapter 14

The History Package

14.1 State




Figure 14.1: VCL Package diagram of package **History**

History's PD (Fig. 14.1) defines **History** as an ensemble package, which imports container packages **CommonPM** and **HistoryPM**.




Figure 14.2: VCL structural diagram of package **History**

14.2 Overall Behaviour




Figure 14.3: VCL behavioural diagram of package **History**

The BD Fig. 14.3 defines several local operations:

- **New** of blob **HistoryEvent** (Fig. 14.4) an abstract operation that factors the commonality of **HistoryEvent** instances.
- **New** of blob **AtrialEvent** (Fig. 14.5) to create a new **AtrialEvent** instance.
- **New** of blob **VentricularEvent** (Fig. 14.6) to create a new **VentricularEvent** instance.

The following global operations are defined in the BD of Fig. 14.3:

- **CreateHistoryAEv** (Fig. 14.7) is the global operation to create new atrial history events.
- **CreateHistoryVEv** (Fig. 14.8) is the global operation to create new ventricular history events.

14.3 Local operations of blob HistoryEvent




Figure 14.4: Contract diagram of operation `New` of blob `HistoryEvent`

14.4 Local operations of blob AtrialEvent




Figure 14.5: Contract diagram of operation `New` of blob `AtrialEvent`

14.5 Local operations of blob VentricularEvent




Figure 14.6: Contract diagram of operation `New` of blob `VentricularEvent`

14.6 Global Operations




Figure 14.7: Contract diagram of operation `CreateHistoryAEv` of package `History`




Figure 14.8: Contract diagram of operation `CreateHistoryVEv` of package `History`

Chapter 15

The PacingWithHistoryI package

The package `PacingWithHistoryI` extends the package `Pacing` (chapter 7) with an interface to enable composition with the package `History` (chapter 14).

15.1 State



Figure 15.1: VCL Package diagram of package `PacingWithHistoryI`

The PD of `PacingWithHistoryI` is given in Fig. 15.1. It defines `PacingWithHistory` as an ensemble package, which imports ensemble package `Pacing`, and container packages `HistoryCommon` and `CommonPM`.

15.2 Overall Behaviour




Figure 15.2: VCL behavioural diagram of package `PacingWithHistory`

BD of package `PacingWithHistoryI` is given in Fig. 15.2. It includes one integral extension to bring in operations from package `Pacing`, which says that all operations of `Pacing` except `TimedAPace` and `TimedVPace` are imported integrally.

In addition, BD of Fig. 15.2 defines two operations defined in **Pacing** that are not defined integrally: **TimedAPace** (Fig. 15.4(a)) and **TimedVPace** (Fig. 15.4(b)). These operations extend **TimedAPace** and **TimedVPace** of **Pacing** with the creation of the proper history markers for those events.

15.3 Global Operations




Figure 15.3: Assertion diagrams of operations **SetEvMarkersA** and **SetEvMarkersV** of package **PacingWithHistory**




Figure 15.4: Contract diagrams of operations `TimedAPace` and `TimedVPACE` of package `PacingWithHistory`

Chapter 16

The ResponseWithHistoryI package

The package **ResponseWithHistoryI** extends the package **Response** (chapter 12) with an interface to enable composition with the package **History** (chapter 14).

16.1 State



Figure 16.1: VCL Package diagram of package **ResponseWithHistoryI**

The PD of **ResponseWithHistoryI** is given in Fig. 16.1. It defines **ResponseWithHistoryI** as an ensemble package, which imports ensemble package **Response**, and container packages **HistoryCommon** and **CommonPM**.

16.2 Overall Behaviour




Figure 16.2: VCL behavioural diagram of package **ResponseWithHistory**

BD of package **ResponseWithHistory** is given in Fig. 16.2. It includes one integral extension to bring in operations from package **Response**, which says that all operations of **Response** except

ASenseRefractory, **VSenseRefractory**, **ASenseOnly**, **VSenseOnly**, **TriggerAPace**, **TriggerVPace**, **InhibitAPace**, **InhibitVPace**, **TrackASense** and **TrackVSense** are imported integrally.

In addition, BD of Fig. 15.2 defines the operations defined in **Response** that are not defined integrally:

- **ASenseRefractory** (Fig. 16.6(a)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersA** (Fig. 16.3).
- **VSenseRefractory** (Fig. 16.6(b)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersV** (Fig. 16.5).
- **ASenseOnly** (Fig. 16.7(a)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersA** (Fig. 16.3).
- **VSenseOnly** (Fig. 16.7(b)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersV** (Fig. 16.5).
- **TriggerAPace** (Fig. 16.8(a)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersA** (Fig. 16.3).
- **TriggerVPace** (Fig. 16.8(b)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersV** (Fig. 16.5).
- **InhibitAPace** (Fig. 16.9(a)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersA** (Fig. 16.3).
- **InhibitVPace** (Fig. 16.9(b)) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersV** (Fig. 16.5).
- **TrackASense** (Fig. 16.10) adds the appropriate event markers for this operation of package **Response**; it uses operation **SetEvMarkersA** (Fig. 16.3).

16.3 Global Operations




Figure 16.3: Contract diagram of operation `SetEvMarkersA` of package `ResponseWithHistoryI`




Figure 16.4: Assertion diagrams of operations **ATMarker** and **ASMarker** of package **ResponseWithHistoryI**





Figure 16.5: Assertion diagrams of operation `SetEvMarkersV` of package `ResponseWithHistoryI`



(a) `ASenseRefractory`

(b) `VSenseRefractory`

Figure 16.6: Contract diagrams of operations `ASenseRefractory` and `VSenseRefractory` in package `ResponseWithHistoryI`




Figure 16.7: Contract diagrams of operations `ASenseOnly` and `VSenseOnly` in package `ResponseWithHistoryI`




Figure 16.8: Contract diagrams of operations `TriggerAPace` and `TriggerVPace` in package `ResponseWithHistoryI`




Figure 16.9: Contract diagrams of operations `InhibitAPace` and `InhibitVPace` in package `ResponseWithHistoryI`




Figure 16.10: Contract diagram of operations `TrackASense` in package `ResponseWithHistoryI`

Chapter 17

The PacemakerWithoutHistory package

The package `PacemakerWithoutHistory` puts together sensing, pacing, response, rate modulation and battery management. This will build a package that puts together all concerns except history.

17.1 State




Figure 17.1: VCL Package diagram of package `PacemakerWithoutHistory`

The PD of `PacemakerWithoutHistory` is given in Fig. 17.1. It defines `PacemakerWithoutHistory` as an ensemble package, which imports ensemble packages `Sensing`, `Battery`, `Pacing`, `RateModulation`, `ResponseWithHistoryI` and `PacingWithHistoryI`.

17.2 Overall Behaviour

BD of package `PacemakerWithoutHistory` is given in Fig. 17.2. It includes one merge box to merge all operations from packages `PacingWithHistoryI`, `Sensing`, `ResponseWithHistoryI`, `Battery` and `RateModulation`.




Figure 17.2: VCL behavioural diagram of package `PacemakerWithoutHistory`

Chapter 18

The Pacemaker package

The package **Pacemaker** puts together the concerns of sensing, pacing, response, rate modulation and battery management already composed in package **PacemakerWithoutHistory** with the history concern isolated in the **History** package.

18.1 State




Figure 18.1: VCL Package diagram of package **Pacemaker**

The PD of **Pacemaker** is given in Fig. 18.1. It defines **Pacemaker** as an ensemble package, which imports ensemble packages **PacemakerWithoutHistory** and **History**.

18.2 Overall Behaviour

BD of package **Pacemaker** is given in Fig. 18.2. It includes one integral extension to bring operations from package **PacemakerWithoutHistory**; this says that all operations from package **PacemakerWithoutHistory** are to be imported integrally except the operations mentioned on the list associated with the **All** operation.

The BD of Fig. 18.2 includes two join boxes to add history to the operations that have not been imported integrally. One join box joins history for the atrial operations and the other for the ventricular operations. These operations are joined through join contracts **HistoryAOp** (Fig. 18.3) and **HistoryVOp** (Fig. 18.4).




Figure 18.2: VCL behavioural diagram of package **Pacemaker**

18.3 Global Operations




Figure 18.3: Contract diagram of operation **HistoryAOp** of package **Pacemaker**




Figure 18.4: Contract diagram of operation `HistoryVOp` of package `Pacing`

Chapter 19

Snapshot Analysis

This chapter presents the snapshots that were used to analyse the pacemaker model presented in the previous chapter. This is done based on the snapshot analysis technique developed in [ASP04, Amá07].




Figure 19.1: Sample snapshots used in for model validation

Snapshot-pair of Fig. 19.1(a) simulates a change from the initial PM mode of DDD (dual chamber sensing and pacing) to VOO (asynchronous ventricular pacing), which is done through the operation `SetBradyCardiaMode`. The snapshot proofs for this pair are provable in the Z/Eves theorem prover, which means that the transition is valid as expected.

The snapshot sequence of Fig. 19.1(b) simulates a change of mode from DDD to DDDR, and then a situation when the battery status changes to ERT (Elective replacement is near), which means that the pacemaker needs to enter a degraded mode. The snapshot proofs for this snapshot sequence are provable in the Z/Eves theorem prover, which means that the sequence is valid as expected.

The snapshot sequence of Fig 19.2 simulates two tracked senses leading to the detection of a premature ventricular contraction (PVC). In the first snapshot(top), the `pvp` property is set to




Figure 19.2: Snapshot sequence simulating two tracked senses leading to a premature ventricular contraction

`false` in **Response**, meaning that PVC has not been detected. After the first **TrackVSense**, we get the second snapshot where the `pvc` flag is set to `true` because there has been a sense and the last sense happened after a pace(see AD `IsPVC` in Fig. 11.12; `vp.lastEvTm > ap.lastEvTm` in the first snapshot, which means that a pace did not follow a sense, as is usually expected). After the third **TrackVSense**, the flag is set to `false` as prescribed by the requirements.




Figure 19.3: Snapshot sequence simulating the detection of atrial tachycardia




Figure 19.4: Snapshot sequence simulating rate smoothing

References

- [AGK11] Nuno Amálio, Christian Glodt, and Pierre Kelsen. Building VCL models and automatically generating Z specifications from them. In *FM 2011*, LNCS. Springer, 2011.
- [AK10] Nuno Amálio and Pierre Kelsen. Modular design by contract visually and formally using VCL. In *VL/HCC 2010*, 2010.
- [AKMG10] Nuno Amálio, Pierre Kelsen, Qin Ma, and Christian Glodt. Using VCL as an aspect-oriented approach to requirements modelling. *TAOSD*, VII:151–199, 2010.
- [Amá07] Nuno Amálio. *Generative frameworks for rigorous model-driven development*. PhD thesis, Dept. Computer Science, Univ. of York, 2007.
- [APS05] Nuno Amálio, Fiona Polack, and Susan Stepney. An object-oriented structuring for Z based on views. In *ZB 2005*, volume 3455 of *LNCS*, pages 262–278. Springer, 2005.
- [ASP04] Nuno Amálio, Susan Stepney, and Fiona Polack. Formal proof from UML models. In *Proc. ICFEM 2004*, volume 3308 of *LNCS*, pages 418–433. Springer, 2004.
- [Bos] Boston Scientific. *ALTRUA 50 and ALTRUA 60: Multiprogrammable pacemakers*.
- [Bos07a] Boston Scientific. *Advanced Timing Cycles*, 2007.
- [Bos07b] Boston Scientific. Pacemaker system specification, 2007. <http://bit.ly/xWbNRM>.
- [Bos07c] Boston Scientific. *Pacing Codes and Modes Concepts*, 2007.
- [Bos07d] Boston Scientific. *Timing Cycles*, 2007.
- [BSS10] S. Serge Barold, Roland X. Stroobandt, and Alfons F. Sinnaeve. *Cardiac Pacemakers and Resynchronization Step-by-Step: an illustrated guide*. Wiley-Blackwell, 2010.
- [JIJ06] Raoul Jetley, S. Purushothaman Iyer, and Paul L. Jones. A formal methods approach to medical device review. *IEEE Computer*, 39(4), 2006.
- [JJA10] Paul Jones, Raoul Jetley, and Jay Abraham. A formal methods-based verification approach to medical device software analysis. *EE Times, Design*, 2010. <http://bit.ly/w5guUS>.
- [Lee11] Jérôme Leemans. Model a pacemaker system using VCL. Master’s thesis, FSTC, Univ of Luxembourg, 2011.
- [MM07] H. Weston Moses and James C. Mullin. *A practical guide to cardiac pacing*. Wolters Kluwer, 2007.

- [TLMB08] Jonathon Timperley, Paul Leeson, Andrew Mitchell, and Timothy Betts. *Cardiac Pacemakers and ICDs*. Oxford University Press, 2008.
- [WLBF09] Jim Woodcock, Peter Gorm Larsen, Juan Bicarregui, and John Fitzgerald. Formal methods: Practice and experience. *ACM Computing Surveys*, 41(4):19–36, 2009.

Appendix A

Requirements Issues

The requirements of the cardiac pacemaker device modelled here are defined informally in English in [Bos07b, p. 29]. During modelling, we found many issues with the requirements, namely, omission, ambiguities and inconsistencies. These issues are summarised in table A.1.

Table A.1: Issues found in the requirements document of the Pacemaker device [Bos07b].

ID	Description
RI1	The requirements document mentions programmable minimum and maximum AV delay [Bos07b, p. 29], but the “programmable parameters table” only defines the minimum “Dynamic AV Delay” and a ”Fixed AV Delay”; it does not define an interval for the allowed values of “Maximum AV Delay”.
RI2	The requirements document distinguishes between paced and sensed AV Delay [Bos07b, p. 29], but there is no distinction in “programmable parameters table” [Bos07b, p. 34].
RI3	The requirements document does not provide much information about a ”factor stored in device memory to create the dynamic AV delay” [Bos07b, p. 29]. Nowhere it says the value for this factor.
RI4	The requirements document misses many details regarding blanking. The requirements document only mentions “ventricular blanking” in the parameters table [Bos07b, p. 29] (it is the only mention of blanking in the document), but the requirements auxiliary documentation ¹ mention both atrial and ventricular blanking.
RI5	The requirements document is not very clear on how it gets the body motion from the accelerometer, which is required for rate adaptive pacing. The requirements [Bos07b, p. 32] just say: “the device shall have the ability to adjust the cardiac cycle in response to metabolic need as measured from body motion using an accelerometer.”
RI6	The requirements document is not very clear on how the rate adaptive pacing rate is calculated. No equations are provided; just a brief natural language description based on the <i>response factor</i> [Bos07b, p. 32], which is not clear to us.
RI7	The requirements document does not say what is the initial battery status. It identifies four possible status [Bos07b, p. 22], but it does not say which one is the initial one.

¹ Available at http://sqrl.mcmaster.ca/pacemaker_docs.htm

RI8	It is not clear what happens to the battery after it becomes depleted. If it is possible to charge it somehow, or that once depleted if it is possible that we get a status update to say that the battery becomes full again.
RI9	It is not clear how the system gets information about the battery status level of the device. The requirements [Bos07b, p. 32] just say: “the battery status for the device shall be used to predict the following battery status levels”. The battery status levels are: BOL (Beginning of Life), Elective Replacement Near (ERN), Elective Replacement Time (ERT) and Elective Replacement Past (ERP). No indication is provided on how the devices reaches these levels.
RI10	It is not clear how the device senses. All the requirements documents says about this is [Bos07b, p. 16]: “Rate Sensing shall be accomplished using bipolar electrodes and sensing circuits. All rate detection decisions shall be based in the measured cardiac cycle lengths of the sensed rhythm. Rate shall be evaluated on an interval-by-interval basis.” It does not say, therefore, how sensing signals are received by the device.
RI11	The atrial tachycardia (AT) detection algorithm of [Bos07b, p. 31] is not clear. There is no information on the number of intervals between atrial senses that are faster than the URL that need to be considered for starting AT treatment (“AT onset” in [Bos07b, p. 31]). There is no information on the number of intervals between atrial senses that are slower than the URL that need to be considered for ceasing AT treatment (“AT cessation” in [Bos07b, p. 31]). There is no information on the “detection period” for starting or ceasing AT treatment [Bos07b, p. 31].
RI12	The requirements document says that the pacemaker device should switch to VVIR mode when it enters the fallback mode of atrial tachycardia detection [Bos07b, p. 33, Sec. 5.6.3]. However, this is contradictory with the requirement that says that fallback should be terminated immediately if the device detects that atrial tachycardia is over [Bos07b, p. 33, Sec. 5.6.3]. It is contradictory because cessation of atrial tachycardia is detected over the atrial channel, which is turned off in the VVIR mode. This means that the device will not be able to exit the fallback mode: it is a dead end. The consequences of this is that the patient may no longer be suffering a tachycardia, but the pacemaker device keeps responding like the patient had one.

Appendix B

Requirements Assumptions

Appendix A documents the issues that we found in the requirements document of the cardiac pacemaker device [Bos07b] modelled here. This chapter documents the requirements assumptions that we had to take concerning these issues. The assumptions are listed in table B.1. The assumptions presented here are based on the challenge's support documentation released by Boston Scientific [Bos07c, Bos07d, Bos07a], on the specification of a device (similar to the one modelled here) produced by Boston Scientific [Bos], and information on DDD pacemakers that is available in the literature [BSS10, MM07, TLMB08].

Table B.1: Requirements assumptions that resolve issues with the requirements. First column gives assumption identifier; second column gives identifier of the requirements issue; and third column describes the assumption

A1	RI1	As the requirements document does define not an interval of allowed values for the maximum AV delay, we assume that it is the same as the "fixed AV Delay" defined in the "programmable parameters table" [Bos07b, p. 34]. The rationale for this is that the most logical thing is that the maximum AV delay programmable parameter must be there (DDD pacemakers seem to have it [BSS10] and another DDD pacemaker of Boston Scientific has it [Bos]) because it would not be possible to calculate dynamic AV delays without it. We take, therefore, the values of the maximum AV delay to be the values of 'fixed AV Delay'; from [Bos] these values seem reasonable.
A2	RI2	To resolve the ambiguity between the paced and sensed AV delay, we consider that there is just one AV Delay, which is updated whenever there is an atrial event (either sensed or paced).
A3	RI3	Since we don't know the value of the "factor stored in the device memory to create the dynamic AV delay" [Bos07b, p. 29], we create in the model a natural number constant to represent this factor, where the value of this constant is left undefined.
A4	RI4	To resolve the lack of information on blanking, we decided to follow the problems's auxiliary documentation ¹ , which considers them. We assume that the interval of allowed values for both atrial and ventricular blanking is as defined by "ventricular blanking" in the "programmable parameters table" [Bos07b, p. 34], which are programmable by the physician.

¹http://sqrl.mcmaster.ca/pacemaker_docs.htm

A5	RI5	To resolve the lack of information regarding how the accelerometer gets the body motion, we assume that the sensor will give a level of workload as defined by the activity threshold parameter of the “programmable parameters table” [Bos07b, p. 34].
A6	RI6	To resolve the lack of information regarding how rate adaptive pacing is calculated in [Bos07b, p. 34], we assume a linear equation based on LRL, MSR and response time when the rate needs to go up, and the LRL, MSR and recovery time when the rate needs to go down.
A7	RI7	We assume that the initial battery status is BOL (Beginning of Life). The device is initialized when the physician programs it; it would not make sense to put a battery in the device that is not in the state BOL. If the status is not correct then the operation to update the battery status is called to correct it.
A8	RI8	We assume that the battery updates tend towards degradation. We also assume that once the device is set to the special bradycardia mode after the battery reaches ERT (Elective Replacement Near), then it is not possible to change the mode unless the device is reprogrammed and the battery’s level is higher than ERT.
A9	RI9	We assume that the battery sensor indicates the battery status level when it changes. The battery status level provided by the sensor is as defined in the table “battery status and therapy availability” [Bos07b, p. 22, table 3]: BOL, ERN, ERT and ERP.
A10	RI10	We assume that sense signals are amplitudes. This is consistent with pacing signals, which are also amplitudes and with what the literature says on pacemakers [BSS10]; the requirements document also say [Bos07b, p. 19]: “The device shall be capable of storing pacing thresholds and P and R wave amplitude information.”
A11	RI11	We assume the programmable parameter ‘Entry Count’, which defines the number of atrial cycles above the URL that are required to initiate atrial tachycardia treatment. This is consistent with other Boston Scientific pacemakers that include a similar parameter, such as the Insignia [BSS10] and the ALTRUA [Bos]. The programmable values range between 1 and 8 cycles, with a nominal value of 8 (following the Insignia model [BSS10]).
A12	RI11	We assume the programmable parameter ‘Exit Count’, which defines the number of atrial cycles below the URL that are required to terminate atrial tachycardia treatment. This is consistent with other Boston Scientific pacemakers that include a similar parameter, such as the Insignia [BSS10] and the ALTRUA [Bos]. The programmable values range between 1 and 8 cycles, with a nominal value of 8 (following the Insignia model [BSS10]).
A13	RI11	Since the detection period for atrial tachycardia is not defined in [Bos07b], we introduce a constant (in milliseconds) in the model to represent this period, where the constant’s value is left undefined.
A14	RI12	Having VVIR as the fallback mode means that the fallback mode will be a dead end: once it enters fallback it will not be able to exit. To solve this, we correct the fallback mode to VDIR, which is the nominal fallback mode in the Insignia [Bos07b] and ALTRUA models [Bos] (in these models the fallback mode is a programmable parameter).

Appendix C

Generated Z Specification

C.1 Preamble

```
section Model_Preamble parents ZOO_Toolkit

CLASS ::= ChamberPacingCl | VentricularEventCl | AtrialEventCl | HistoryEventCl

subCl : CLASS ↔ CLASS
subCl = {(VentricularEventCl, HistoryEventCl), (AtrialEventCl, HistoryEventCl) }

∅ : CLASS → ℙ1 OBJ
∅x : CLASS → ℙ1 OBJ
disjoint ∅x

∀ cl : CLASS •
  ∅ cl = ∅x cl ∪ ∪(∅x ∥ (subCl+)~ (cl)) )

∀ cl, cl' : CLASS | cl ↦ cl' ∈ subCl •
  ∅ cl ⊆ ∅ cl'
```

C.2 Package CommonPM

```
section CommonPM parents ZOO_Toolkit, Model_Preamble

BradycardiaMode ::= O | R | A | V | D | T | I

Bool ::= True | False

Switch ::= On | Off
```

ChambersPaced == {*O, A, V, D*}

MicroSec == \mathbb{N}

MilliSec == \mathbb{N}

PPM == \mathbb{N}

MilliVolt == \mathbb{N}

MicroVolt == \mathbb{N}

ChambersSensed == {*O, A, V, D*}

ResponseToSensing == {*O, T, I, D*}

BatteryStatus ::= *BOL* | *ERN* | *ERT* | *ERP*

RateModulation == {*O, R*}

C.3 Package Battery

section *Battery* **parents** *ZOO_Toolkit*, *Model_Preamble*, *CommonPM*

BatteryGblSt

bs : *BatteryStatus*

BatteryGbl

BatteryGblSt

BatteryInit

BatteryGbl

bs = *BOL*

BatteryIsBatteryLevelHigh _____

<i>BatteryGbl</i>
<i>bs = BOL</i> \vee <i>bs = ERN</i>

BatteryUpdateBatteryStatus _____

Δ <i>BatteryGbl</i>
<i>bs? : BatteryStatus</i>
<i>bs' = bs?</i>

BatteryUpdateBatteryStatusUFr == *BatteryGbl \ (bs)*

BatteryUpdateBatteryStatusFr == Δ *BatteryGbl* \wedge \exists *BatteryUpdateBatteryStatusUFr*

BatteryVUpdateBatteryStatus _____

<i>BatteryUpdateBatteryStatusFr</i>
<i>BatteryUpdateBatteryStatus</i>

BatterySetBradycardiaMode _____

<i>BatteryGbl</i>
<i>BatteryIsBatteryLevelHigh</i>

BatteryVSetBradycardiaMode _____

<i>BatterySetBradycardiaMode</i>

C.4 Package RateSmoothing

section *RateSmoothing* **parents** *ZOO_Toolkit*, *Model_Preamble*, *CommonPM*

RateSmoothingValue == {*o* : \mathbb{N} | *o* \leq 25}

RateSmoothingGblSt _____

<i>rsUp : RateSmoothingValue</i>
<i>rsDown : RateSmoothingValue</i>

RateSmoothingGbl _____
 |
 | *RateSmoothingGblSt*

 | _____
 | *RateSmoothingInit* _____
 | | *RateSmoothingGblSt*
 | |
 | | *rsUp* = 0
 | | *rsDown* = 0

 | _____
 | *RateSmoothingIncRSDownMedium* _____
 | | *ΔRateSmoothingGblSt*

 | |
 | | *rsDown* ≤ 18
 | | *rsDown'* = *rsDown* + 3

 | _____
 | *RateSmoothingIncRSDownHigh* _____
 | | *ΔRateSmoothingGblSt*

 | |
 | | *rsDown* ≥ 21
 | | *rsDown'* = 25

 | _____
 | *RateSmoothingIncRSDown* _____
 | | *ΔRateSmoothingGbl*

 | | *RateSmoothingIncRSDownMedium* ∨ *RateSmoothingIncRSDownHigh*

RateSmoothingIncRSDownUFr == *RateSmoothingGbl* \ (*rsDown*)

RateSmoothingIncRSDownFr == *ΔRateSmoothingGbl* ∧ *ΞRateSmoothingIncRSDownUFr*

 | _____
 | *RateSmoothingVIncRSDown* _____
 | | *RateSmoothingIncRSDownFr*
 | | *RateSmoothingIncRSDown*

RateSmoothingDecRSDownHigh _____

$\Delta RateSmoothingGblSt$

$rsDown = 25$

$rsDown' = 21$

RateSmoothingDecRSDownMedium _____

$\Delta RateSmoothingGblSt$

$rsDown > 0 \wedge rsDown \leq 21$

$rsDown' = rsDown - 3$

RateSmoothingDecRSDownLow _____

$\Delta RateSmoothingGblSt$

$rsDown = 0$

$rsDown' = 0$

RateSmoothingDecRSDown _____

$\Delta RateSmoothingGbl$

$RateSmoothingDecRSDownHigh \vee RateSmoothingDecRSDownMedium \vee RateSmoothingDecRSDownLow$

$RateSmoothingDecRSDownUFr == RateSmoothingGbl \setminus (rsDown)$

$RateSmoothingDecRSDownFr == \Delta RateSmoothingGbl \wedge \exists RateSmoothingDecRSDownUFr$

RateSmoothingVDecRSDown _____

$RateSmoothingDecRSDownFr$

$RateSmoothingDecRSDown$

RateSmoothingIncRSUpMedium _____

$\Delta RateSmoothingGblSt$

$rsUp \leq 18$

$rsUp' = rsUp + 3$

$\text{RateSmoothingIncRSUpHigh}$

$\Delta \text{RateSmoothingGblSt}$

$rsUp \geq 21$

$rsUp' = 25$

$\text{RateSmoothingIncRSUp}$

$\Delta \text{RateSmoothingGbl}$

$\text{RateSmoothingIncRSUpMedium} \vee \text{RateSmoothingIncRSUpHigh}$

$\text{RateSmoothingIncRSUpUFr} == \text{RateSmoothingGbl} \setminus (rsUp)$

$\text{RateSmoothingIncRSUpFr} == \Delta \text{RateSmoothingGbl} \wedge \exists \text{RateSmoothingIncRSUpUFr}$

$\text{RateSmoothingVIncRSUp}$

$\text{RateSmoothingIncRSUpFr}$

$\text{RateSmoothingIncRSUp}$

$\text{RateSmoothingDecRSUpHigh}$

$\Delta \text{RateSmoothingGblSt}$

$rsUp = 25$

$rsUp' = 21$

$\text{RateSmoothingDecRSUpMedium}$

$\Delta \text{RateSmoothingGblSt}$

$rsUp > 0 \wedge rsUp \leq 21$

$rsUp' = rsUp - 3$

$\text{RateSmoothingDecRSUpLow}$

$\Delta \text{RateSmoothingGblSt}$

$rsUp = 0$

$rsUp' = 0$

$\text{RateSmoothingDecRSUp}$

$\Delta \text{RateSmoothingGbl}$

$\text{RateSmoothingDecRSUpHigh} \vee \text{RateSmoothingDecRSUpMedium} \vee \text{RateSmoothingDecRSUpLow}$

$\text{RateSmoothingDecRSUpUFr} == \text{RateSmoothingGbl} \setminus (\text{rsUp})$
 $\text{RateSmoothingDecRSUpFr} == \Delta \text{RateSmoothingGbl} \wedge \exists \text{RateSmoothingDecRSUpUFr}$

$\text{RateSmoothingVDecRSUp}$ $\text{RateSmoothingDecRSUpFr}$ $\text{RateSmoothingDecRSUp}$
--

$\text{RateSmoothingCalcUBSmoothedCCI}$ $\text{RateSmoothingGblSt}$ ubCCI! : MilliSec $\text{currentCCI? : MilliSec}$ $\text{ubCCI!} = (\text{currentCCI?} * (100 + \text{rsDown}) \text{ div } 100)$
--

$\text{RateSmoothingCalcLBSmoothedCCI}$ $\text{RateSmoothingGblSt}$ lbCCI! : MilliSec $\text{currentCCI? : MilliSec}$ $\text{lbCCI!} = (\text{currentCCI?} * (100 - \text{rsUp}) \text{ div } 100)$
--

$\text{RateSmoothingCalcSmoothedCCI0}$ RateSmoothingGbl $\text{targetCCI? : MilliSec}$ $\text{currentCCI? : MilliSec}$ $\text{newCCI! : MilliSec}$ $\text{RateSmoothingCalcUBSmoothedCCI}$ $\text{RateSmoothingCalcLBSmoothedCCI}$ $(\text{targetCCI?} > \text{currentCCI?} \wedge \text{targetCCI?} > \text{ubCCI!} \wedge \text{rsDown} > 0) \Rightarrow \text{newCCI!} = \text{ubCCI!}$ $(\text{targetCCI?} < \text{currentCCI?} \wedge \text{targetCCI?} < \text{lbCCI!} \wedge \text{rsUp} > 0) \Rightarrow \text{newCCI!} = \text{lbCCI!}$ $(\text{targetCCI?} > \text{currentCCI?} \wedge (\text{rsDown} = 0 \vee \text{targetCCI?} \leq \text{ubCCI!})) \Rightarrow \text{newCCI!} = \text{targetCCI?}$ $(\text{targetCCI?} \leq \text{currentCCI?} \wedge (\text{rsUp} = 0 \vee \text{targetCCI?} \geq \text{lbCCI!})) \Rightarrow \text{newCCI!} = \text{targetCCI?}$
--

$\text{RateSmoothingCalcSmoothedCCI} == \text{RateSmoothingCalcSmoothedCCI0} \setminus (\text{lbCCI!}, \text{ubCCI!})$

C.5 Package AVDelay

section *AVDelay parents ZOO_Toolkit, Model_Preamble, CommonPM*

| $dynDelayFactor : \mathbb{N}$

$MaxDynAVDelay == \{o : MilliSec \mid o \geq 70 \wedge o \leq 300\}$

$MinDynAVDelay == \{o : MilliSec \mid o = 30 \wedge o = 100\}$

$SensedAVDelayOffset == \{o : MilliSec \mid o \leq 100\}$

$AVDelayGblSt$

$dynMode : Switch$
 $fixedAVI : MaxDynAVDelay$
 $minDynDelay : MinDynAVDelay$
 $savDelayOffset : SensedAVDelayOffset$
 $maxDynDelay : MaxDynAVDelay$

$AVDelayGbl$

$AVDelayGblSt$

$AVDelayInit$

$AVDelayGbl$
 $fixedAVI = 150$
 $dynMode = Off$
 $minDynDelay = 50$
 $savDelayOffset = 0$
 $maxDynDelay = 150$

$AVDelayCalcDynDelayValue$

$AVDelayGblSt$
 $cdynAVI! : MilliSec$
 $cci? : MilliSec$
 $cdynAVI! = cci? \text{ div } dynDelayFactor$

AVDelayCalcDynDelay0

AVDelayGblSt
dynAVI! : MilliSec
cci? : MilliSec
AVDelayCalcDynDelayValue

cdynAVI! \geq *minDynDelay* \wedge *cdynAVI!* \leq *maxDynDelay* \Rightarrow *dynAVI!* = *cdynAVI!*
cdynAVI! < *minDynDelay* \Rightarrow *dynAVI!* = *minDynDelay*
cdynAVI! > *maxDynDelay* \Rightarrow *dynAVI!* = *maxDynDelay*

$$AVDelayCalcDynDelay == AVDelayCalcDynDelay0 \setminus (cdynAVI!)$$

AVDelayGetAVDelay0

AVDelayGbl
avi! : MilliSec
cci? : MilliSec
AVDelayCalcDynDelay

dynMode = *On* \Rightarrow *avi!* = *dynAVI!*
dynMode = *Off* \Rightarrow *avi!* = *fixedAVI*

$$AVDelayGetAVDelay == AVDelayGetAVDelay0 \setminus (dynAVI!)$$

AVDelayDynamicDelayOn

Δ *AVDelayGbl*

dynMode = *Off*
dynMode' = *On*

$$AVDelayDynamicDelayOnUFr == AVDelayGbl \setminus (dynMode)$$

$$AVDelayDynamicDelayOnFr == \Delta AVDelayGbl \wedge \exists AVDelayDynamicDelayOnUFr$$

AVDelayVDynamicDelayOn

AVDelayDynamicDelayOnFr
AVDelayDynamicDelayOn

$AVDelayDynamicDelayOff$ _____

$\Delta AVDelayGbl$

$dynMode = On$

$dynMode' = Off$

$AVDelayDynamicDelayOffUFr == AVDelayGbl \setminus (dynMode)$

$AVDelayDynamicDelayOffFr == \Delta AVDelayGbl \wedge \exists AVDelayDynamicDelayOffUFr$

$AVDelayVDynamicDelayOff$ _____

$AVDelayDynamicDelayOffFr$

$AVDelayDynamicDelayOff$

$AVDelayIncFixedAVIIInc$ _____

$\Delta AVDelayGblSt$

$fixedAVI \leq 290$

$fixedAVI' = fixedAVI + 10$

$AVDelayIncFixedAVISame$ _____

$\Delta AVDelayGblSt$

$fixedAVI > 290$

$fixedAVI' = 300$

$AVDelayIncFixedAVI$ _____

$\Delta AVDelayGbl$

$AVDelayIncFixedAVIIInc \vee AVDelayIncFixedAVISame$

$AVDelayIncFixedAVIUFr == AVDelayGbl \setminus (fixedAVI)$

$AVDelayIncFixedAVIFr == \Delta AVDelayGbl \wedge \exists AVDelayIncFixedAVIUFr$

$AVDelayVIncFixedAVI$ _____

$AVDelayIncFixedAVIFr$

$AVDelayIncFixedAVI$

$AVDelayDecFixedAVIDec$

$\Delta AVDelayGblSt$

$fixedAVI \geq 80$

$fixedAVI' = fixedAVI - 10$

$AVDelayDecFixedAVISame$

$\Delta AVDelayGblSt$

$fixedAVI < 80$

$fixedAVI' = 70$

$AVDelayDecFixedAVI$

$\Delta AVDelayGbl$

$AVDelayDecFixedAVIDec \vee AVDelayDecFixedAVISame$

$AVDelayDecFixedAVIUFr == AVDelayGbl \setminus (fixedAVI)$

$AVDelayDecFixedAVIFr == \Delta AVDelayGbl \wedge \exists AVDelayDecFixedAVIUFr$

$AVDelayVDecFixedAVI$

$AVDelayDecFixedAVIFr$

$AVDelayDecFixedAVI$

$AVDelayIncMinDynAVDelayInc$

$\Delta AVDelayGblSt$

$minDynDelay \leq 90$

$minDynDelay' = minDynDelay + 10$

$AVDelayIncMinDynAVDelaySame$

$\Delta AVDelayGblSt$

$minDynDelay > 90$

$minDynDelay' = 100$

$\overline{AVDelayIncMinDynAVDelay}$	$\Delta AVDelayGbl$
	$AVDelayIncMinDynAVDelayInc \vee AVDelayIncMinDynAVDelaySame$

$$AVDelayIncMinDynAVDelayUFr == AVDelayGbl \setminus (minDynDelay)$$

$$AVDelayIncMinDynAVDelayFr == \Delta AVDelayGbl \wedge \exists AVDelayIncMinDynAVDelayUFr$$

$\overline{AVDelayVIncMinDynAVDelay}$	$AVDelayIncMinDynAVDelayFr$
	$AVDelayIncMinDynAVDelay$

$\overline{AVDelayDecMinDynAVDelayDec}$	$\Delta AVDelayGblSt$
	$minDynDelay \geq 40$
	$minDynDelay' = minDynDelay - 10$

$\overline{AVDelayDecMinDynAVDelaySame}$	$\Delta AVDelayGblSt$
	$minDynDelay < 40$
	$minDynDelay' = 30$

$\overline{AVDelayDecMinDynAVDelay}$	$\Delta AVDelayGbl$
	$AVDelayDecMinDynAVDelayDec \vee AVDelayDecMinDynAVDelaySame$

$$AVDelayDecMinDynAVDelayUFr == AVDelayGbl \setminus (minDynDelay)$$

$$AVDelayDecMinDynAVDelayFr == \Delta AVDelayGbl \wedge \exists AVDelayDecMinDynAVDelayUFr$$

$\overline{AVDelayVDecMinDynAVDelay}$	$AVDelayDecMinDynAVDelayFr$
	$AVDelayDecMinDynAVDelay$

$\overline{AVDelayIncMaxDynAVDelayInc}$

$\Delta AVDelayGblSt$

$maxDynDelay \leq 290$

$maxDynDelay' = maxDynDelay + 10$

$\overline{AVDelayIncMaxDynAVDelaySame}$

$\Delta AVDelayGblSt$

$maxDynDelay > 290$

$maxDynDelay' = 300$

$\overline{AVDelayIncMaxDynAVDelay}$

$\Delta AVDelayGbl$

$AVDelayIncMaxDynAVDelayInc \vee AVDelayIncMaxDynAVDelaySame$

$AVDelayIncMaxDynAVDelayUFr == AVDelayGbl \setminus (maxDynDelay)$

$AVDelayIncMaxDynAVDelayFr == \Delta AVDelayGbl \wedge \exists AVDelayIncMaxDynAVDelayUFr$

$\overline{AVDelayVIncMaxDynAVDelay}$

$AVDelayIncMaxDynAVDelayFr$

$AVDelayIncMaxDynAVDelay$

$\overline{AVDelayDecMaxDynAVDelayDec}$

$\Delta AVDelayGblSt$

$maxDynDelay \geq 80$

$maxDynDelay' = maxDynDelay - 10$

$\overline{AVDelayDecMaxDynAVDelaySame}$

$\Delta AVDelayGblSt$

$maxDynDelay > 80$

$maxDynDelay' = 70$

$\overline{AVDelayDecMaxDynAVDelay}$

$\Delta AVDelayGbl$

$AVDelayDecMaxDynAVDelayDec \vee AVDelayDecMaxDynAVDelaySame$

$AVDelayDecMaxDynAVDelayUFr == AVDelayGbl \setminus (maxDynDelay)$
 $AVDelayDecMaxDynAVDelayFr == \Delta AVDelayGbl \wedge \exists AVDelayDecMaxDynAVDelayUFr$

$AVDelayVDecMaxDynAVDelay$ $\Delta AVDelayDecMaxDynAVDelayFr$ $AVDelayDecMaxDynAVDelay$

$AVDelayIncSAVDelayOffsetInc$ $\Delta AVDelayGblSt$ $savDelayOffset \leq 90$ $savDelayOffset' = savDelayOffset + 10$

$AVDelayIncSAVDelayOffsetSame$ $\Delta AVDelayGblSt$ $savDelayOffset > 90$ $savDelayOffset' = 100$

$AVDelayIncSAVDelayOffset$ $\Delta AVDelayGbl$ $AVDelayIncSAVDelayOffsetInc \vee AVDelayIncSAVDelayOffsetSame$
--

$AVDelayIncSAVDelayOffsetUFr == AVDelayGbl \setminus (savDelayOffset)$
 $AVDelayIncSAVDelayOffsetFr == \Delta AVDelayGbl \wedge \exists AVDelayIncSAVDelayOffsetUFr$

$AVDelayVIncSAVDelayOffset$ $\Delta AVDelayIncSAVDelayOffsetFr$ $AVDelayIncSAVDelayOffset$
--

$AVDelayDecSAVDelayOffsetDec$ $\Delta AVDelayGblSt$ $savDelayOffset \geq 10$ $savDelayOffset' = savDelayOffset - 10$

$AVDelayDecSAVDelayOffsetSame$	_____
$\Delta AVDelayGblSt$	_____
$savDelayOffset < 10$	_____
$savDelayOffset' = 0$	_____

$AVDelayDecSAVDelayOffset$	_____
$\Delta AVDelayGbl$	_____
$AVDelayDecSAVDelayOffsetDec \vee AVDelayDecSAVDelayOffsetSame$	_____

$AVDelayDecSAVDelayOffsetUFr == AVDelayGbl \setminus (savDelayOffset)$
 $AVDelayDecSAVDelayOffsetFr == \Delta AVDelayGbl \wedge \exists AVDelayDecSAVDelayOffsetUFr$

$AVDelayVDecSAVDelayOffset$	_____
$AVDelayDecSAVDelayOffsetFr$	_____
$AVDelayDecSAVDelayOffset$	_____

C.6 Package RateModulation

section *RateModulation* **parents** *ZOO_Toolkit*, *Model_Preamble*, *CommonPM*

MaximumSensorRate == { $o : PPM \mid o \geq 50 \wedge o \leq 175$ }

ResponseFactor == { $o : \mathbb{N} \mid o \geq 1 \wedge o \leq 16$ }

ReactionTime == { $o : MilliSec \mid o \geq 10000 \wedge o \leq 50000$ }

RecoveryTime == { $o : MilliSec \mid o \geq 120000 \wedge o \leq 960000$ }

Activity == { $o : \mathbb{N} \mid o \geq 1 \wedge o \leq 7$ }

LowerRateLimit == { $o : PPM \mid o \geq 30 \wedge o \leq 175$ }

RateModulationGblSt _____
 threshold : Activity
 msr : MaximumSensorRate
 rm : RateModulationMode
 responseF : ResponseFactor
 reactionTm : ReactionTime
 recoveryTm : RecoveryTime
 lastSwitchTm : MilliSec
 lastActivity : Activity
 prm : RateModulationMode
 targetCCI : MilliSec
 lrl : LowerRateLimit

RateModulationGbl _____
RateModulationGblSt

RateModulationInit _____
RateModulationGbl
 now? : MilliSec

 msr = 120
 responseF = 8
 lastSwitchTm = now?
 threshold = 4 \wedge lastActivity = 1
 reactionTm = 30000 \wedge recoveryTm = 600000
 rm = O \wedge prm = O
 lrl = 60 \wedge targetCCI = 60000 div lrl

RateModulationSetBradycardiaMode _____
 Δ *RateModulationGbl*
 rm? : RateModulationMode

 $rm' = rm?$
 $prm' = rm$

RateModulationSetBradycardiaModeUFr == *RateModulationGbl* \ (rm, prm)

RateModulationSetBradycardiaModeFr == Δ *RateModulationGbl* \wedge Ξ *RateModulationSetBradycardiaModeUFr*

RateModulationVSetBradycardiaMode

RateModulationSetBradycardiaModeFr

RateModulationSetBradycardiaMode

RateModulationIsAboveThreshold

RateModulationGbl

activity? : Activity

activity? > threshold

RateModulationCalcTargetCCISwitchingUp

RateModulationGbl

activity? : Activity

targetCCI! : MilliSec

RateModulationIsAboveThreshold \wedge \neg *RateModulationIsAboveThreshold*[*lastActivity/activity?*]

targetCCI! = 60000 div lrl

RateModulationCalcTargetCCISwitchingDown

RateModulationGbl

targetCCI! : MilliSec

activity? : Activity

\neg *RateModulationIsAboveThreshold* \wedge *RateModulationIsAboveThreshold*[*lastActivity/activity?*]

targetCCI! = 60000 div msr

RateModulationSetTargetCCISwitching0

Δ *RateModulationGbl*

now? : MilliSec

activity? : Activity

targetCCI! : MilliSec

RateModulationCalcTargetCCISwitchingUp \vee *RateModulationCalcTargetCCISwitchingDown*

lastSwitchTm' = now?

targetCCI' = targetCCI!

RateModulationSetTargetCCISwitching == *RateModulationSetTargetCCISwitching0* \ (*targetCCI!*)

RateModulationCalcTargetCCIIIncreasing _____

RateModulationGbl

activity? : Activity

now? : MilliSec

targetCCI! : MilliSec

$lastSwitchTm \leq now? - reactionTm \Rightarrow targetCCI! = 60000 \text{ div } msr$

$lastSwitchTm > now? - reactionTm \Rightarrow targetCCI! = (60000 \text{ div } lrl) + ((60000 * (now? - lastSwitchTm)) \text{ div } msr)$

RateModulationSetTargetCCIIIncreasing0 _____

$\Delta RateModulationGbl$

now? : MilliSec

activity? : Activity

targetCCI! : MilliSec

RateModulationCalcTargetCCIIIncreasing

targetCCI' = targetCCI!

lastSwitchTm' = lastSwitchTm

RateModulationSetTargetCCIIIncreasing == RateModulationSetTargetCCIIIncreasing0 \ (targetCCI!)

RateModulationCalcTargetCCIDecreasing _____

RateModulationGbl

activity? : Activity

now? : MilliSec

targetCCI! : MilliSec

$lastSwitchTm \leq now? - reactionTm \Rightarrow targetCCI! = 60000 \text{ div } lrl$

$lastSwitchTm > now? - recoveryTm \Rightarrow targetCCI! = (60000 \text{ div } msr) - ((60000 * (now? - lastSwitchTm)) \text{ div } msr)$

RateModulationSetTargetCCIDecreasing0 _____

$\Delta RateModulationGbl$

now? : MilliSec

activity? : Activity

targetCCI! : MilliSec

RateModulationCalcTargetCCIDecreasing

targetCCI' = targetCCI!

lastSwitchTm' = lastSwitchTm

$\text{RateModulationSetTargetCCIDecreasing} == \text{RateModulationSetTargetCCIDecreasing0} \setminus (\text{targetCCI}!)$

$\text{RateModulationSetTargetCCI}$ _____

$\Delta \text{RateModulationGbl}$
 $\text{now?} : \text{MilliSec}$
 $\text{activity?} : \text{Activity}$

$\text{RateModulationSetTargetCCISwitching} \vee \text{RateModulationSetTargetCCIIncreasing} \vee \text{RateModulationSetTargetCCIDecreasing}$

$\text{rm} = R$

$\text{lastActivity}' = \text{activity?}$

$\text{RateModulationSetTargetCCIUFr} == \text{RateModulationGbl} \setminus (\text{lastActivity}, \text{targetCCI}, \text{lastSwitchTm})$

$\text{RateModulationSetTargetCCIFr} == \Delta \text{RateModulationGbl} \wedge \exists \text{RateModulationSetTargetCCIUFr}$

$\text{RateModulationVSetTargetCCI}$ _____

$\text{RateModulationSetTargetCCIFr}$
 $\text{RateModulationSetTargetCCI}$

$\text{RateModulationSetTargetCCIToLRL}$ _____

$\Delta \text{RateModulationGbl}$

$\text{targetCCI}' = 60000 \text{ div } lrl$

$\text{RateModulationSetTargetCCIToLRLUFr} == \text{RateModulationGbl} \setminus (\text{targetCCI})$

$\text{RateModulationSetTargetCCIToLRLFr} == \Delta \text{RateModulationGbl} \wedge \exists \text{RateModulationSetTargetCCIToLRLUFr}$

$\text{RateModulationVSetTargetCCIToLRL}$ _____

$\text{RateModulationSetTargetCCIToLRLFr}$
 $\text{RateModulationSetTargetCCIToLRL}$

$\text{RateModulationIsLRLInc5}$ _____

$\text{RateModulationGblSt}$

$lrl < 50 \vee lrl \geq 90 \wedge lrl < 175$

RateModulationIncLRL5

$\Delta RateModulationGbl$

RateModulationIsLRLInc5

$lrl' = lrl + 5$

RateModulationIncLRL1

$\Delta RateModulationGbl$

$lrl \geq 50 \wedge lrl < 90$

$lrl' = lrl + 1$

RateModulationIncLRLSame

$\Delta RateModulationGbl$

$lrl = 175$

$lrl' = lrl$

RateModulationIncLRL

$\Delta RateModulationGbl$

RateModulationIncLRL5 \vee *RateModulationIncLRL1* \vee *RateModulationIncLRLSame*

RateModulationIncLRLUFr == *RateModulationGbl* \setminus (lrl)

RateModulationIncLRLFr == $\Delta RateModulationGbl \wedge \exists RateModulationIncLRLUFr$

RateModulationVIncLRL

RateModulationIncLRLFr

RateModulationIncLRL

RateModulationIsLRLDec5

RateModulationGblSt

$lrl > 90 \vee lrl \leq 50 \wedge lrl > 30$

RateModulationDeclRL5 _____

$\Delta RateModulationGbl$

RateModulationIsRLDec5

$lrl' = lrl - 5$

RateModulationDeclRL1 _____

$\Delta RateModulationGbl$

$lrl \leq 90 \wedge lrl > 50$

$lrl' = lrl - 1$

RateModulationDeclRLSame _____

$\Delta RateModulationGbl$

$lrl = 30$

$lrl' = lrl$

RateModulationDeclRL _____

$\Delta RateModulationGbl$

RateModulationDeclRL5 \vee *RateModulationDeclRL1* \vee *RateModulationDeclRLSame*

RateModulationDeclRLUFr == *RateModulationGbl* \ (lrl)

RateModulationDeclRLFr == $\Delta RateModulationGbl \wedge \exists RateModulationDeclRLUFr$

RateModulationVDeclRL _____

RateModulationDeclRLFr

RateModulationDeclRL

RateModulationIncMSRInc _____

$\Delta RateModulationGbl$

$msr < 175$

$msr' = msr + 5$

RateModulationIncMSRSame _____

$\Delta RateModulationGbl$

$msr = 175$

$msr' = msr$

$RateModulationIncMSR == RateModulationIncMSRInc$
 $\vee RateModulationIncMSRSame$

$RateModulationIncMSRUFr == RateModulationGbl \setminus (msr)$
 $RateModulationIncMSRFr == \Delta RateModulationGbl \wedge \exists RateModulationIncMSRUFr$

$RateModulationVIncMSR$ _____

 $RateModulationIncMSRFr$
 $RateModulationIncMSR$

$RateModulationDecMSRDec$ _____

 $\Delta RateModulationGbl$

 $msr > 50$

 $msr' = msr - 5$

$RateModulationDecMSRSame$ _____

 $\Delta RateModulationGbl$

 $msr = 50$

 $msr' = msr$

$RateModulationDecMSR == RateModulationDecMSRDec$
 $\vee RateModulationDecMSRSame$

$RateModulationDecMSRUFr == RateModulationGbl \setminus (msr)$
 $RateModulationDecMSRFr == \Delta RateModulationGbl \wedge \exists RateModulationDecMSRUFr$

$RateModulationVDecMSR$ _____

 $RateModulationDecMSRFr$
 $RateModulationDecMSR$

RateModulationIncRFInc _____

$\Delta RateModulationGbl$

$responseF < 16$

$responseF' = responseF + 1$

RateModulationIncRFSame _____

$\Delta RateModulationGbl$

$responseF = 16$

$responseF' = responseF$

RateModulationIncRF == RateModulationIncRFInc

$\vee RateModulationIncRFSame$

RateModulationIncRFUFr == RateModulationGbl \ (responseF)

RateModulationIncRFFr == ΔRateModulationGbl \ ∃ RateModulationIncRFUFr

RateModulationVIncRF _____

RateModulationIncRFFr

RateModulationIncRF

RateModulationDecRFDec _____

$\Delta RateModulationGbl$

$responseF > 1$

$responseF' = responseF - 1$

RateModulationDecRFSame _____

$\Delta RateModulationGbl$

$responseF = 1$

$responseF' = responseF$

RateModulationDecRF == RateModulationDecRFDec

$\vee RateModulationDecRFSame$

$\text{RateModulationDecRFUFr} == \text{RateModulationGbl} \setminus (\text{responseF})$
 $\text{RateModulationDecRFFr} == \Delta \text{RateModulationGbl} \wedge \exists \text{RateModulationDecRFUFr}$

$\text{RateModulationVDecRF} \quad \dots$
 $\text{RateModulationDecRFFr}$
 $\text{RateModulationDecRF}$

$\text{RateModulationIncReactionTmInc} \quad \dots$
 $\Delta \text{RateModulationGbl}$
 $\text{reactionTm} < 50000$
 $\text{reactionTm}' = \text{reactionTm} + 10000$

$\text{RateModulationIncReactionTmSame} \quad \dots$
 $\Delta \text{RateModulationGbl}$
 $\text{reactionTm} = 50000$
 $\text{reactionTm}' = \text{reactionTm}$

$\text{RateModulationIncReactionTm} == \text{RateModulationIncReactionTmInc}$
 $\vee \text{RateModulationIncReactionTmSame}$

$\text{RateModulationIncReactionTmUFr} == \text{RateModulationGbl} \setminus (\text{reactionTm})$
 $\text{RateModulationIncReactionTmFr} == \Delta \text{RateModulationGbl} \wedge \exists \text{RateModulationIncReactionTmUFr}$

$\text{RateModulationVIncReactionTm} \quad \dots$
 $\text{RateModulationIncReactionTmFr}$
 $\text{RateModulationIncReactionTm}$

$\text{RateModulationDecReactionTmDec} \quad \dots$
 $\Delta \text{RateModulationGbl}$
 $\text{reactionTm} > 10000$
 $\text{reactionTm}' = \text{reactionTm} - 10000$

$\neg RateModulationDecReactionTmSame$ _____

$\Delta RateModulationGbl$

$reactionTm = 10000$

$reactionTm' = reactionTm$

$RateModulationDecReactionTm == RateModulationDecReactionTmDec$

$\vee RateModulationDecReactionTmSame$

$RateModulationDecReactionTmUFr == RateModulationGbl \setminus (reactionTm)$

$RateModulationDecReactionTmFr == \Delta RateModulationGbl \wedge \exists RateModulationDecReactionTmUFr$

$\neg RateModulationVDecReactionTm$ _____

$RateModulationDecReactionTmFr$

$RateModulationDecReactionTm$

$\neg RateModulationIncRecoveryTmInc$ _____

$\Delta RateModulationGbl$

$recoveryTm < 960000$

$recoveryTm' = recoveryTm + 60000$

$\neg RateModulationIncRecoveryTmSame$ _____

$\Delta RateModulationGbl$

$recoveryTm = 960000$

$recoveryTm' = recoveryTm$

$RateModulationIncRecoveryTm == RateModulationIncRecoveryTmInc$

$\vee RateModulationIncRecoveryTmSame$

$RateModulationIncRecoveryTmUFr == RateModulationGbl \setminus (recoveryTm)$

$RateModulationIncRecoveryTmFr == \Delta RateModulationGbl \wedge \exists RateModulationIncRecoveryTmUFr$

RateModulationVIncRecoveryTm _____
 ┌── *RateModulationIncRecoveryTmFr*
 └── *RateModulationIncRecoveryTm*

RateModulationDecRecoveryTmDec _____
 ┌── $\Delta \text{RateModulationGbl}$
 └── $\text{recoveryTm} > 120000$
 $\text{recoveryTm}' = \text{recoveryTm} - 60000$

RateModulationDecRecoveryTmSame _____
 ┌── $\Delta \text{RateModulationGbl}$
 └── $\text{recoveryTm} = 120000$
 $\text{recoveryTm}' = \text{recoveryTm}$

RateModulationDecRecoveryTm == *RateModulationDecRecoveryTmDec*
 \vee *RateModulationDecRecoveryTmSame*

RateModulationDecRecoveryTmUFr == *RateModulationGbl* \ (*recoveryTm*)
RateModulationDecRecoveryTmFr == $\Delta \text{RateModulationGbl} \wedge \exists \text{RateModulationDecRecoveryTmUFr}$

RateModulationVDecRecoveryTm _____
 ┌── *RateModulationDecRecoveryTmFr*
 └── *RateModulationDecRecoveryTm*

RateModulationIncActThresholdInc _____
 ┌── $\Delta \text{RateModulationGbl}$
 └── $\text{threshold} < 7$
 $\text{threshold}' = \text{threshold} + 1$

$\neg RateModulationIncActThresholdSame$ _____

$\Delta RateModulationGbl$

$\overline{threshold} = 7$

$threshold' = threshold$

$RateModulationIncActThreshold == RateModulationIncActThresholdInc$

$\vee RateModulationIncActThresholdSame$

$RateModulationIncActThresholdUFr == RateModulationGbl \setminus (threshold)$

$RateModulationIncActThresholdFr == \Delta RateModulationGbl \wedge \exists RateModulationIncActThresholdUFr$

$\neg RateModulationVIncActThreshold$ _____

$RateModulationIncActThresholdFr$

$RateModulationIncActThreshold$

$\neg RateModulationDecActThresholdDec$ _____

$\Delta RateModulationGbl$

$\overline{threshold} > 1$

$threshold' = threshold - 1$

$\neg RateModulationDecActThresholdSame$ _____

$\Delta RateModulationGbl$

$\overline{threshold} = 1$

$threshold' = threshold$

$RateModulationDecActThreshold == RateModulationDecActThresholdDec$

$\vee RateModulationDecActThresholdSame$

$RateModulationDecActThresholdUFr == RateModulationGbl \setminus (threshold)$

$RateModulationDecActThresholdFr == \Delta RateModulationGbl \wedge \exists RateModulationDecActThresholdUFr$

$\neg RateModulationVDecActThreshold$ _____

$RateModulationDecActThresholdFr$

$RateModulationDecActThreshold$

C.7 Package Pacing

section *Pacing* **parents** *ZOO_Toolkit*, *Model_Preamble*, *CommonPM*, *RateSmoothing*, *AVDelay*, *RateModulator*

PulseWidth == {*o* : *MicroSec* | *o* ≥ 50 ∧ *o* ≤ 1900}

PulseAmplitudeMode ::= *Regulated* | *Unregulated*

PulseAmplitudeUnregulated == {*o* : *MilliVolt* | *o* ≥ 1250 ∧ *o* = 5000}

PulseAmplitudeRegulated == {*o* : *MilliVolt* | *o* ≥ 500 ∧ *o* ≤ 7000}

ChamberPacing0
—
amplMode : *PulseAmplitudeMode*
pulseWidth : *PulseWidth*
amplReg : *PulseAmplitudeRegulated*
amplUReg : *PulseAmplitudeUnregulated*
lastEvTm : *MilliSec*

ChamberPacing
—
ChamberPacing0

SChamberPacing
—
sChamberPacing : $\mathbb{P}(\mathbb{O} \text{ ChamberPacingCl})$
stChamberPacing : $\mathbb{O} \text{ ChamberPacingCl} \rightarrow \text{ChamberPacing}$
dom *stChamberPacing* = *sChamberPacing*

SChamberPacingInit
—
SChamberPacing
—
sChamberPacing = \emptyset ∧ *stChamberPacing* = \emptyset

ChamberPacingNew _____

Δ *ChamberPacing*
now? : MilliSec

amplMode = *Unregulated*
amplReg = 3500
amplUReg = 3750
pulseWidth = 400
lastEvTm = *now?*

ChamberPacingSetLastEventTm _____

Δ *ChamberPacing*
now? : MilliSec

amplMode' = *amplMode'*
amplReg' = *amplReg'*
amplUReg' = *amplUReg'*
pulseWidth' = *pulseWidth*
lastEvTm' = *now?*

ChamberPacingIncPulseAmplRegLow _____

Δ *ChamberPacing*

amplReg < 3500
amplReg' = *amplReg* + 100
amplMode' = *amplMode*
amplUReg' = *amplUReg*
pulseWidth' = *pulseWidth*
lastEvTm' = *lastEvTm*

ChamberPacingIncPulseAmplRegMedium

$\Delta_{ChamberPacing}$

$amplReg \geq 3500 \wedge amplReg < 7000$

$amplReg' = amplReg + 500$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseAmplRegSame

$\Delta_{ChamberPacing}$

$amplReg = 7000$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseAmplReg

$\Delta_{ChamberPacing}$

ChamberPacingIncPulseAmplRegLow \vee *ChamberPacingIncPulseAmplRegMedium* \vee *ChamberPacingIncPulseAmplRegHigh*

ChamberPacingDecPulseAmplRegLow

$\Delta_{ChamberPacing}$

$amplReg > 500 \wedge amplReg \leq 3500$

$amplReg' = amplReg - 100$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseAmplRegMedium

$\Delta_{ChamberPacing}$

$amplReg > 3500 \wedge amplReg \leq 7000$

$amplReg' = amplReg - 500$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseAmplRegSame

$\Delta_{ChamberPacing}$

$amplReg = 500$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseAmplReg

$\Delta_{ChamberPacing}$

$ChamberPacingDecPulseAmplRegLow \vee ChamberPacingDecPulseAmplRegMedium \vee ChamberPacingDecPulseAmplRegHigh$

ChamberPacingIncPulseAmplURegInc

$\Delta_{ChamberPacing}$

$amplUReg < 5000$

$amplUReg' = amplUReg + 1250$

$amplMode' = amplMode$

$amplReg' = amplReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseAmplURegSame _____

$\Delta_{ChamberPacing}$

$amplUReg = 5000$

$amplUReg' = amplUReg$

$amplMode' = amplMode$

$amplReg' = amplReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseAmplUReg _____

$\Delta_{ChamberPacing}$

ChamberPacingIncPulseAmplURegInc \vee *ChamberPacingIncPulseAmplURegSame*

ChamberPacingDecPulseAmplURegDec _____

$\Delta_{ChamberPacing}$

$amplUReg > 1250$

$amplUReg' = amplUReg - 1250$

$amplMode' = amplMode$

$amplReg' = amplReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseAmplURegSame _____

$\Delta_{ChamberPacing}$

$amplUReg = 1250$

$amplUReg' = amplUReg$

$amplMode' = amplMode$

$amplReg' = amplReg$

$pulseWidth' = pulseWidth$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseAmplUReg _____

$\Delta_{ChamberPacing}$

ChamberPacingDecPulseAmplURegDec \vee *ChamberPacingDecPulseAmplURegSame*

ChamberPacingIncPulseWidthLow _____

$\Delta_{ChamberPacing}$

$pulseWidth < 100$

$pulseWidth' = pulseWidth + 50$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseWidthMedium _____

$\Delta_{ChamberPacing}$

$pulseWidth \geq 100 \wedge pulseWidth < 1900$

$pulseWidth' = pulseWidth + 100$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseWidthSame _____

$\Delta_{ChamberPacing}$

$pulseWidth = 1900$

$pulseWidth' = pulseWidth$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$lastEvTm' = lastEvTm$

ChamberPacingIncPulseWidth _____

$\Delta_{ChamberPacing}$

ChamberPacingIncPulseWidthLow \vee *ChamberPacingIncPulseWidthMedium* \vee *ChamberPacingIncPulseWidthSame*

ChamberPacingDecPulseWidthHigh _____

$\Delta_{ChamberPacing}$

$pulseWidth \leq 1900 \wedge pulseWidth > 100$

$pulseWidth' = pulseWidth - 100$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseWidthLow _____

$\Delta_{ChamberPacing}$

$pulseWidth \leq 100$

$pulseWidth' = pulseWidth - 50$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseWidthSame _____

$\Delta_{ChamberPacing}$

$pulseWidth = 50$

$pulseWidth' = pulseWidth$

$amplReg' = amplReg$

$amplMode' = amplMode$

$amplUReg' = amplUReg$

$lastEvTm' = lastEvTm$

ChamberPacingDecPulseWidth _____

$\Delta_{ChamberPacing}$

$ChamberPacingDecPulseWidthHigh \vee ChamberPacingDecPulseWidthLow \vee ChamberPacingDecPulseWidthSame$

ChamberPacingGetPulseAmpl _____

ChamberPacing
pulseAmpl! : *MilliVolt*

amplMode = *Regulated* \Rightarrow *pulseAmpl!* = *amplReg*

amplMode = *Unregulated* \Rightarrow *pulseAmpl!* = *amplUReg*

ChamberPacingGetPulseWidth _____

ChamberPacing
pulseWidth! : *MicroSec*

pulseWidth! = *pulseWidth*

ChamberPacingGetLastEventTm _____

ChamberPacing
lastEvTm! : *MilliSec*

lastEvTm! = *lastEvTm*

SChamberPacingNF _____

$\Delta S_{ChamberPacing}$

ChamberPacing

o! : $\odot ChamberPacingCl$

o! $\in \odot x_{ChamberPacingCl} \setminus s_{ChamberPacing}$

sChamberPacing' = *sChamberPacing'* $\cup \{o!\}$

stChamberPacing' = *stChamberPacing'* $\cup \{(o!, \theta ChamberPacing)\}$

SChamberPacingUF _____

$\Delta S_{ChamberPacing}$

$\Delta ChamberPacing$

o? : $\odot x ChamberPacingCl$

o? $\in s_{ChamberPacing}$

$\theta ChamberPacing = st_{ChamberPacing} o?$

sChamberPacing' = *sChamberPacing*

stChamberPacing' = *stChamberPacing'* $\oplus \{(o?, \theta ChamberPacing')\}$

<i>SChamberPacingOF</i>	_____
<i>SChamberPacing</i>	
<i>ChamberPacing</i>	
<i>o? : O ChamberPacingCl</i>	
<i>o? ∈ sChamberPacing</i>	
<i>θ ChamberPacing = stChamberPacing o?</i>	

SChamberPacingNew == $\exists \text{ ChamberPacing} \bullet$
SChamberPacingNF[cp!/o!] ∧ ChamberPacingNew

SChamberPacingSetLastEventTm == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingSetLastEventTm

SChamberPacingIncPulseAmplReg == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingIncPulseAmplReg

SChamberPacingDecPulseAmplReg == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingDecPulseAmplReg

SChamberPacingIncPulseAmplUReg == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingIncPulseAmplUReg

SChamberPacingDecPulseAmplUReg == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingDecPulseAmplUReg

SChamberPacingIncPulseWidth == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingIncPulseWidth

SChamberPacingDecPulseWidth == $\exists \Delta \text{ChamberPacing} \bullet$
SChamberPacingUF ∧ ChamberPacingDecPulseWidth

SChamberPacingGetLastEventTm == $\exists \text{ ChamberPacing} \bullet$
SChamberPacingOF ∧ ChamberPacingGetLastEventTm

SChamberPacingGetPulseWidth == $\exists \text{ ChamberPacing} \bullet$
SChamberPacingOF ∧ ChamberPacingGetPulseWidth

SChamberPacingGetPulseAmpl == $\exists \text{ ChamberPacing} \bullet$
SChamberPacingOF ∧ ChamberPacingGetPulseAmpl

PacingPkgSt

hystMode : *Switch*
hrl : *LowerRateLimit*
aP : \textcircled{O} *ChamberPacingCl*
vP : \textcircled{O} *ChamberPacingCl*
cp : *ChambersPaced*
currentCCI : *MilliSec*
pcp : *ChambersPaced*

PacingOnlyGblSt

PacingPkgSt
SChamberPacing

PacingGblSt

PacingOnlyGblSt
AVDelayGbl
RateSmoothingGbl
RateModulationGbl

PacingTargetCCIIInv

PacingGblSt

$rm = O \Rightarrow targetCCI \leq 60000 \text{ div } lrl$
 $(rm = O \wedge hystMode = On) \Rightarrow targetCCI = 60000 \text{ div } hrl$
 $(rm = R \wedge hystMode = Off) \Rightarrow targetCCI = 60000 \text{ div } lrl$

PacingGbl

PacingGblSt

PacingTargetCCIIInv

PacingIntExtRateSmoothingFr == $\Delta PacingGbl \wedge \exists PacingOnlyGblSt \wedge \exists AVDelayGbl \wedge \exists RateModulationGbl$

PacingVIncRSDown

PacingIntExtRateSmoothingFr
RateSmoothingVIncRSDown

PacingVDecRSDown

PacingIntExtRateSmoothingFr
RateSmoothingVDecRSDown

PacingVIIncRSUp

PacingIntExtRateSmoothingFr
RateSmoothingVIIncRSUp

PacingVDecRSUp

PacingIntExtRateSmoothingFr
RateSmoothingVIIncRSUp

PacingIntExtAvDelayFr == ΔPacingGbl ∧ ∃PacingOnlyGblSt ∧ ∃RateSmoothingGbl ∧ ∃RateModulationGbl

PacingVDynamicDelayOn

PacingIntExtAvDelayFr
AVDelayVDynamicDelayOn

PacingVDynamicDelayOff

PacingIntExtAvDelayFr
AVDelayVDynamicDelayOff

PacingVIIncFixedAVI

PacingIntExtAvDelayFr
AVDelayIncFixedAVI

PacingVDecFixedAVI

PacingIntExtAvDelayFr
AVDelayVDecFixedAVI

$PacingVI_{Inc}MinDynAVDelay$ _____
 └── $PacingIntExtAvDelayFr$
 $AVDelayVI_{Inc}MinDynAVDelay$

$PacingVDecMinDynAVDelay$ _____
 └── $PacingIntExtAvDelayFr$
 $AVDelayVDecMinDynAVDelay$

$PacingVI_{Inc}MaxDynAVDelay$ _____
 └── $PacingIntExtAvDelayFr$
 $AVDelayVI_{Inc}MaxDynAVDelay$

$PacingVDecMaxDynAVDelay$ _____
 └── $PacingIntExtAvDelayFr$
 $AVDelayVDecMaxDynAVDelay$

$PacingVI_{Inc}SAVDelayOffset$ _____
 └── $PacingIntExtAvDelayFr$
 $AVDelayVI_{Inc}SAVDelayOffset$

$PacingVDecSAVDelayOffset$ _____
 └── $PacingIntExtAvDelayFr$
 $AVDelayVDecSAVDelayOffset$

$PacingIntExtRateModulationFr == \Delta PacingGbl \wedge \exists PacingOnlyGblSt \wedge \exists RateSmoothingGbl \wedge \exists AVDelayGbl$

$PacingVSetTargetCCI$ _____
 └── $PacingIntExtRateModulationFr$
 $RateModulationVSetTargetCCI$

$PacingVI_{Inc}LRL$ _____
 └── $PacingIntExtRateModulationFr$
 $RateModulationVI_{Inc}LRL$

PacingVDecLRL

PacingIntExtRateModulationFr

RateModulationVDecLRL

PacingVIIncMSR

PacingIntExtRateModulationFr

RateModulationVIIncMSR

PacingVDecMSR

PacingIntExtRateModulationFr

RateModulationVDecMSR

PacingVIIncRF

PacingIntExtRateModulationFr

RateModulationVIIncRF

PacingVDecRF

PacingIntExtRateModulationFr

RateModulationVDecRF

PacingVIIncReactionTm

PacingIntExtRateModulationFr

RateModulationVIIncReactionTm

PacingVDecReactionTm

PacingIntExtRateModulationFr

RateModulationVDecReactionTm

PacingVIIncRecoveryTm

PacingIntExtRateModulationFr

RateModulationVIIncRecoveryTm

PacingVDecRecoveryTm
PacingIntExtRateModulationFr
RateModulationVDecRecoveryTm

PacingVIIncActThreshold
PacingIntExtRateModulationFr
RateModulationVIIncActThreshold

PacingVDecActThreshold
PacingIntExtRateModulationFr
RateModulationVDecActThreshold

PacingInit0
PacingGbl
SChamberPacingInit

PacingInit
 Δ PacingGbl
PacingInit0
AVDelayInit
RateSmoothingInit
RateModulationInit
now? : MilliSec
SChamberPacingNew[cpa!/cp!] ; SChamberPacingNew[cpv!/cp!]
 $hystMode' = Off$
 $hrl' = 60$
 $currentCCI' = targetCCI$
 $cp' = D$
 $aP' = cpa!$
 $vP' = cpv!$

PacingSetBradycardiaMode —————

$\Delta PacingGbl$

$cp? : ChambersPaced$

$RateModulationSetBradycardiaMode$

$cp' = cp?$

$pcp' = pcp$

$PacingSetBradycardiaModeUFr == PacingGbl \setminus (cp, pcp, rm, prm)$

$PacingSetBradycardiaModeFr == \Delta PacingGbl \wedge \exists PacingSetBradycardiaModeUFr$

PacingVSetBradycardiaMode —————

$PacingSetBradycardiaModeFr$

$PacingSetBradycardiaMode$

PacingIsBatteryDegraded —————

$\Delta PacingGblSt$

$bs? : BatteryStatus$

$bs? = ERT \vee bs? = ERP$

PacingUpdateBatteryChangeMode —————

$\Delta PacingGblSt$

$PacingIsBatteryDegraded$

$PacingSetBradycardiaMode[V / cp?, R / rm?]$

PacingUpdateBatterySame —————

$\Delta PacingGblSt$

$cp' = cp$

$rm' = rm$

$prm' = prm$

$pcp' = pcp$

PacingUpdateBatteryStatus —————

$\Delta PacingGbl$

$bs? : BatteryStatus$

$PacingUpdateBatteryChangeMode \vee PacingUpdateBatterySame$

PacingUpdateBatteryStatusUFr == PacingGbl \ (cp, pcp, rm, prm)

PacingUpdateBatteryStatusFr == ΔPacingGbl ∧ ∃PacingUpdateBatteryStatusUFr

PacingVUpdateBatteryStatus

PacingUpdateBatteryStatusFr

PacingUpdateBatteryStatus

PacingIncAPulseAmplReg

ΔPacingGbl

SChamberPacingIncPulseAmplReg[aP/o?]

PacingVIncAPulseAmplReg

ΔPacingGbl

ΞPacingPkgSt

ΞAVDelayGbl

ΞRateSmoothingGbl

PacingIncAPulseAmplReg

PacingDecAPulseAmplReg

ΔPacingGbl

SChamberPacingDecPulseAmplReg[aP/o?]

PacingVDecAPulseAmplReg

ΔPacingGbl

ΞPacingOnlyGblSt

ΞAVDelayGbl

ΞRateSmoothingGbl

PacingDecAPulseAmplReg

PacingIncVPulseAmplReg

ΔPacingGbl

SChamberPacingIncPulseAmplReg[vP/o?]

*PacingVI*_{nc}*VPulseAmplReg* —

Δ *PacingGbl*
 Ξ *PacingPkgSt*
 Ξ *AVDelayGbl*
 Ξ *RateSmoothingGbl*
PacingIncVPulseAmplReg

PacingDecVPulseAmplReg —

Δ *PacingGbl*
SChamberPacingDecPulseAmplReg[*vP/o?*]

PacingVDecVPulseAmplReg —

Δ *PacingGbl*
 Ξ *PacingPkgSt*
 Ξ *AVDelayGbl*
 Ξ *RateSmoothingGbl*
PacingDecVPulseAmplReg

PacingIncAPulseAmplUReg —

Δ *PacingGbl*
SChamberPacingIncPulseAmplUReg[*aP/o?*]

*PacingVI*_{nc}*APulseAmplUReg* —

Δ *PacingGbl*
 Ξ *PacingPkgSt*
 Ξ *AVDelayGbl*
 Ξ *RateSmoothingGbl*
PacingIncAPulseAmplUReg

PacingDecAPulseAmplUReg —

Δ *PacingGbl*
SChamberPacingDecPulseAmplUReg[*aP/o?*]

PacingVDecAPulseAmplUReg —

Δ *PacingGbl*
 Ξ *PacingPkgSt*
 Ξ *AVDelayGbl*
 Ξ *RateSmoothingGbl*
PacingDecAPulseAmplUReg

PacingIncVPulseAmplUReg —

$\Delta PacingGbl$

$SChamberPacingIncPulseAmplUReg[vP/o?]$

PacingVIncVPulseAmplUReg —

$\Delta PacingGbl$

$\Xi PacingPkgSt$

$\Xi AVDelayGbl$

$\Xi RateSmoothingGbl$

PacingIncVPulseAmplUReg

PacingDecVPulseAmplUReg —

$\Delta PacingGbl$

$SChamberPacingDecPulseAmplUReg[vP/o?]$

PacingVDecVPulseAmplUReg —

$\Delta PacingGbl$

$\Xi PacingPkgSt$

$\Xi AVDelayGbl$

$\Xi RateSmoothingGbl$

PacingDecVPulseAmplUReg

PacingIncAPulseWidth —

$\Delta PacingGbl$

$SChamberPacingIncPulseWidth[aP/o?]$

PacingVIncAPulseWidth —

$\Delta PacingGbl$

$\Xi PacingPkgSt$

$\Xi AVDelayGbl$

$\Xi RateSmoothingGbl$

PacingIncAPulseWidth

PacingDecAPulseWidth —

$\Delta PacingGbl$

$SChamberPacingDecPulseWidth[aP/o?]$

PacingVDecAPulseWidth _____

$\Delta PacingGbl$
 $\Xi PacingPkgSt$
 $\Xi AVDelayGbl$
 $\Xi RateSmoothingGbl$
PacingDecAPulseWidth

PacingIncVPulseWidth _____

$\Delta PacingGbl$
SChamberPacingIncPulseWidth [$vP/o?$]

PacingVIIncVPulseWidth _____

$\Delta PacingGbl$
 $\Xi PacingPkgSt$
 $\Xi AVDelayGbl$
 $\Xi RateSmoothingGbl$
PacingIncVPulseWidth

PacingDecVPulseWidth _____

$\Delta PacingGbl$
SChamberPacingDecPulseWidth [$vP/o?$]

PacingVDecVPulseWidth _____

$\Delta PacingGbl$
 $\Xi PacingPkgSt$
 $\Xi AVDelayGbl$
 $\Xi RateSmoothingGbl$
PacingDecVPulseWidth

PacingIsHRLInc5 _____

PacingGblSt
 $hrl < 50 \vee hrl \geq 90 \wedge hrl < 175$

PacingIncHRL5 _____

$\Delta PacingGbl$
PacingIsHRLInc5
 $hrl' = hrl + 5$

PacingIncHRL1

$\Delta PacingGbl$

$hrl \geq 50 \wedge hrl < 90$

$hrl' = hrl + 1$

PacingIncHRLSame

$\Delta PacingGbl$

$hrl = 175$

$hrl' = hrl$

PacingIncHRL

$\Delta PacingGbl$

$PacingIncHRL5 \vee PacingIncHRL1 \vee PacingIncHRLSame$

$PacingIncHRLUFr == PacingGbl \setminus (hrl)$

$PacingIncHRLFr == \Delta PacingGbl \wedge \exists PacingIncHRLUFr$

PacingVIncHRL

$PacingIncHRLFr$

$PacingIncHRL$

PacingIsHRLDec5

$PacingGblSt$

$hrl > 90 \vee hrl \leq 50 \wedge hrl > 30$

PacingDecHRL5

$\Delta PacingGbl$

$PacingIsHRLDec5$

$hrl' = hrl - 5$

PacingDecHRL1

$\Delta PacingGbl$

$hrl \leq 90 \wedge hrl > 50$

$hrl' = hrl - 1$

PacingDecHRLSame

$\Delta PacingGbl$

$hrl = 30$

$hrl' = hrl$

PacingDecHRL

$\Delta PacingGbl$

$PacingDecHRL5 \vee PacingDecHRL1 \vee PacingDecHRLSame$

$PacingDecHRLUFr == PacingGbl \setminus (hrl)$

$PacingDecHRLFr == \Delta PacingGbl \wedge \exists PacingDecHRLUFr$

PacingVDecHRL

$PacingDecHRLFr$

$PacingDecHRL$

PacingUpdateCCI0

$\Delta PacingGbl$

$RateSmoothingCalcSmoothedCCI[targetCCI/targetCCI?, currentCCI/currentCCI?]$

$currentCCI' = newCCI!$

$PacingUpdateCCI == \exists newCCI! : MilliSec \bullet PacingUpdateCCI0$

PacingCurrentCCISame

$\Delta PacingGbl$

$currentCCI' = currentCCI$

PacingAPaceEnabled

$PacingGbl$

$cp = A \vee cp = D$

PacingDoesARateSmoothing

PacingGbl

PacingAPaceEnabled

$\neg (cp = D)$

PacingUpdateCCIIfDoesRateSmoothing

$\Delta PacingGbl$

$(PacingDoesARateSmoothing \wedge PacingUpdateCCI) \vee PacingCurrentCCISame$

PacingAPaceCommon

$\Delta PacingGbl$

now? : MilliSec

PacingAPaceEnabled

SChamberPacingSetLastEventTm[aP/o?]

PacingUpdateCCIIfDoesRateSmoothing

PacingVPaceEnabled

PacingGblSt

$cp = V \vee cp = D$

PacingVPaceCommon

$\Delta PacingGbl$

now? : MilliSec

PacingVPaceEnabled

SChamberPacingSetLastEventTm[vP/o?]

PacingUpdateCCI

PacingTriggerAPace

$\Delta PacingGbl$

now? : MilliSec

aPulseAmp! : MilliVolt

aPulseWidth! : MicroSec

PacingAPaceCommon

SChamberPacingGetPulseWidth[aP/o?, aPulseWidth!/pulseWidth!]

SChamberPacingGetPulseAmp![aP/o?, aPulseAmp!/pulseAmp!]

$PacingTriggerAPaceUFr == PacingGbl \setminus (currentCCI, sChamberPacing, stChamberPacing)$
 $PacingTriggerAPaceFr == \Delta PacingGbl \wedge \exists PacingTriggerAPaceUFr$

<i>PacingVTriggerAPace</i>	_____
<i>PacingTriggerAPaceFr</i>	
<i>PacingTriggerAPace</i>	

<i>PacingTriggerVPace</i>	_____
$\Delta PacingGbl$	
<i>now? : MilliSec</i>	
<i>vPulseAmpl! : MilliVolt</i>	
<i>vPulseWidth! : MicroSec</i>	
<i>PacingVPaceCommon</i>	
<i>SChamberPacingGetPulseWidth[vP/o?, vPulseWidth!/pulseWidth!]</i>	
<i>SChamberPacingGetPulseAmpl[vP/o?, vPulseAmpl!/pulseAmpl!]</i>	

$PacingTriggerVPaceUFr == PacingGbl \setminus (currentCCI, sChamberPacing, stChamberPacing)$
 $PacingTriggerVPaceFr == \Delta PacingGbl \wedge \exists PacingTriggerVPaceUFr$

<i>PacingVTriggerVPace</i>	_____
<i>PacingTriggerVPaceFr</i>	
<i>PacingTriggerVPace</i>	

<i>PacingInhibitAPace</i>	_____
$\Delta PacingGbl$	
<i>now? : MilliSec</i>	
<i>PacingAPaceCommon</i>	

$PacingInhibitAPaceUFr == PacingGbl \setminus (currentCCI, sChamberPacing, stChamberPacing)$
 $PacingInhibitAPaceFr == \Delta PacingGbl \wedge \exists PacingInhibitAPaceUFr$

$PacingVIInhibitAPace$ _____
 $PacingInhibitAPaceFr$
 $PacingInhibitAPace$

$PacingInhibitVPace$ _____
 $\Delta PacingGbl$
 $now? : MilliSec$
 $PacingVPaceCommon$

$PacingInhibitVPaceUFr == PacingGbl \setminus (currentCCI, sChamberPacing, stChamberPacing)$
 $PacingInhibitVPaceFr == \Delta PacingGbl \wedge \exists PacingInhibitVPaceUFr$

$PacingVIInhibitVPace$ _____
 $PacingInhibitVPaceFr$
 $PacingInhibitVPace$

$PacingTimeToAPace0$ _____
 $PacingGblSt$
 $now? : MilliSec$
 $SChamberPacingGetLastEventTm[aP/o?]$
 $lastEvTm! = now? - currentCCI$

$PacingTimeToAPace == \exists lastEvTm! : \mathbb{Z} \bullet$
 $PacingTimeToAPace0$

$PacingTimedAPace$ _____
 $\Delta PacingGbl$
 $PacingTimeToAPace$
 $PacingTriggerAPace$

$PacingTimedAPaceUFr == PacingGbl \setminus (currentCCI, sChamberPacing, stChamberPacing)$
 $PacingTimedAPaceFr == \Delta PacingGbl \wedge \exists PacingTimedAPaceUFr$

PacingVTimedAPace

PacingTimedAPaceFr
PacingTimedAPace

PacingTimeToVPace0

PacingGblSt
now? : MilliSec
SChamberPacingGetLastEventTm[vP/o?]

lastEvTm! = now? - currentCCI

PacingTimeToVPace == $\exists \text{lastEvTm!} : \mathbb{Z}$ •
PacingTimeToVPace0

PacingTimeToAVDelay

PacingGblSt
now? : MilliSec
AVDelayGetAVDelay
SChamberPacingGetLastEventTm[aP/o?]

lastEvTm! = now? - avi!

lastEvTm! > now? - avi!

cci? = currentCCI

PacingTimedVPace0

Δ *PacingGbl*
PacingTimeToVPace
PacingTriggerVPace
avi! : MilliSec
cci? : MilliSec
lastEvTm! : MilliSec

PacingTimeToAPace \vee *PacingTimeToAVDelay*

PacingTimedVPace == $\exists \text{avi!} : \mathbb{Z}; \text{cci} : \mathbb{Z}; \text{lastEvTm!} : \mathbb{Z}$ •
PacingTimedVPace0

PacingTimedVPaceUFr == PacingGbl \ (currentCCI, sChamberPacing, stChamberPacing)

PacingTimedVPaceFr == \Delta PacingGbl \ \Xi PacingTimedVPaceUFr

PacingVTimedVPace

PacingTimedVPaceFr

PacingTimedVPace

PacingTrackASense

$\Delta PacingGbl$

PacingTimeToVPace

PacingTriggerVPace

now? : MilliSec

SChamberPacingSetLastEventTm[aP/o?]

PacingTrackASenseUFr == PacingGbl \ (currentCCI, sChamberPacing, stChamberPacing)

PacingTrackASenseFr == \Delta PacingGbl \wedge \exists PacingTrackASenseUFr

PacingVTrackASense

PacingTrackASenseFr

PacingTrackASense

PacingExecuteEventATRFallback

$\Delta PacingGbl$

ev? : PMEvent

PacingSetBradycardiaMode[V / cp?, R / rm?]

ev? = EvATRFallback

targetCCI' = targetCCI

PacingExecuteEventTerminateATR

$\Delta PacingGbl$

ev? : PMEvent

ev? = EvTerminateATR

cp' = pcp

pcp' = cp

rm' = prm

prm' = rm

targetCCI' = targetCCI

PacingExecuteEventFallbackTmElapsed _____

$\Delta PacingGbl$
 $ev? : PMEvent$
 $RateModulationVSetTargetCCIToLRL$

$ev? = EvFallbackTmElapsed$

$cp' = cp$

$pcp' = pcp$

$rm' = rm$

$prm' = prm$

PacingExecuteEventNone _____

$\Delta PacingGbl$
 $ev? : PMEvent$

$ev? = EvNone$

$cp' = cp$

$pcp' = pcp$

$rm' = rm$

$prm' = prm$

$targetCCI' = targetCCI$

PacingExecuteEvent == PacingExecuteEventATRFallback

$\vee PacingExecuteEventTerminateATR$
 $\vee PacingExecuteEventFallbackTmElapsed$
 $\vee PacingExecuteEventNone$

PacingExecuteEventUFr == PacingGbl \ (cp, pcp, rm, prm, targetCCI)

PacingExecuteEventFr == $\Delta PacingGbl \wedge \exists PacingExecuteEventUFr$

PacingVEexecuteEvent _____

PacingExecuteEventFr
PacingExecuteEvent

C.8 Package Sensing

section *Sensing parents ZOO_Toolkit, Model_Preamble, CommonPM*

Sensitivity == { $o : MicroVolt \mid o \geq 250 \wedge o \leq 10000$ }

Blanking == { $o : MilliSec \mid o \geq 30 \wedge o \leq 60$ }

SensingGblSt _____
 aSensitivity : *Sensitivity*
 vSensitivity : *Sensitivity*
 cs : *ChambersSensed*
 aBlanking : *Blanking*
 vBlanking : *Blanking*
 pcs : *ChambersSensed*

SensingGbl _____
 SensingGblSt

SensingInit _____
 SensingGbl
 aSensitivity = 750
 vSensitivity = 2500
 cs = D
 aBlanking = 40
 vBlanking = 40
 pcs = D

SensingSetBradycardiaMode _____
 Δ *SensingGbl*
 cs? : *ChambersSensed*
 cs' = cs?
 pcs' = cs

SensingSetBradycardiaModeUFr == *SensingGbl* $\setminus (cs, pcs)$

SensingSetBradycardiaModeFr == Δ *SensingGbl* $\wedge \exists$ *SensingSetBradycardiaModeUFr*

$\text{SensingVSetBradycardiaMode}$
 $\text{SensingSetBradycardiaModeFr}$
 $\text{SensingSetBradycardiaMode}$

$\text{SensingIsBatteryDegraded}$
 $\Delta \text{SensingGblSt}$

$bs? : \text{BatteryStatus}$

$bs? = ERT \vee bs? = ERP$

$\text{SensingUpdateBatteryChangeMode}$
 $\Delta \text{SensingGblSt}$

$\text{SensingIsBatteryDegraded}$

$\text{SensingSetBradycardiaMode}[V / cs?]$

$\text{SensingUpdateBatterySame}$
 $\Delta \text{SensingGblSt}$

$cs' = cs$

$pcs' = pcs$

$\text{SensingUpdateBatteryStatus}$
 $\Delta \text{SensingGbl}$

$bs? : \text{BatteryStatus}$

$\text{SensingUpdateBatteryChangeMode} \vee \text{SensingUpdateBatterySame}$

$\text{SensingUpdateBatteryStatusUfr} == \text{SensingGbl} \setminus (cs, pcs)$

$\text{SensingUpdateBatteryStatusFr} == \Delta \text{SensingGbl} \wedge \exists \text{SensingUpdateBatteryStatusUfr}$

$\text{SensingVUpdateBatteryStatus}$
 $\text{SensingUpdateBatteryStatusFr}$

$\text{SensingUpdateBatteryStatus}$

SensingExecuteEventATRFallback _____

$\Delta \text{SensingGbl}$
 $\text{ev?} : \text{PMEvent}$
 $\text{SensingSetBradycardiaMode}[D / cs?]$
 $\text{ev?} = \text{EvATRFallback}$

SensingExecuteEventTerminateATR _____

$\Delta \text{SensingGbl}$
 $\text{ev?} : \text{PMEvent}$
 $\text{ev?} = \text{EvTerminateATR}$
 $\text{cs}' = \text{pcs}$
 $\text{pcs}' = \text{cs}$

SensingEvsDoingNothing _____

SensingGbl
 $\text{ev?} : \text{PMEvent}$
 $\text{ev?} = \text{EvFallbackTmElapsed} \vee \text{ev?} = \text{EvNone}$

SensingExecuteEventDoingNothing _____

$\Delta \text{SensingGbl}$
 $\text{SensingEvsDoingNothing}$
 $\text{ev?} : \text{PMEvent}$
 $\text{cs}' = \text{pcs}$
 $\text{pcs}' = \text{cs}$

$\text{SensingExecuteEvent} == \text{SensingExecuteEventATRFallback}$

$\vee \text{SensingExecuteEventTerminateATR}$
 $\vee \text{SensingExecuteEventDoingNothing}$

$\text{SensingExecuteEventUFr} == \text{SensingGbl} \setminus (\text{cs}, \text{pcs})$

$\text{SensingExecuteEventFr} == \Delta \text{SensingGbl} \wedge \exists \text{SensingExecuteEventUFr}$

SensingVExecuteEvent _____

$\text{SensingExecuteEventFr}$
 $\text{SensingExecuteEvent}$

SensingASenseOnly

SensingGblSt
sAAmpl? : MicroVolt
now? : MilliSec
lastAEvTm? : MilliSec

sAAmpl? > aSensitivity
aBlanking < now? - lastAEvTm?
cs = A ∨ cs = D

SensingVASenseOnly

SensingASenseOnly

SensingTriggerAPace

SensingGbl
SensingASenseOnly

SensingVTriggerAPace

SensingTriggerAPace

SensingInhibitAPace

SensingGbl
SensingASenseOnly

SensingVIhibitAPace

SensingInhibitAPace

SensingASenseRefractory

SensingGbl
SensingASenseOnly

SensingVASenseRefractory

SensingASenseRefractory

SensingVSenseOnly

SensingGblSt

sVAmpl? : MicroVolt

now? : MilliSec

lastVEvTm? : MilliSec

sVAmpl? > vSensitivity

vBlanking < now? - lastVEvTm?

cs = V ∨ cs = D

SensingVVSenseOnly

SensingGbl

SensingVSenseOnly

SensingVSenseRefractory

SensingGbl

SensingVSenseOnly

SensingVVSenseRefractory

SensingVSenseRefractory

SensingTriggerVPace

SensingGbl

SensingVSenseOnly

SensingVTriggerVPace

SensingTriggerVPace

SensingInhibitVPace

SensingGbl

SensingVSenseOnly

SensingVInhibitVPace _____
 ┌── *SensingInhibitVPace*

SensingTrackASense _____
 ┌── *SensingGbl*
 ┌── *SensingASenseOnly*

SensingVTrackASense _____
 ┌── *SensingTrackASense*

SensingIncASensitivityLow _____
 ┌── *ΔSensingGbl*
 ┌── *aSensitivity < 100*
 ┌── *aSensitivity' = aSensitivity + 250*

SensingIncASensitivityMedium _____
 ┌── *ΔSensingGbl*
 ┌── *aSensitivity ≥ 1000 ∧ aSensitivity < 10000*
 ┌── *aSensitivity' = aSensitivity + 500*

SensingIncASensitivitySame _____
 ┌── *ΔSensingGbl*
 ┌── *aSensitivity = 1000*
 ┌── *aSensitivity' = aSensitivity*

SensingIncASensitivity == SensingIncASensitivityLow
 \vee *SensingIncASensitivityMedium* \vee *SensingIncASensitivitySame*

SensingIncASensitivityUFr == SensingGbl \ (aSensitivity)
SensingIncASensitivityFr == ΔSensingGbl \ \exists SensingIncASensitivityUFr

SensingVI \in *ASensitivity*

SensingIncA \in *ASensitivityFr*

SensingIncA \in *ASensitivity*

SensingDecA \in *ASensitivityLow*

Δ *SensingGbl*

$aSensitivity \leq 1000 \wedge aSensitivity > 250$

$aSensitivity' = aSensitivity - 250$

SensingDecA \in *ASensitivityMedium*

Δ *SensingGbl*

$aSensitivity > 1000$

$aSensitivity' = aSensitivity - 500$

SensingDecA \in *ASensitivitySame*

Δ *SensingGbl*

$aSensitivity = 250$

$aSensitivity' = aSensitivity$

SensingDecA \in *ASensitivity == SensingDecA* \in *ASensitivityLow*

\vee *SensingDecA* \in *ASensitivityMedium*

\vee *SensingDecA* \in *ASensitivitySame*

SensingDecA \in *ASensitivityU* \neq *SensingGbl* \setminus ($aSensitivity$)

SensingDecA \in *ASensitivityFr* \equiv Δ *SensingGbl* \wedge \exists *SensingDecA* \in *ASensitivityU* \neq

SensingVI \in *DecASensitivity*

SensingDecA \in *ASensitivityFr*

SensingDecA \in *ASensitivity*

SensingIncV \in *VSensitivityLow*

Δ *SensingGbl*

$vSensitivity < 100$

$vSensitivity' = aSensitivity + 250$

SensingIncVSensitivityMedium _____

$\Delta S \text{ensing} Gbl$

$vSensitivity \geq 1000 \wedge vSensitivity < 10000$

$vSensitivity' = vSensitivity + 500$

SensingIncVSensitivitySame _____

$\Delta S \text{ensing} Gbl$

$vSensitivity = 1000$

$vSensitivity' = vSensitivity$

SensingIncVSensitivity == SensingIncVSensitivityLow

$\vee S \text{sensing} Inc VSensitivity Medium$

$\vee S \text{sensing} Inc VSensitivity Same$

SensingIncVSensitivityUFr == SensingGbl \setminus (vSensitivity)

SensingIncVSensitivityFr == \Delta S \text{ensing} Gbl \wedge \exists S \text{sensing} Inc VSensitivity UFr

SensingVIcVSensitivity _____

$S \text{sensing} Inc VSensitivity Fr$

$S \text{sensing} Inc VSensitivity$

SensingDecVSensitivityLow _____

$\Delta S \text{ensing} Gbl$

$vSensitivity \leq 1000 \wedge vSensitivity > 250$

$vSensitivity' = vSensitivity - 250$

SensingDecVSensitivityMedium _____

$\Delta S \text{ensing} Gbl$

$vSensitivity > 1000$

$vSensitivity' = vSensitivity - 500$

SensingDecVSensitivitySame _____

$\Delta S_{\text{SensingGbl}}$

$v_{\text{Sensitivity}} = 250$

$v_{\text{Sensitivity}'} = v_{\text{Sensitivity}}$

SensingDecVSensitivity == SensingDecVSensitivityLow

$\vee S_{\text{SensingDecVSensitivityMedium}}$

$\vee S_{\text{SensingDecVSensitivitySame}}$

SensingDecVSensitivityUFr == S_{\text{SensingGbl}} \setminus (v_{\text{Sensitivity}})

SensingDecVSensitivityFr == \Delta S_{\text{SensingGbl}} \wedge \exists S_{\text{SensitivityUFr}}

SensingVDecVSensitivity _____

SensingDecVSensitivityFr

SensingDecVSensitivity

SensingIncABlankingInc _____

$\Delta S_{\text{SensingGbl}}$

$a_{\text{Blanking}} < 60$

$a_{\text{Blanking}'} = a_{\text{Blanking}} + 10$

SensingIncABlankingSame _____

$\Delta S_{\text{SensingGbl}}$

$a_{\text{Blanking}} = 60$

$a_{\text{Blanking}'} = a_{\text{Blanking}}$

SensingIncABlanking == SensingIncABlankingInc

$\vee S_{\text{SensingIncABlankingSame}}$

SensingIncABlankingUFr == S_{\text{SensingGbl}} \setminus (a_{\text{Blanking}})

SensingIncABlankingFr == \Delta S_{\text{SensingGbl}} \wedge \exists S_{\text{SensitivityUFr}}

SensingVIncABlanking _____

SensingIncABlankingFr

SensingIncABlanking

SensingDecABlankingDec _____

$\Delta \text{SensingGbl}$

$a\text{Blanking} > 30$

$a\text{Blanking}' = a\text{Blanking} - 10$

SensingDecABlankingSame _____

$\Delta \text{SensingGbl}$

$a\text{Blanking} = 30$

$a\text{Blanking}' = a\text{Blanking}$

SensingDecABlanking == SensingDecABlankingDec

$\vee \text{SensingDecABlankingSame}$

SensingDecABlankingUFr == SensingGbl \setminus (a\text{Blanking})

SensingDecABlankingFr == $\Delta \text{SensingGbl} \wedge \exists \text{SensingDecABlankingUFr}$

SensingVDecABlanking _____

SensingDecABlankingFr

SensingDecABlanking

SensingIncVBlankingInc _____

$\Delta \text{SensingGbl}$

$v\text{Blanking} < 60$

$v\text{Blanking}' = v\text{Blanking} + 10$

SensingIncVBlankingSame _____

$\Delta S \text{ensing} G \text{bl}$

$v \text{Blanking} = 60$

$v \text{Blanking}' = v \text{Blanking}$

SensingIncVBlanking == SensingIncVBlankingInc

$\vee S \text{ensing} I \text{nc} V \text{Blanking} S \text{ame}$

SensingIncVBlankingUFr == SensingGbl \setminus (v \text{Blanking})

SensingIncVBlankingFr == \Delta S \text{ensing} G \text{bl} \wedge \exists S \text{ensing} Inc V \text{Blanking} U \text{Fr}

SensingVIncVBlanking _____

SensingIncVBlankingFr

SensingIncVBlanking

SensingDecVBlankingDec _____

$\Delta S \text{ensing} G \text{bl}$

$v \text{Blanking} > 30$

$v \text{Blanking}' = v \text{Blanking} - 10$

SensingDecVBlankingSame _____

$\Delta S \text{ensing} G \text{bl}$

$v \text{Blanking} = 30$

$v \text{Blanking}' = v \text{Blanking}$

SensingDecVBlanking == SensingDecVBlankingDec

$\vee S \text{ensing} Dec V \text{Blanking} S \text{ame}$

SensingDecVBlankingUFr == SensingGbl \setminus (v \text{Blanking})

SensingDecVBlankingFr == \Delta S \text{ensing} G \text{bl} \wedge \exists S \text{ensing} Dec V \text{Blanking} U \text{Fr}

SensingVDecVBlanking _____

SensingDecVBlankingFr

SensingDecVBlanking

C.9 Package ResponseCommon

```
section ResponseCommon parents ZOO_Toolkit, Model_Preamble, CommonPM
```

```
RefractoryPeriod == {o : MilliSec | o ≥ 150 ∧ o ≤ 500}
```

C.10 Package AtrialResponse

```
section AtrialResponse parents ZOO_Toolkit, Model_Preamble, CommonPM, ResponseCommon
```

```
| detectionPeriod : MilliSec
```

```
UpperRateLimit == {o : PPM | o ≥ 50 ∧ o ≤ 175}
```

```
ATRDuration == {o : N | o ≥ 10 ∧ o ≤ 2000}
```

```
ATRFallbackTime == {o : MilliSec | o ≥ 60000 ∧ o ≤ 300000}
```

```
ATRDetectionCount == {o : N | o ≥ 1 ∧ o ≤ 8}
```

```
ATRStatus ::= NoATR | Duration | Fallback | FallbackTmElapsed
```

```
AtrialResponseGblSt _____  
url : UpperRateLimit  
arp : RefractoryPeriod  
atrMode : Switch  
atrDuration : ATRDuration  
atrFbTm : ATRFallbackTime  
atrEntryC : ATRDetectionCount  
atrExitC : ATRDetectionCount  
atrDetC : N  
atrStartDetTm : MilliSec  
atrStatus : ATRStatus  
atrDurC : N  
atrStartFbTm : MilliSec
```

AtrialResponseGbl _____
| *AtrialResponseGblSt*

— *AtrialResponseInit* _____
| *AtrialResponseGbl*

| *url* = 120
| *arp* = 250
| *atrMode* = *Off*
| *atrDuration* = 20 \wedge *atrFbTm* = 60000
| *atrEntryC* = 8 \wedge *atrExitC* = 8
| *atrStartDetTm* = 0 \wedge *atrStartFbTm* = 0
| *atrStatus* = *NoATR* \wedge (*atrDetC* = 0 \wedge *atrDurC* = 0)

— *AtrialResponseInATRFallBackMode* _____
| *AtrialResponseGbl*

| *atrStatus* = *Fallback* \vee *atrStatus* = *FallbackTmElapsed*

— *AtrialResponseVInATRFallBackMode* _____
| *AtrialResponseInATRFallBackMode*

— *AtrialResponseDoNoATR* _____
| Δ *AtrialResponseGbl*
| *ev!* : *PMEEvent*
| *atrMode* = *Off*
| *ev!* = *EvNone*
| *atrStartDetTm'* = *atrStartDetTm*
| *atrDetC'* = *atrDetC*
| *atrStatus'* = *atrStatus*
| *atrDurC'* = *atrDurC*

AtrialResponseSetPMEv _____

<i>AtrialResponseGbl</i>
<i>ev!</i> : PMEvent
<i>evDetect?</i> : PMEvent
<i>evDurationFallback?</i> : PMEvent

$(evDetect? \neq EvNone \wedge ev! = evDetect?) \vee ev! = evDurationFallback?$

AtrialResponseGetURLI _____

<i>AtrialResponseGbl</i>
<i>urlI!</i> : MilliSec

$urlI! = 60000 \text{ div } url$

AtrialResponseATRDetectEntryReset _____

$\Delta AtrialResponseGbl$
<i>now?</i> : MilliSec
<i>laet?</i> : MilliSec
<i>evDetect!</i> : PMEvent
<i>AtrialResponseGetURLI</i>

$detectionPeriod < now? - atrStartDetTm$

$urlI! < now? - laet?$

$atrStartDetTm' = now?$

$atrDetC' = 1$

$evDetect! = EvNone$

$atrDurC' = atrDurC$

$atrStatus' = atrStatus$

AtrialResponseATRDetectEntryInc _____

$\Delta AtrialResponseGbl$
 $now? : MilliSec$
 $laet? : MilliSec$
 $evDetect! : PMEvent$
 $AtrialResponseGetURLI$

$urlI! < now? - laet?$
 $detectionPeriod \geq now? - atrStartDetTm$
 $atrEntryC > atrDetC + 1$
 $atrDetC' = atrDetC + 1$
 $evDetect! = EvNone$
 $atrStartDetTm' = atrStartDetTm$
 $atrDurC' = atrDurC$
 $atrStatus' = atrStatus$

AtrialResponseATRDetectEntryToDuration _____

$\Delta AtrialResponseGbl$
 $now? : MilliSec$
 $laet? : MilliSec$
 $evDetect! : PMEvent$
 $AtrialResponseGetURLI$

$urlI! < now? - laet?$
 $detectionPeriod \geq now? - atrStartDetTm$
 $atrEntryC = atrDetC + 1$
 $atrStatus' = Duration$
 $atrDurC' = 0$
 $evDetect! = EvATRDuration$
 $atrStartDetTm' = atrStartDetTm$
 $atrDetC' = atrDetC$

AtrialResponseATRDetectEntry0 == AtrialResponseATRDetectEntryReset

$\vee AtrialResponseATRDetectEntryInc$
 $\vee AtrialResponseATRDetectEntryToDuration$

AtrialResponseATRDetectEntry _____

$\Delta AtrialResponseGbl$

AtrialResponseATRDetectEntry0

$atrStatus = NoATR$

AtrialResponseATRDetectExitReset _____

$\Delta AtrialResponseGbl$

$now? : MilliSec$

$laet? : MilliSec$

$evDetect! : PMEvent$

AtrialResponseGetURLI

$detectionPeriod < now? - atrStartDetTm$

$urlI! \geq now? - laet?$

$atrStartDetTm' = now?$

$atrDetC' = 1$

$evDetect! = EvNone$

$atrDurC' = atrDurC$

$atrStatus' = atrStatus$

AtrialResponseATRDetectExitInc _____

$\Delta AtrialResponseGbl$

$now? : MilliSec$

$laet? : MilliSec$

$evDetect! : PMEvent$

AtrialResponseGetURLI

$urlI! \geq now? - laet?$

$detectionPeriod \geq now? - atrStartDetTm$

$atrExitC > atrDetC + 1$

$atrDetC' = atrDetC + 1$

$evDetect! = EvNone$

$atrStartDetTm' = atrStartDetTm$

$atrDurC' = atrDurC$

$atrStatus' = atrStatus$

AtrialResponseATRDetectExitSetPMEEvent

AtrialResponseGbl
evDetect! : PMEvent

atrStatus = Fallback \Rightarrow *evDetect!* = EvTerminateATR

atrStatus \neq Fallback \Rightarrow *evDetect!* = EvNone

AtrialResponseATRDetectEntryToNoATR

Δ *AtrialResponseGbl*

now? : MilliSec

laet? : MilliSec

AtrialResponseGetURLI

AtrialResponseATRDetectExitSetPMEEvent

urlI! \geq *now?* – *laet?*

detectionPeriod \geq *now?* – *atrStartDetTm*

atrExitC = *atrDetC* + 1

atrStatus' = NoATR

atrDurC' = *atrDurC*

atrStartDetTm' = *atrStartDetTm*

atrDetC' = *atrDetC*

AtrialResponseATRDetectExit0 == *AtrialResponseATRDetectExitReset*

\vee *AtrialResponseATRDetectExitInc*

\vee *AtrialResponseATRDetectEntryToNoATR*

AtrialResponseATRDetectExit

Δ *AtrialResponseGbl*

AtrialResponseATRDetectExit0

atrStatus \neq NoATR

AtrialResponseATRDetect == *AtrialResponseATRDetectEntry*

\vee *AtrialResponseATRDetectExit*

AtrialResponseATRDoDurationInc _____

$\Delta AtrialResponseGbl$
 $evDurationFallback! : PMEvent$
 $atrDuration > atrDurC + 1$
 $atrDurC' = atrDurC + 1$
 $evDurationFallback! = EvNone$
 $atrStatus' = atrStatus$
 $atrStartFbTm' = atrStartFbTm$

AtrialResponseATRDoDurationToFallback _____

$\Delta AtrialResponseGbl$
 $now? : MilliSec$
 $evDurationFallback! : PMEvent$
 $atrDuration = atrDurC + 1$
 $atrStatus' = Fallback$
 $atrStartFbTm' = now?$
 $evDurationFallback! = EvATRFallback$
 $atrDurC' = atrDurC$

AtrialResponseATRDoDuration0 == AtrialResponseATRDoDurationInc

$\vee AtrialResponseATRDoDurationToFallback$

AtrialResponseATRDoDuration _____

$\Delta AtrialResponseGbl$
 $AtrialResponseATRDoDuration0$
 $atrStatus = Duration$

AtrialResponseATRFallbackTmNotElapsed _____

$\Delta AtrialResponseGbl$
 $now? : MilliSec$
 $evDurationFallback! : PMEvent$
 $atrFbTm < now? - atrStartFbTm$
 $evDurationFallback! = EvNone$
 $atrStatus' = atrStatus$

AtrialResponseATRFallbackTmElapsed _____

$\Delta AtrialResponseGbl$
 $now? : MilliSec$
 $evDurationFallback! : PMEvent$

$atrFbTm \geq now? - atrStartFbTm$
 $evDurationFallback! = EvFallbackTmElapsed$
 $atrStatus' = FallbackTmElapsed$

AtrialResponseATRFallback0 == AtrialResponseATRFallbackTmNotElapsed

$\vee AtrialResponseATRFallbackTmElapsed$

AtrialResponseATRFallback _____

$\Delta AtrialResponseGbl$
AtrialResponseATRFallback0

$atrStatus = Fallback$
 $atrStartFbTm' = atrStartFbTm$
 $atrDurC' = atrDurC$

AtrialResponseATRDoNothing _____

$\Delta AtrialResponseGbl$
 $evDurationFallback! : PMEvent$

$atrStatus \neq Duration \wedge atrStatus \neq Fallback$
 $evDurationFallback! = EvNone$
 $atrStatus' = atrStatus$
 $atrStartFbTm' = atrStartFbTm$
 $atrDurC' = atrDurC$

AtrialResponseDoATRDisj0 == AtrialResponseATRDoDuration

$\vee AtrialResponseATRFallback$
 $\vee AtrialResponseATRDoNothing$

AtrialResponseDoATRDisj _____

$\Delta AtrialResponseGbl$
AtrialResponseDoATRDisj0

$atrDetC' = atrDetC$
 $atrStartDetTm' = atrStartDetTm$

$AtrialResponseDoATRSemi0 == AtrialResponseATRDetect$
 $\quad \quad \quad \dagger AtrialResponseDoATRDisj$

$AtrialResponseDoATR0$ _____

$\Delta AtrialResponseGbl$
$now? : MilliSec$
$ev! : PMEvent$
$laet? : MilliSec$
$AtrialResponseDoATRSemi0$
$AtrialResponseSetPMEv[evDetect!/evDetect?, evDurationFallback!/evDurationFallback?]$
$atrMode = On$

$AtrialResponseDoATR == AtrialResponseDoATR0 \setminus (evDetect!, evDurationFallback!)$

$AtrialResponsePerformATR == AtrialResponseDoATR \vee AtrialResponseDoNoATR$

$AtrialResponsePerformATRUFr == AtrialResponseGbl \setminus (atrDurC, atrStartFbTm, atrStatus, atrDetC, atrStartFbTm)$
 $AtrialResponsePerformATRFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponsePerformATRUFr$

$AtrialResponseVPerformATR$ _____

$AtrialResponsePerformATRFr$
$AtrialResponsePerformATR$

$AtrialResponseIsInRefractory$ _____

$AtrialResponseGbl$
$now? : MilliSec$
$lastAEvTm? : MilliSec$
$lastVEvTm? : MilliSec$
$pvarp? : MilliSec$
$arp \leq now? - lastAEvTm? \vee pvarp? \leq now? - lastVEvTm?$

AtrialResponseASenseOnly

$\Delta AtrialResponseGbl$
now? : MilliSec
lastAEvTm? : MilliSec
lastVEvTm? : MilliSec
pvarp? : MilliSec
AtrialResponsePerformATR

$\neg AtrialResponseIsInRefractory$

$AtrialResponseASenseOnlyUFr == AtrialResponseGbl \setminus (atrDurC, atrStartFbTm, atrStatus, atrDetC, atrStart)$
 $AtrialResponseASenseOnlyFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseASenseOnlyUFr$

AtrialResponseVASenseOnly

AtrialResponseASenseOnlyFr
AtrialResponseASenseOnly

AtrialResponseASenseRefractory

$\Delta AtrialResponseGbl$
now? : MilliSec
lastAEvTm? : MilliSec
lastVEvTm? : MilliSec
pvarp? : MilliSec
AtrialResponsePerformATR

AtrialResponseIsInRefractory

$AtrialResponseASenseRefractoryUFr == AtrialResponseGbl \setminus (atrDurC, atrStartFbTm, atrStatus, atrDetC, atrStart)$
 $AtrialResponseASenseRefractoryFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseASenseRefractoryUFr$

AtrialResponseVASenseRefractory

AtrialResponseASenseRefractoryFr
AtrialResponseASenseRefractory

AtrialResponseIncURLInc

$\Delta AtrialResponseGbl$

url < 175

url' = url + 5

AtrialResponseIncURLSame

$\Delta AtrialResponseGbl$

$url = 175$

$url' = url$

$AtrialResponseIncURLUFr == AtrialResponseGbl \setminus (url)$

$AtrialResponseIncURLFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseIncURLUFr$

$AtrialResponseIncURL == AtrialResponseIncURLInc$

$\vee AtrialResponseIncURLSame$

AtrialResponseVIncURL

$AtrialResponseIncURLFr$

$AtrialResponseIncURL$

AtrialResponseDecURLDec

$\Delta AtrialResponseGbl$

$url > 50$

$url' = url - 5$

AtrialResponseDecURLSame

$\Delta AtrialResponseGbl$

$url = 50$

$url' = url$

$AtrialResponseDecURL == AtrialResponseDecURLDec$

$\vee AtrialResponseDecURLSame$

$AtrialResponseDecURLUFr == AtrialResponseGbl \setminus (url)$

$AtrialResponseDecURLFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseDecURLUFr$

$\neg AtrialResponseVDecURL$

$\Delta AtrialResponseDecURLFr$

$AtrialResponseDecURL$

$AtrialResponseIncARPInc$

$\Delta AtrialResponseGbl$

$arp < 500$

$arp' = arp + 10$

$AtrialResponseIncARPSame$

$\Delta AtrialResponseGbl$

$arp = 500$

$arp' = arp$

$AtrialResponseIncCARPUFr == AtrialResponseGbl \setminus (arp)$

$AtrialResponseIncARPFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseIncCARPUFr$

$AtrialResponseIncARP == AtrialResponseIncARPInc$

$\vee AtrialResponseIncARPSame$

$\neg AtrialResponseVIncARP$

$\Delta AtrialResponseIncARPFr$

$AtrialResponseIncARP$

$\neg AtrialResponseDecARPDec$

$\Delta AtrialResponseGbl$

$arp > 150$

$arp' = arp - 10$

$\neg AtrialResponseDecARPSame$

$\Delta AtrialResponseGbl$

$arp = 150$

$arp' = arp$

$AtrialResponseDecARP == AtrialResponseDecARPDec$
 $\vee AtrialResponseDecARPSame$

$AtrialResponseDecARPUFr == AtrialResponseGbl \setminus (arp)$
 $AtrialResponseDecARPFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseDecARPUFr$

$\neg AtrialResponseVDecARP$ _____
AtrialResponseDecARP
AtrialResponseDecARP

$AtrialResponseATROn$ _____
 $\Delta AtrialResponseGbl$
atrMode = Off
atrMode' = On

$AtrialResponseATROnUFr == AtrialResponseGbl \setminus (atrMode)$
 $AtrialResponseATROnFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseATROnUFr$

$\neg AtrialResponseVATROn$ _____
AtrialResponseATROnFr
AtrialResponseATROn

$AtrialResponseATROff$ _____
 $\Delta AtrialResponseGbl$
atrMode = On
atrMode' = Off

$AtrialResponseATROffUFr == AtrialResponseGbl \setminus (atrMode)$
 $AtrialResponseATROffFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseATROnUFr$

AtrialResponseVATROff

$\Delta AtrialResponseGbl$

$AtrialResponseATROffFr$

$AtrialResponseATROff$

AtrialResponseIncATRDurLow

$\Delta AtrialResponseGbl$

$atrDuration < 20$

$atrDuration' = atrDuration + 10$

AtrialResponseIncATRDurMedium

$\Delta AtrialResponseGbl$

$atrDuration \geq 20 \wedge atrDuration \leq 80$

$atrDuration' = atrDuration + 20$

AtrialResponseIncATRDurHigh

$\Delta AtrialResponseGbl$

$atrDuration \geq 100 \wedge atrDuration < 2000$

$atrDuration' = atrDuration + 100$

AtrialResponseIncATRDurSame

$\Delta AtrialResponseGbl$

$atrDuration = 2000$

$atrDuration' = atrDuration$

$AtrialResponseIncATRDuration == AtrialResponseIncATRDurLow$

$\vee AtrialResponseIncATRDurMedium \vee AtrialResponseIncATRDurHigh \vee AtrialResponseIncATRDurSame$

$AtrialResponseIncATRDurationUFr == AtrialResponseGbl \setminus (atrDuration)$

$AtrialResponseIncATRDurationFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseIncATRDurationUFr$

AtrialResponseVIncATRDur

$\Delta AtrialResponseIncATRDur$

$AtrialResponseIncATRDur$

AtrialResponseDecATRDurLow

$\Delta AtrialResponseGbl$

$atrDuration \leq 20$

$atrDuration' = 10$

AtrialResponseDecATRDurMedium

$\Delta AtrialResponseGbl$

$atrDuration > 20 \wedge atrDuration \leq 80$

$atrDuration' = atrDuration - 20$

AtrialResponseDecATRDurHigh

$\Delta AtrialResponseGbl$

$atrDuration \geq 100$

$atrDuration' = atrDuration - 100$

AtrialResponseDecATRDur == AtrialResponseDecATRDurLow

$\vee AtrialResponseDecATRDurMedium$

$\vee AtrialResponseDecATRDurHigh$

AtrialResponseDecATRDurUFr == AtrialResponseGbl \ (atrDuration)

AtrialResponseDecATRDurFr == $\Delta AtrialResponseGbl \wedge \exists AtrialResponseDecATRDurUFr$

AtrialResponseVDecATRDur

$\Delta AtrialResponseDecATRDur$

$AtrialResponseDecATRDur$

AtrialResponseIncATRFbTmInc

$\Delta AtrialResponseGbl$

$atrFbTm < 5$

$atrFbTm' = atrFbTm + 1$

AtrialResponseIncATRFbTmSame _____

$\Delta AtrialResponseGbl$

$atrFbTm = 5$

$atrFbTm' = atrFbTm$

AtrialResponseIncATRFbTmUFr == AtrialResponseGbl \ (atrFbTm)

AtrialResponseIncATRFbTmFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseIncATRFbTmUFr

AtrialResponseIncATRFbTm == AtrialResponseIncATRFbTmInc

$\vee AtrialResponseIncATRFbTmSame$

AtrialResponseVIncATRFbTm _____

AtrialResponseIncATRFbTmFr

AtrialResponseIncATRFbTm

AtrialResponseDecATRFbTmDec _____

$\Delta AtrialResponseGbl$

$atrFbTm > 1$

$atrFbTm' = atrFbTm - 1$

AtrialResponseDecATRFbTmSame _____

$\Delta AtrialResponseGbl$

$atrFbTm = 1$

$atrFbTm' = atrFbTm$

AtrialResponseDecATRFbTmUFr == AtrialResponseGbl \ (atrFbTm)

AtrialResponseDecATRFbTmFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseDecATRFbTmUFr

AtrialResponseDecATRFbTm == AtrialResponseDecATRFbTmDec

$\vee AtrialResponseDecATRFbTmSame$

$AtrialResponseVDecATRFbTm$ _____
 $\Delta AtrialResponseDecATRFbTmFr$
 $AtrialResponseDecATRFbTm$

$AtrialResponseIncATREntryCInc$ _____
 $\Delta AtrialResponseGbl$
 $atrEntryC < 8$
 $atrEntryC' = atrEntryC + 1$

$AtrialResponseIncATREntryCSame$ _____
 $\Delta AtrialResponseGbl$
 $atrEntryC = 8$
 $atrEntryC' = atrEntryC$

$AtrialResponseIncATREntryCUFr == AtrialResponseGbl \setminus (atrEntryC)$
 $AtrialResponseIncATREntryCFr == \Delta AtrialResponseGbl \wedge \exists AtrialResponseIncATREntryCUFr$

$AtrialResponseIncATREntryC == AtrialResponseIncATREntryCInc$
 $\vee AtrialResponseIncATREntryCSame$

$AtrialResponseVIncATREntryC$ _____
 $\Delta AtrialResponseIncATREntryCFr$
 $AtrialResponseIncATREntryC$

$AtrialResponseDecATREntryCDec$ _____
 $\Delta AtrialResponseGbl$
 $atrEntryC > 1$
 $atrEntryC' = atrEntryC - 1$

AtrialResponseDecATREntryCSame _____

$\Delta AtrialResponseGbl$

$atrEntryC = 1$

$atrEntryC' = atrEntryC$

AtrialResponseDecATREntryCUFr == AtrialResponseGbl \ (atrEntryC)

AtrialResponseDecATREntryCFr == $\Delta AtrialResponseGbl \wedge \exists AtrialResponseDecATREntryCUFr$

AtrialResponseDecATREntryC == AtrialResponseDecATREntryCDec

$\vee AtrialResponseDecATREntryCSame$

AtrialResponseVDecATREntryC _____

AtrialResponseDecATREntryCFr

AtrialResponseDecATREntryC

AtrialResponseIncATRExitCInc _____

$\Delta AtrialResponseGbl$

$atrExitC < 8$

$atrExitC' = atrExitC + 1$

AtrialResponseIncATRExitCSame _____

$\Delta AtrialResponseGbl$

$atrExitC = 8$

$atrExitC' = atrExitC$

AtrialResponseIncATRExitCUFr == AtrialResponseGbl \ (atrExitC)

AtrialResponseIncATRExitCFr == $\Delta AtrialResponseGbl \wedge \exists AtrialResponseIncATRExitCUFr$

AtrialResponseIncATRExitC == AtrialResponseIncATRExitCInc

$\vee AtrialResponseIncATRExitCSame$

AtrialResponseVI *IncATRExitC*

AtrialResponseIncATRExitCFr

AtrialResponseIncATRExitC

AtrialResponseDecATRExitCDec

$\Delta AtrialResponseGbl$

$atrExitC > 1$

$atrExitC' = atrExitC - 1$

AtrialResponseDecATRExitCSame

$\Delta AtrialResponseGbl$

$atrExitC = 1$

$atrExitC' = atrExitC$

AtrialResponseDecATRExitCUFr == *AtrialResponseGbl* \ ($atrExitC$)

AtrialResponseDecATRExitCFr == $\Delta AtrialResponseGbl \wedge \exists AtrialResponseDecATRExitCUFr$

AtrialResponseDecATRExitC == *AtrialResponseDecATRExitCDec*
 $\vee AtrialResponseDecATRExitCSame$

AtrialResponseVDecATRExitC

AtrialResponseDecATRExitCFr

AtrialResponseDecATRExitC

C.11 Package VentricularResponse

section *VentricularResponse* **parents** *ZOO_Toolkit*, *Model_Preamble*, *CommonPM*, *ResponseCommon*

PVARPExtension == { $o : MilliSec \mid o \geq 0 \wedge o \leq 400$ }

VentricularResponseGblSt

vRP : *RefractoryPeriod*

pvarp : *RefractoryPeriod*

pvarpExt : *PVARPExtension*

pvc : *Bool*

```
  └── VentricularResponseGbl _____  
    └── VentricularResponseGblSt
```

```
  └── VentricularResponseInit _____  
    └── VentricularResponseGbl  
      └── vrp = 320  
      └── pvarp = 250  
      └── pvarpExt = 0  
      └── pvc = False
```

```
  └── VentricularResponseGetPVARP _____  
    └── VentricularResponseGbl  
      └── pvarp! : RefractoryPeriod  
      └── pvc = False ⇒ pvarp! = pvarp  
      └── pvc = True ⇒ pvarp! = pvarp + pvarpExt
```

```
  └── VentricularResponseVGetPVARP _____  
    └── VentricularResponseGetPVARP
```

```
  └── VentricularResponseIsPVC _____  
    └── VentricularResponseGbl  
      └── cs? : ChambersSensed  
      └── lastAEvTm? : MilliSec  
      └── lastVEvTm? : MilliSec  
      └── cs? = D ∧ lastVEvTm? > lastAEvTm?
```

```
  └── VentricularResponseVIspVC _____  
    └── VentricularResponseIsPVC
```

```
  └── VentricularResponseIncVRPInc _____  
    └── Δ VentricularResponseGbl  
      └── vrp < 500  
      └── vrp' = vrp + 10
```

$\neg \text{VentricularResponseIncVRPSame}$ _____

$\Delta \text{VentricularResponseGbl}$

$vRP = 500$

$vRP' = vRP$

$\text{VentricularResponseIncVRPUFr} == \text{VentricularResponseGbl} \setminus (vRP)$

$\text{VentricularResponseIncVRPFr} == \Delta \text{VentricularResponseGbl} \wedge \exists \text{VentricularResponseIncVRPUFr}$

$\text{VentricularResponseIncVRP} == \text{VentricularResponseIncVRPInc}$

$\vee \text{VentricularResponseIncVRPSame}$

$\neg \text{VentricularResponseVIncVRP}$ _____

$\text{VentricularResponseIncVRPFr}$

$\text{VentricularResponseIncVRP}$

$\neg \text{VentricularResponseDecVRPDec}$ _____

$\Delta \text{VentricularResponseGbl}$

$vRP > 150$

$vRP' = vRP - 10$

$\neg \text{VentricularResponseDecVRPSame}$ _____

$\Delta \text{VentricularResponseGbl}$

$vRP = 150$

$vRP' = vRP$

$\text{VentricularResponseDecVRPUFr} == \text{VentricularResponseGbl} \setminus (vRP)$

$\text{VentricularResponseDecVRPFr} == \Delta \text{VentricularResponseGbl} \wedge \exists \text{VentricularResponseDecVRPUFr}$

$\text{VentricularResponseDecVRP} == \text{VentricularResponseDecVRPDec}$

$\vee \text{VentricularResponseDecVRPSame}$

$\text{VentricularResponseVDecVRP}$
 $\text{VentricularResponseDecVRPFr}$
 $\text{VentricularResponseDecVRP}$

$\text{VentricularResponseIncPVARPInc}$
 $\Delta \text{VentricularResponseGbl}$
 $pvarp < 500$
 $pvarp' = pvarp + 10$

$\text{VentricularResponseIncPVARPSame}$
 $\Delta \text{VentricularResponseGbl}$
 $pvarp = 500$
 $pvarp' = pvarp$

$\text{VentricularResponseIncPVARPUFr} == \text{VentricularResponseGbl} \setminus (pvarp)$

$\text{VentricularResponseIncPVARPFr} == \Delta \text{VentricularResponseGbl} \wedge \exists \text{VentricularResponseIncPVARPUFr}$

$\text{VentricularResponseIncPVARP} == \text{VentricularResponseIncPVARPInc}$
 $\vee \text{VentricularResponseIncPVARPSame}$

$\text{VentricularResponseVIncPVARP}$
 $\text{VentricularResponseIncPVARPFr}$
 $\text{VentricularResponseIncPVARP}$

$\text{VentricularResponseDecPVARPDec}$
 $\Delta \text{VentricularResponseGbl}$
 $pvarp > 150$
 $pvarp' = pvarp - 10$

$\text{VentricularResponseDecPVARPSame} \quad \dots$

$\Delta \text{VentricularResponseGbl}$

$pvarp = 150$

$pvarp' = pvarp$

$\text{VentricularResponseDecPVARPUFr} == \text{VentricularResponseGbl} \setminus (pvarp)$

$\text{VentricularResponseDecPVARPFr} == \Delta \text{VentricularResponseGbl} \wedge \exists \text{VentricularResponseDecPVARPUFr}$

$\text{VentricularResponseDecPVARP} == \text{VentricularResponseDecPVARPDec}$

$\vee \text{VentricularResponseDecPVARPSame}$

$\text{VentricularResponseVDecPVARP} \quad \dots$

$\text{VentricularResponseDecPVARPFr}$

$\text{VentricularResponseDecPVARP}$

$\text{VentricularResponseIncPVARPExtInc} \quad \dots$

$\Delta \text{VentricularResponseGbl}$

$pvarpExt < 400$

$pvarpExt' = pvarpExt + 50$

$\text{VentricularResponseIncPVARPExtSame} \quad \dots$

$\Delta \text{VentricularResponseGbl}$

$pvarpExt = 400$

$pvarpExt' = pvarpExt$

$\text{VentricularResponseIncPVARPExtUFr} == \text{VentricularResponseGbl} \setminus (pvarpExt)$

$\text{VentricularResponseIncPVARPExtFr} == \Delta \text{VentricularResponseGbl} \wedge \exists \text{VentricularResponseIncPVARPExtUFr}$

$\text{VentricularResponseIncPVARPExt} == \text{VentricularResponseIncPVARPExtInc}$

$\vee \text{VentricularResponseIncPVARPExtSame}$

$\text{VentricularResponseVI}_{\text{ncPVARPEst}}$ _____
 $\Delta \text{VentricularResponseIncPVARPEstFr}$
 $\text{VentricularResponseIncPVARPEst}$

$\text{VentricularResponseDecPVARPEstDec}$ _____
 $\Delta \text{VentricularResponseGbl}$
 $pvarpExt > 0$
 $pvarpExt' = pvarpExt - 50$

$\text{VentricularResponseDecPVARPEstSame}$ _____
 $\Delta \text{VentricularResponseGbl}$
 $pvarpExt = 0$
 $pvarpExt' = pvarpExt$

$\text{VentricularResponseDecPVARPEstUfr} == \text{VentricularResponseGbl} \setminus (pvarpExt)$
 $\text{VentricularResponseDecPVARPEstFr} == \Delta \text{VentricularResponseGbl} \wedge \exists \text{VentricularResponseDecPVARPEstUfr}$

$\text{VentricularResponseDecPVARPEst} == \text{VentricularResponseDecPVARPEstDec}$
 $\vee \text{VentricularResponseDecPVARPEstSame}$

$\text{VentricularResponseVDecPVARPEst}$ _____
 $\Delta \text{VentricularResponseDecPVARPEstFr}$
 $\text{VentricularResponseDecPVARPEst}$

$\text{VentricularResponsePVCDetected}$ _____
 $\Delta \text{VentricularResponseGbl}$
 $\text{VentricularResponseIsPVC}$
 $pvc = \text{False}$
 $pvc' = \text{True}$

$\text{VentricularResponseDisablePVC}$ _____
 $\Delta \text{VentricularResponseGbl}$
 $pvc' = \text{False}$

$VentricularResponsePVCUPdate == VentricularResponsePVCDetected$
 $\vee VentricularResponseDisablePVC$

$VentricularResponseVSenseOnly$ _____

$\Delta VentricularResponseGbl$
 $VentricularResponsePVCUPdate$
 $lastVEvTm? : MilliSec$
 $now? : MilliSec$

$vRP < now? - lastVEvTm?$

$VentricularResponseVSenseOnlyUFr == VentricularResponseGbl \setminus (pvc)$

$VentricularResponseVSenseOnlyFr == \Delta VentricularResponseGbl \wedge \exists VentricularResponseVSenseOnlyUFr$

$VentricularResponseVVSenseOnly$ _____

$VentricularResponseVSenseOnlyFr$
 $VentricularResponseVSenseOnly$

$VentricularResponseVSenseRefractory$ _____

$\Delta VentricularResponseGbl$
 $VentricularResponsePVCUPdate$
 $lastVEvTm? : MilliSec$
 $now? : MilliSec$

$vRP \geq now? - lastVEvTm?$

$VentricularResponseVSenseRefractoryUFr == VentricularResponseGbl \setminus (pvc)$

$VentricularResponseVSenseRefractoryFr == \Delta VentricularResponseGbl \wedge \exists VentricularResponseVSenseRefractoryUFr$

$VentricularResponseVVSenseRefractory$ _____

$VentricularResponseVSenseRefractoryFr$
 $VentricularResponseVSenseRefractory$

C.12 Package Response

section *Response* **parents** *ZOO_Toolkit*, *Model_Preamble*, *CommonPM*, *VentricularResponse*, *AtrialResponse*

ResponsePkgSt _____
 |rts : *ResponseToSensing*
 |prts : *ResponseToSensing*

ResponseOnlyGblSt _____
 |*ResponsePkgSt*

ResponseGblSt _____
 |*ResponseOnlyGblSt*
 |*VentricularResponseGbl*
 |*AtrialResponseGbl*

ResponseGbl _____
 |*ResponseGblSt*

ResponseInit _____
 |*AtrialResponseInit*
 |*VentricularResponseInit*
 |*ResponseGbl*

 |rts = D
 |prts = D

ResponseIntExtAtrialResponseFr == Δ *ResponseGbl* \wedge \exists *ResponseOnlyGblSt* \wedge \exists *VentricularResponseGbl*

ResponseVInATRFallBackMode _____
 |*ResponseGbl*
 |*AtrialResponseVInATRFallBackMode*

ResponseVIIncURL

ResponseIntExtAtrialResponseFr
AtrialResponseVIIncURL

ResponseVDecURL

ResponseIntExtAtrialResponseFr
AtrialResponseVDecURL

ResponseVIIncARP

ResponseIntExtAtrialResponseFr
AtrialResponseVIIncARP

ResponseVDecARP

ResponseIntExtAtrialResponseFr
AtrialResponseVDecARP

ResponseVATROn

ResponseIntExtAtrialResponseFr
AtrialResponseVATROn

ResponseVATROff

ResponseIntExtAtrialResponseFr
AtrialResponseVATROff

ResponseVIIncATRDURATION

ResponseIntExtAtrialResponseFr
AtrialResponseVIIncATRDURATION

ResponseVDecATRDURATION

ResponseIntExtAtrialResponseFr
AtrialResponseVDecATRDURATION

Response VIncATRFbTm _____
| *ResponseIntExtAtrialResponseFr*
| *AtrialResponseVIncATRDuratio*n
| _____

Response VDecATRFbTm _____
| *ResponseIntExtAtrialResponseFr*
| *AtrialResponseVDecATRFbTm*
| _____

Response VIncATREntryC _____
| *ResponseIntExtAtrialResponseFr*
| *AtrialResponseVIncATREntryC*
| _____

Response VDecATREntryC _____
| *ResponseIntExtAtrialResponseFr*
| *AtrialResponseVDecATREntryC*
| _____

Response VIncATRExitC _____
| *ResponseIntExtAtrialResponseFr*
| *AtrialResponseVIncATRExitC*
| _____

Response VDecATRExitC _____
| *ResponseIntExtAtrialResponseFr*
| *AtrialResponseVDecATRExitC*
| _____

ResponseIntExtVentricularResponseFr == $\Delta ResponseGbl \wedge \exists ResponseOnlyGblSt \wedge \exists AtrialResponseGbl$

Response VIsPVC _____
| *ResponseGbl*
| *VentricularResponseVIsPVC*
| _____

Response VIncVRP _____
| *ResponseIntExtVentricularResponseFr*
| *VentricularResponseVIncVRP*
| _____

Response VDecVRP

ResponseIntExtVentricularResponseFr
VentricularResponseVIncVRP

Response VIncPVARP

ResponseIntExtVentricularResponseFr
VentricularResponseVIncPVARP

Response VDecPVARP

ResponseIntExtVentricularResponseFr
VentricularResponseVIncPVARP

Response VIncPVARPExt

ResponseIntExtVentricularResponseFr
VentricularResponseVIncPVARPExt

Response VDecPVARPExt

ResponseIntExtVentricularResponseFr
VentricularResponseVIncPVARPExt

Response VVSenseRefractory

ResponseIntExtVentricularResponseFr
VentricularResponseVVSenseRefractory

ResponseSetBradycardiaMode

Δ *ResponseGbl*

rts? : ResponseToSensing

rts' = rts?

prts' = rts

ResponseSetBradycardiaModeUFr == ResponseGbl \ (rts, prts)

ResponseSetBradycardiaModeFr == \Delta ResponseGbl \wedge \exists ResponseSetBradycardiaModeUFr

Response VSetBradycardiaMode _____
 ┌── *ResponseSetBradycardiaModeFr*
 ┌── *ResponseSetBradycardiaMode*

ResponseIsBatteryDegraded _____
 ┌── *ResponseGblSt*
 ┌── *bs? : BatteryStatus*
 ┌── *bs? = ERT* \vee *bs? = ERP*

ResponseUpdateBatteryChangeMode _____
 ┌── Δ *ResponseGblSt*
 ┌── *ResponseIsBatteryDegraded*
 ┌── *ResponseSetBradycardiaMode[I/rts?]*

ResponseUpdateBatterySame _____
 ┌── Δ *ResponseGblSt*
 ┌── *rts' = rts*
 ┌── *prts' = prts*

ResponseUpdateBatteryStatus _____
 ┌── Δ *ResponseGbl*
 ┌── *bs? : BatteryStatus*
 ┌── *ResponseUpdateBatteryChangeMode* \vee *ResponseUpdateBatterySame*

ResponseUpdateBatteryStatusUFr == ResponseGbl \ (rts, prts)
ResponseUpdateBatteryStatusFr == \Delta ResponseGbl \ \exists ResponseUpdateBatteryStatusUFr

Response VUpdateBatteryStatus _____
 ┌── *ResponseUpdateBatteryStatusFr*
 ┌── *ResponseUpdateBatteryStatus*

ResponseExecuteEventATRFallback _____

$\Delta ResponseGbl$
 $ev? : PMEvent$
 $ResponseSetBradycardiaMode[I/rts?]$
 $ev? = EvATRFallback$

ResponseExecuteEventTerminateATR _____

$\Delta ResponseGbl$
 $ev? : PMEvent$
 $ev? = EvTerminateATR$
 $rts' = prts$
 $prts' = rts$

ResponseEvsDoingNothing _____

$ResponseGbl$
 $ev? : PMEvent$
 $ev? = EvFallbackTmElapsed \vee ev? = EvNone$

ResponseExecuteEventDoingNothing _____

$\Delta ResponseGbl$
ResponseEvsDoingNothing
 $ev? : PMEvent$
 $rts' = rts$
 $rts' = rts$

ResponseExecuteEvent == ResponseExecuteEventATRFallback

$\vee ResponseExecuteEventTerminateATR$
 $\vee ResponseExecuteEventDoingNothing$

ResponseExecuteEventUFr == ResponseGbl \setminus (rts, prts)

ResponseExecuteEventFr == \Delta ResponseGbl \wedge \exists ResponseExecuteEventUFr

ExecuteEvent _____

ResponseExecuteEventFr
ResponseExecuteEvent

ResponseASenseRefractory

$\Delta ResponseGbl$
 $now? : MilliSec$
 $lastAEvTm? : MilliSec$
 $lastVEvTm? : MilliSec$
 $VentricularResponseVGetPVARP$
 $AtrialResponseVASenseRefractory[pvarp!/pvarp?]$

$ResponseASenseRefractoryUFr == ResponsePkgSt \wedge VentricularResponseGbl$
 $ResponseASenseRefractoryFr == \Delta ResponseGbl \wedge \exists ResponseExecuteEventUFr$

ResponseVASenseRefractory

$ResponseASenseRefractoryFr$
 $ResponseASenseRefractory$

ResponseASenseOnly

$\Delta ResponseGbl$
 $AtrialResponseVASenseOnly$

$rts = O$

$ResponseASenseOnlyUFr == ResponsePkgSt \wedge VentricularResponseGbl$
 $ResponseASenseOnlyFr == \Delta ResponseGbl \wedge \exists ResponseASenseOnlyUFr$

ResponseVASenseOnly

$ResponseASenseOnlyFr$
 $ResponseASenseOnly$

ResponseVSenseOnly

$\Delta ResponseGbl$
 $VentricularResponseVVSenseOnly$

$rts = O$

$ResponseVSenseOnlyUFr == ResponsePkgSt \wedge AtrialResponseGbl$
 $ResponseVSenseOnlyFr == \Delta ResponseGbl \wedge \exists ResponseVSenseOnlyUFr$

Response VVSenseOnly _____
 ┌── *Response VSenseOnlyFr*
 ┌── *Response VSenseOnly*

Response TrackASense _____
 ┌── Δ *ResponseGbl*
 ┌── *AtrialResponse VASenseOnly*
 ┌── *rts = D*

Response TrackASenseUFr == ResponsePkgSt \wedge *VentricularResponseGbl*
Response TrackASenseFr == Δ ResponseGbl \wedge \exists *Response TrackASenseUFr*

Response VTrackASense _____
 ┌── *Response TrackASenseFr*
 ┌── *Response TrackASense*

Response TriggerAPace _____
 ┌── Δ *ResponseGbl*
 ┌── *AtrialResponse VASenseOnly*
 ┌── *rts = T*

Response TriggerAPaceUFr == ResponsePkgSt \wedge *VentricularResponseGbl*
Response TriggerAPaceFr == Δ ResponseGbl \wedge \exists *Response TriggerAPaceUFr*

Response VTriggerAPace _____
 ┌── *Response TriggerAPaceFr*
 ┌── *Response TriggerAPace*

Response TriggerVPace _____
 ┌── Δ *ResponseGbl*
 ┌── *VentricularResponse VVSenseOnly*
 ┌── *rts = T*

$\text{ResponseTriggerVPaceUFr} == \text{ResponsePkgSt} \wedge \text{AtrialResponseGbl}$
 $\text{ResponseTriggerVPaceFr} == \Delta \text{ResponseGbl} \wedge \exists \text{ResponseTriggerVPaceUFr}$

$\text{ResponseVTriggerVPace} \ldots$

$\text{ResponseTriggerVPaceFr}$
$\text{ResponseTriggerVPace}$

$\text{ResponseInhibitAPace} \ldots$

$\Delta \text{ResponseGbl}$
$\text{AtrialResponseVASenseOnly}$

$rts = I$

$\text{ResponseInhibitAPaceUFr} == \text{ResponsePkgSt} \wedge \text{VentricularResponseGbl}$
 $\text{ResponseInhibitAPaceFr} == \Delta \text{ResponseGbl} \wedge \exists \text{ResponseInhibitAPaceUFr}$

$\text{ResponseVIInhibitAPace} \ldots$

$\text{ResponseInhibitAPaceFr}$
$\text{ResponseInhibitAPace}$

$\text{ResponseInhibitedOrTracked} \ldots$

ResponseGbl

$rts = I \vee rts = T$

$\text{ResponseInhibitVPace} \ldots$

$\Delta \text{ResponseGbl}$
$\text{ResponseInhibitedOrTracked}$
$\text{VentricularResponseVVSenseOnly}$

$\text{ResponseInhibitVPaceUFr} == \text{ResponsePkgSt} \wedge \text{AtrialResponseGbl}$
 $\text{ResponseInhibitVPaceFr} == \Delta \text{ResponseGbl} \wedge \exists \text{ResponseInhibitVPaceUFr}$

$\text{ResponseVIInhibitVPace} \ldots$

$\text{ResponseInhibitVPaceFr}$
$\text{ResponseInhibitVPace}$

C.13 Package HistoryCommon

```
section HistoryCommon parents ZOO_Toolkit, Model_Preamble
```

```
AugmentationMarker ::= ATRDur | ATRFB | ATREnd | PVP | AugMNone
```

```
VentricularMarker ::= VS | VP | PVC | VTN
```

```
AtrialMarker ::= AS | AP | AT | ATN
```

C.14 Package History

```
section History parents ZOO_Toolkit, Model_Preamble, HistoryCommon, CommonPM
```

```
HistoryEvent0
```

```
time : MilliSec  
augMarker : AugmentationMarker  
mRe : Bool  
mHy : Bool  
mSr : Bool  
mUp : Bool  
mDown : Bool
```

```
HistoryEvent
```

```
HistoryEvent0
```

```
SHistoryEvent
```

```
sHistoryEvent : P(∅ HistoryEventCl)  
stHistoryEvent : ∅ HistoryEventCl → HistoryEvent
```

```
dom stHistoryEvent = sHistoryEvent  
sHistoryEvent ∩ ∅x HistoryEventCl = ∅
```

```
SHistoryEventInit
```

```
SHistoryEvent
```

```
sHistoryEvent = ∅ ∧ stHistoryEvent = ∅
```

HistoryEventNew _____

HistoryEvent
now? : *MilliSec*
mre? : *Bool*
mhy? : *Bool*
msr? : *Bool*
mUp? : *Bool*
mDown? : *Bool*

time = *now?*
mRe = *mre?*
mHy = *mhy?*
mSr = *msr?*
mUp = *mUp?*
mDown = *mDown?*

SHistoryEventNF _____

$\Delta S\text{HistoryEvent}$
HistoryEvent
o! : $\textcircled{\textcircled{}} \text{HistoryEventCl}$

$s\text{HistoryEvent}' = s\text{HistoryEvent}' \cup \{ o! \}$
 $st\text{HistoryEvent}' = st\text{HistoryEvent}' \cup \{ (o!, \theta \text{ HistoryEvent}) \}$

AtrialEvent0 _____

HistoryEvent
aMarker : *AtrialMarker*

AtrialEvent _____

AtrialEvent0

SAtrialEvent _____

$s\text{AtrialEvent} : \mathbb{P}(\textcircled{\textcircled{}} \text{AtrialEventCl})$
 $st\text{AtrialEvent} : \textcircled{\textcircled{}} \text{AtrialEventCl} \rightarrow \text{AtrialEvent}$

$\text{dom } st\text{AtrialEvent} = s\text{AtrialEvent}$

*S*AtrialEventInit _____
 ┌── *S*AtrialEvent
 └── $sAtrialEvent = \emptyset \wedge stAtrialEvent = \emptyset$

*S*AtrialEventIsHistoryEvent _____
 ┌── *S*HistoryEvent; *S*AtrialEvent
 └── $sAtrialEvent \subseteq sHistoryEvent$
 $\forall o : sAtrialEvent \bullet$
 $(\lambda AtrialEvent \bullet \theta HistoryEvent)(stAtrialEvent o) = (stHistoryEvent o)$

AtrialEventNew _____
 ┌── *AtrialEvent*
 ┌── *HistoryEventNew*
 ┌── *am?* : *AtrialMarker*
 └── $aMarker = am?$

*S*AtrialEventNF _____
 ┌── ΔS AtrialEvent
 ┌── *S*HistoryEventNF
 ┌── *AtrialEvent*
 ┌── $o! : \bigcirc AtrialEventCl$
 └── $o! \in \bigcirc x AtrialEventCl \setminus sAtrialEvent$
 $sAtrialEvent' = sAtrialEvent' \cup \{o!\}$
 $stAtrialEvent' = stAtrialEvent' \cup \{(o!, \theta AtrialEvent)\}$

S AtrialEventNew == $\exists AtrialEvent \bullet$
 S AtrialEventNF[ae!/o!] \wedge AtrialEventNew

VentricularEvent0 _____
 ┌── *HistoryEvent*
 ┌── *vMarker* : *VentricularMarker*

VentricularEvent _____
VentricularEvent0

— *SVentricularEvent* _____
sVentricularEvent : $\mathbb{P}(\mathbb{O} \text{ VentricularEventCl})$
stVentricularEvent : $\mathbb{O} \text{ VentricularEventCl} \rightarrow \text{VentricularEvent}$
 dom *stVentricularEvent* = *sVentricularEvent*

SVentricularEventInit _____
SVentricularEvent
sVentricularEvent = $\emptyset \wedge \text{stVentricularEvent} = \emptyset$

SVentricularEventIsHistoryEvent _____
SHistoryEvent; SVentricularEvent
sVentricularEvent \subseteq *sHistoryEvent*
 $\forall o : sVentricularEvent \bullet$
 $(\lambda \text{ VentricularEvent } \bullet \theta \text{ HistoryEvent})(\text{stVentricularEvent } o) = (\text{sHistoryEvent } o)$

VentricularEventNew _____
VentricularEvent
HistoryEventNew
vm? : VentricularMarker
vMarker = *vm?*

SVentricularEventNF _____
 $\Delta \text{SVentricularEvent}$
SHistoryEventNF
VentricularEvent
 $o! : \mathbb{O} \text{ VentricularEventCl}$
 $o! \in \mathbb{O} \text{ VentricularEventCl} \setminus sVentricularEvent$
 $sVentricularEvent' = sVentricularEvent' \cup \{o!\}$
 $\text{stVentricularEvent}' = \text{stVentricularEvent}' \cup \{(o!, \theta \text{ VentricularEvent})\}$

SVentricularEventNew == $\exists \text{ VentricularEvent} \bullet$
 $SVentricularEventNF[ve!/o!] \wedge \text{VentricularEventNew}$

HistoryGblSt _____
| *SVentricularEvent*
| *SAtrialEvent*
| *SHistoryEvent*
| _____

HistoryGbl _____
| *HistoryGblSt*
| *SAtrialEventIsHistoryEvent*
| *SVentricularEventIsHistoryEvent*
| _____

HistoryCreateHistoryAEv0 _____
| Δ *HistoryGbl*
| *SAtrialEventNew*
| _____

HistoryCreateHistoryAEv == *HistoryCreateHistoryAEv0* \ (ae!)

HistoryVCreateHistoryAEv _____
| Δ *HistoryGbl*
| Ξ *SVentricularEvent*
| *HistoryCreateHistoryAEv*
| _____

HistoryCreateHistoryVEv0 _____
| Δ *HistoryGbl*
| *SVentricularEventNew*
| _____

HistoryCreateHistoryVEv == *HistoryCreateHistoryVEv0* \ (ve!)

HistoryVCreateHistoryVEv _____
| Δ *HistoryGbl*
| Ξ *SAtrialEvent*
| *HistoryCreateHistoryVEv*
| _____

C.15 Package PacingWithHistoryI

section *PacingWithHistoryI* **parents** *ZOO_Toolkit*, *Model_Preamble*, *HistoryCommon*, *CommonPM*, *Pacing*

— *PacingWithHistoryIGblSt* _____
 PacingGbl

— *PacingWithHistoryIGbl* _____
 PacingWithHistoryIGblSt

*PacingWithHistoryI**IntExtPacingFr* == Δ *PacingWithHistoryIGbl*

— *PacingWithHistoryIInit* _____
 *PacingWithHistoryI**IntExtPacingFr*
 PacingInit

— *PacingWithHistoryIVIncRSDown* _____
 *PacingWithHistoryI**IntExtPacingFr*
 PacingVIncRSDown

— *PacingWithHistoryIVDecRSDown* _____
 *PacingWithHistoryI**IntExtPacingFr*
 PacingVDecRSDown

— *PacingWithHistoryIVIncRSUp* _____
 *PacingWithHistoryI**IntExtPacingFr*
 PacingVIncRSUp

— *PacingWithHistoryIVDecRSUp* _____
 *PacingWithHistoryI**IntExtPacingFr*
 PacingVDecRSUp

PacingWithHistoryIVDynamicDelayOn _____
PacingWithHistoryIIntExtPacingFr
PacingVDynamicDelayOn

PacingWithHistoryIVDynamicDelayOff _____
PacingWithHistoryIIntExtPacingFr
PacingVDynamicDelayOff

PacingWithHistoryIVIncFixedAVI _____
PacingWithHistoryIIntExtPacingFr
PacingVIncFixedAVI

PacingWithHistoryIVDecFixedAVI _____
PacingWithHistoryIIntExtPacingFr
PacingVDecFixedAVI

PacingWithHistoryIVIncMinDynAVDelay _____
PacingWithHistoryIIntExtPacingFr
PacingVIncMinDynAVDelay

PacingWithHistoryIVDecMinDynAVDelay _____
PacingWithHistoryIIntExtPacingFr
PacingVDecMinDynAVDelay

PacingWithHistoryIVIncMaxDynAVDelay _____
PacingWithHistoryIIntExtPacingFr
PacingVIncMaxDynAVDelay

PacingWithHistoryIVDecMaxDynAVDelay _____
PacingWithHistoryIIntExtPacingFr
PacingVDecMaxDynAVDelay

PacingWithHistoryIVIncSAVDelayOffset _____

PacingWithHistoryIIntExtPacingFr

PacingVIIncSAVDelayOffset

PacingWithHistoryIVDecSAVDelayOffset _____

PacingWithHistoryIIntExtPacingFr

PacingVDecSAVDelayOffset

PacingWithHistoryIVSetBradycardiaMode _____

PacingWithHistoryIIntExtPacingFr

PacingVSetBradycardiaMode

PacingWithHistoryIVUpdateBatteryStatus _____

PacingWithHistoryIIntExtPacingFr

PacingVUpdateBatteryStatus

PacingWithHistoryIVIncAPulseAmplReg _____

PacingWithHistoryIIntExtPacingFr

PacingVIIncAPulseAmplReg

PacingWithHistoryIVDecAPulseAmplReg _____

PacingWithHistoryIIntExtPacingFr

PacingVDecAPulseAmplReg

PacingWithHistoryIVIncVPulseAmplReg _____

PacingWithHistoryIIntExtPacingFr

PacingVIIncVPulseAmplReg

PacingWithHistoryIVDecVPulseAmplReg _____

PacingWithHistoryIIntExtPacingFr

PacingVDecVPulseAmplReg

PacingWithHistoryIVIncAPulseAmplUReg

PacingWithHistoryIIntExtPacingFr

PacingVIIncAPulseAmplUReg

PacingWithHistoryIVDecAPulseAmplUReg

PacingWithHistoryIIntExtPacingFr

PacingVDecAPulseAmplUReg

PacingWithHistoryIVIncVPulseAmplUReg

PacingWithHistoryIIntExtPacingFr

PacingVIIncVPulseAmplUReg

PacingWithHistoryIVDecVPulseAmplUReg

PacingWithHistoryIIntExtPacingFr

PacingVDecVPulseAmplUReg

PacingWithHistoryIVIncAPulseWidth

PacingWithHistoryIIntExtPacingFr

PacingVIIncAPulseWidth

PacingWithHistoryIVDecAPulseWidth

PacingWithHistoryIIntExtPacingFr

PacingVDecAPulseWidth

PacingWithHistoryIVIncVPulseWidth

PacingWithHistoryIIntExtPacingFr

PacingVIIncVPulseWidth

PacingWithHistoryIVDecVPulseWidth

PacingWithHistoryIIntExtPacingFr

PacingVDecVPulseWidth

PacingWithHistoryIVIncLRL
PacingWithHistoryIIntExtPacingFr
PacingVIIncLRL

PacingWithHistoryIVDecLRL
PacingWithHistoryIIntExtPacingFr
PacingVDecLRL

PacingWithHistoryIVIncHRL
PacingWithHistoryIIntExtPacingFr
PacingVIIncHRL

PacingWithHistoryIVDecHRL
PacingWithHistoryIIntExtPacingFr
PacingVDecHRL

PacingWithHistoryIVSetTargetCCI
PacingWithHistoryIIntExtPacingFr
PacingVSetTargetCCI

PacingWithHistoryIVTriggerAPace
PacingWithHistoryIIntExtPacingFr
PacingVTriggerAPace

PacingWithHistoryIVTriggerVPace
PacingWithHistoryIIntExtPacingFr
PacingVTriggerVPace

PacingWithHistoryIVInhibitAPace
PacingWithHistoryIIntExtPacingFr
PacingVInhibitAPace

PacingWithHistoryIVInhibitVPace

PacingWithHistoryIIntExtPacingFr

PacingVInhibitVPace

PacingWithHistoryIVTrackASense

PacingWithHistoryIIntExtPacingFr

PacingVTrackASense

PacingWithHistoryIVExecuteEvent

PacingWithHistoryIIntExtPacingFr

PacingVExecuteEvent

PacingWithHistoryISetEvMarkersCommon

PacingWithHistoryIGbl

augm! : AugmentationMarker

mRe! : Bool

mhy! : Bool

msr! : Bool

mUp! : Bool

mDown! : Bool

$(rm = R \wedge targetCCI < 60000 \text{ div } lrl) \Leftrightarrow msr! = True$

$(rm = O \wedge hystMode = On) \Leftrightarrow mhy! = True$

$targetCCI < currentCCI \Leftrightarrow mUp! = True$

$targetCCI > currentCCI \Leftrightarrow mDown! = True$

$augm! = AugMNone \wedge mRe! = False$

PacingWithHistoryISetEvMarkersA

PacingWithHistoryIGbl

PacingWithHistoryISetEvMarkersCommon

am! : AtrialMarker

$am! = AP$

PacingWithHistoryISetEvMarkersV _____
 | *PacingWithHistoryIGbl*
 | *PacingWithHistoryISetEvMarkersCommon*
 | *vm! : VentricularMarker*
 | *vm! = VP*

PacingWithHistoryIVTimedAPace _____
 | *PacingWithHistoryIGbl*
 | *PacingVTimedAPace*
 | *PacingWithHistoryISetEvMarkersA*

PacingWithHistoryIVTimedVPace _____
 | *PacingWithHistoryIGbl*
 | *PacingVTimedVPace*
 | *PacingWithHistoryISetEvMarkersV*

C.16 Package ResponseWithHistoryI

section *ResponseWithHistoryI* **parents** *ZOO_Toolkit*, *Model_Preamble*, *HistoryCommon*, *CommonPM*, *Respo*

ResponseWithHistoryIGblSt _____
 | *ResponseGbl*

ResponseWithHistoryIGbl _____
 | *ResponseWithHistoryIGblSt*

ResponseWithHistoryIIntExtResponseFr == ΔResponseWithHistoryIGbl

ResponseWithHistoryIInit _____
 | *ResponseWithHistoryIIntExtResponseFr*
 | *ResponseInit*

ResponseWithHistoryIVSetBradycardiaMode _____
 | *ResponseWithHistoryIIntExtResponseFr*
 | *ResponseVSetBradycardiaMode*

Response WithHistoryIVUpdateBatteryStatus _____

Response WithHistoryIIntExtResponseFr

Response VUpdateBatteryStatus

Response WithHistoryIVExecuteEvent _____

Response WithHistoryIIntExtResponseFr

Response VExecuteEvent

Response WithHistoryIVInATRFallBackMode _____

Response WithHistoryIGbl

Response VInATRFallBackMode

Response WithHistoryIVIncURL _____

Response WithHistoryIIntExtResponseFr

Response VIncURL

Response WithHistoryIVDecURL _____

Response WithHistoryIIntExtResponseFr

Response VDecURL

Response WithHistoryIVIncARP _____

Response WithHistoryIIntExtResponseFr

Response VIncARP

Response WithHistoryIVDecARP _____

Response WithHistoryIIntExtResponseFr

Response VDecARP

Response WithHistoryIVATROn _____

Response WithHistoryIIntExtResponseFr

Response VATROn

Response WithHistoryIVATROff _____
[*Response WithHistoryIIntExtResponseFr*
Response VATROff]

Response WithHistoryIVIncATRDuration _____
[*Response WithHistoryIIntExtResponseFr*
Response VIncATRDuration]

Response WithHistoryIVDecATRDuration _____
[*Response WithHistoryIIntExtResponseFr*
Response VDecATRDuration]

Response WithHistoryIVIncATRFbTm _____
[*Response WithHistoryIIntExtResponseFr*
Response VIncATRDuration]

Response WithHistoryIVDecATRFbTm _____
[*Response WithHistoryIIntExtResponseFr*
Response VDecATRFbTm]

Response WithHistoryIVIncATREntryC _____
[*Response WithHistoryIIntExtResponseFr*
Response VIncATREntryC]

Response WithHistoryIVDecATREntryC _____
[*Response WithHistoryIIntExtResponseFr*
Response VDecATREntryC]

Response WithHistoryIVIncATRExitC _____
[*Response WithHistoryIIntExtResponseFr*
Response VIncATRExitC]

Response WithHistoryIVDecATRExitC

Response WithHistoryIIntExtResponseFr

Response VDecATRExitC

Response WithHistoryIVIncVRP

Response WithHistoryIIntExtResponseFr

Response VIncVRP

Response WithHistoryIVDecVRP

Response WithHistoryIIntExtResponseFr

Response VDecVRP

Response WithHistoryIVIncPVARP

Response WithHistoryIIntExtResponseFr

Response VIncPVARP

Response WithHistoryIVDecPVARP

Response WithHistoryIIntExtResponseFr

Response VDecPVARP

Response WithHistoryIVIncPVARPExt

Response WithHistoryIIntExtResponseFr

Response VIncPVARPExt

Response WithHistoryIVDecPVARPExt

Response IntExtVentricularResponseFr

Response VDecPVARPExt

Response WithHistoryIASMarker

Response WithHistoryIGbl

am! : AtrialMarker

augm! : AugmentationMarker

ev? : PMEvent

(atrStatus = NoATR ∨ atrStatus = FallbackTmElapsed) ∨ (ev? = EvTerminateATR ∨ ev? = EvFallbackTr

((ev? = EvTerminateATR ∨ ev? = EvFallbackTmElapsed) ∧ augm! = ATREnd) ∨ augm! = AugMNone

am! = AS

ResponseWithHistoryIATMarker _____

ResponseWithHistoryIGbl

am! : *AtrialMarker*

augm! : *AugmentationMarker*

ev? : *PMEvent*

ev? = EvATRDuration \vee *atrStatus* \neq *NoATR*

am! = AT

(ev? = EvATRDuration \wedge *augm! = ATRDur*) \vee *augm! = AugMNone*

ResponseWithHistoryISetEvMarkersADF == ResponseWithHistoryIASMarker \vee *ResponseWithHistoryIATMarker*

ResponseWithHistoryISetEvMarkersA _____

ResponseWithHistoryIGbl

ResponseWithHistoryISetEvMarkersADF

mhy! : *Bool*

msr! : *Bool*

mUp! : *Bool*

mDown! : *Bool*

mhy! = False

msr! = False

mUp! = False

mDown! = False

ResponseWithHistoryIVSetEvMarkersA _____

ResponseWithHistoryIGbl

ResponseWithHistoryISetEvMarkersA

ResponseWithHistoryIAugMWithPVC _____

ResponseWithHistoryIGbl

ResponseVIsPVC

augm! : *AugmentationMarker*

augm! = PVP

ResponseWithHistoryIAugMWithoutPVC

ResponseWithHistoryIGbl

ResponseVIspvc

augm! : AugmentationMarker

augm! = AugMNone

ResponseWithHistoryISetEvMarkersVDF == ResponseWithHistoryIAugMWithPVC ∨ ResponseWithHistoryIAugMWithoutPVC

ResponseWithHistoryISetEvMarkersV

ResponseWithHistoryIGbl

ResponseWithHistoryISetEvMarkersVDF

vm! : VentricularMarker

mhy! : Bool

msr! : Bool

mUp! : Bool

mDown! : Bool

vm! = VS

mhy! = False

msr! = False

mUp! = False

mDown! = False

ResponseWithHistoryIVSetEvMarkersV

ResponseWithHistoryIGbl

ResponseWithHistoryISetEvMarkersV

ResponseWithHistoryIASenseRefractory

ΔResponseWithHistoryIGbl

ResponseVASenseRefractory

ResponseWithHistoryISetEvMarkersA

mre! : Bool

mre! = True

ResponseWithHistoryIASenseRefractoryUFr == ResponseWithHistoryIGbl

ResponseWithHistoryIASenseRefractoryFr == ΔResponseWithHistoryIGbl ∧ ∃ResponseWithHistoryIASenseR

Response WithHistoryIVASenseRefractory _____
 Response WithHistoryIASenseRefractoryFr
 Response WithHistoryIASenseRefractory

Response WithHistoryIVSenseRefractory _____
 Δ *Response WithHistoryIGbl*
 Response VVSenseRefractory
 Response WithHistoryISetEvMarkersV
 mre! : Bool
 mre! = True

Response WithHistoryIVSenseRefractoryUFr == Response WithHistoryIGbl

Response WithHistoryIVSenseRefractoryFr == Δ Response WithHistoryIGbl $\wedge \exists$ Response WithHistoryIVSenseR

Response WithHistoryIVVSenseRefractory _____
 Response WithHistoryIVSenseRefractoryFr
 Response WithHistoryIVSenseRefractory

Response WithHistoryIASenseOnly _____
 Δ *Response WithHistoryIGbl*
 Response VASenseOnly
 Response WithHistoryISetEvMarkersA
 mre! : Bool
 mre! = False

Response WithHistoryIASenseOnlyUFr == Response WithHistoryIGbl

Response WithHistoryIASenseOnlyFr == Δ Response WithHistoryIGbl $\wedge \exists$ Response WithHistoryIASenseOnlyU

Response WithHistoryIVASenseOnly _____
 Response WithHistoryIASenseOnlyFr
 Response WithHistoryIASenseOnly

ResponseWithHistoryIVSenseOnly

$\Delta \text{ResponseWithHistoryIGbl}$

ResponseVVSenseOnly

ResponseWithHistoryISetEvMarkersV

mre! : Bool

$mre! = \text{False}$

ResponseWithHistoryIVSenseOnlyUFr == ResponseWithHistoryIGbl

ResponseWithHistoryIVSenseOnlyFr == $\Delta \text{ResponseWithHistoryIGbl} \wedge \exists \text{ResponseWithHistoryIVSenseOnlyUFr}$

ResponseWithHistoryIVVSenseOnly

ResponseWithHistoryIVSenseOnlyFr

ResponseWithHistoryIVSenseOnly

ResponseWithHistoryITriggerAPace

$\Delta \text{ResponseWithHistoryIGbl}$

ResponseVTriggerAPace

ResponseWithHistoryISetEvMarkersA

mre! : Bool

$mre! = \text{False}$

ResponseWithHistoryITriggerAPaceUFr == ResponseWithHistoryIGbl

ResponseWithHistoryITriggerAPaceFr == $\Delta \text{ResponseWithHistoryIGbl} \wedge \exists \text{ResponseWithHistoryITriggerAPaceUFr}$

ResponseWithHistoryIVTriggerAPace

ResponseWithHistoryITriggerAPaceFr

ResponseWithHistoryITriggerAPace

ResponseWithHistoryITriggerVPace

$\Delta \text{ResponseWithHistoryIGbl}$

ResponseVTriggerVPace

ResponseWithHistoryISetEvMarkersV

mre! : Bool

$mre! = \text{False}$

ResponseWithHistoryITriggerVPaceUFr == ResponseWithHistoryIGbl

ResponseWithHistoryITriggerVPaceFr == ΔResponseWithHistoryIGbl ∧ ∃ResponseWithHistoryITriggerVPace

ResponseWithHistoryIVTriggerVPace

ResponseWithHistoryITriggerVPaceFr

ResponseWithHistoryITriggerVPace

ResponseWithHistoryIIInhibitAPace

ΔResponseWithHistoryIGbl

ResponseVInhibitAPace

ResponseWithHistoryISetEvMarkersA

mre! : Bool

mre! = False

ResponseWithHistoryIIInhibitAPaceUFr == ResponseWithHistoryIGbl

ResponseWithHistoryIIInhibitAPaceFr == ΔResponseWithHistoryIGbl ∧ ∃ResponseWithHistoryIIInhibitAPace

ResponseWithHistoryIVInhibitAPace

ResponseWithHistoryIIInhibitAPaceFr

ResponseWithHistoryIIInhibitAPace

ResponseWithHistoryIIInhibitVPace

ΔResponseWithHistoryIGbl

ResponseVInhibitVPace

ResponseWithHistoryISetEvMarkersV

mre! : Bool

mre! = False

ResponseWithHistoryIIInhibitVPaceUFr == ResponseWithHistoryIGbl

ResponseWithHistoryIIInhibitVPaceFr == ΔResponseWithHistoryIGbl ∧ ∃ResponseWithHistoryIIInhibitVPace

```
ResponseWithHistoryIVInhibitVPace _____
| ResponseWithHistoryIInhibitVPaceFr
| ResponseWithHistoryIInhibitVPace
```

```
ResponseWithHistoryITrackASense _____
| ΔResponseWithHistoryIGbl
| ResponseVTrackASense
| ResponseWithHistoryISetEvMarkersA
| mre! : Bool
| mre! = False
```

ResponseWithHistoryITrackASenseUFr == ResponseWithHistoryIGbl

ResponseWithHistoryITrackASenseFr == ΔResponseWithHistoryIGbl ∧ ∃ResponseWithHistoryITrackASenseU

```
ResponseWithHistoryIVTrackASense _____
| ResponseWithHistoryITrackASenseFr
| ResponseWithHistoryITrackASense
```

C.17 Package PacemakerWithoutHistory

section *PacemakerWithoutHistory* **parents** *ZOO_Toolkit*, *Model_Preamble*, *PacingWithHistoryI*, *Sensing*, *Ba*

```
PacemakerWithoutHistoryGblSt _____
| RateModulationGbl
| BatteryGbl
| SensingGbl
| PacingWithHistoryIGbl
| ResponseWithHistoryIGbl
```

```
PacemakerWithoutHistoryGbl _____
| PacemakerWithoutHistoryGblSt
```

PacemakerWithoutHistoryIntExtPacingWithHistoryIFr == ΔPacemakerWithoutHistoryGbl

PacemakerWithoutHistoryVIncRSDown _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncRSDown

PacemakerWithoutHistoryVDecRSDown _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecRSDown

PacemakerWithoutHistoryVIncRSUp _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncRSUp

PacemakerWithoutHistoryVDecRSUp _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecRSUp

PacemakerWithoutHistoryVDynamicDelayOn _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDynamicDelayOn

PacemakerWithoutHistoryVDynamicDelayOff _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDynamicDelayOff

PacemakerWithoutHistoryVIncFixedAVI _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncFixedAVI

PacemakerWithoutHistoryVDecFixedAVI _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecFixedAVI

PacemakerWithoutHistoryVIncMinDynAVDelay _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncMinDynAVDelay

PacemakerWithoutHistoryVDecMinDynAVDelay _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecMinDynAVDelay

PacemakerWithoutHistoryVIncMaxDynAVDelay _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncMaxDynAVDelay

PacemakerWithoutHistoryVDecMaxDynAVDelay _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecMaxDynAVDelay

PacemakerWithoutHistoryVIncSAVDelayOffset _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncSAVDelayOffset

PacemakerWithoutHistoryVDecSAVDelayOffset _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecSAVDelayOffset

PacemakerWithoutHistoryVIncAPulseAmplReg _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncAPulseAmplReg

PacemakerWithoutHistoryVDecAPulseAmplReg _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecAPulseAmplReg

PacemakerWithoutHistoryVIIncVPulseAmplReg
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncVPulseAmplReg

PacemakerWithoutHistoryVDecVPulseAmplReg
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecVPulseAmplReg

PacemakerWithoutHistoryVIIncAPulseAmplUReg
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncAPulseAmplUReg

PacemakerWithoutHistoryVDecAPulseAmplUReg
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecAPulseAmplUReg

PacemakerWithoutHistoryVIIncVPulseAmplUReg
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncVPulseAmplUReg

PacemakerWithoutHistoryVDecVPulseAmplUReg
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecVPulseAmplUReg

PacemakerWithoutHistoryVIIncAPulseWidth
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncAPulseWidth

PacemakerWithoutHistoryVDecAPulseWidth
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecAPulseWidth

PacemakerWithoutHistoryVIIncVPulseWidth _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncVPulseWidth

PacemakerWithoutHistoryVDecVPulseWidth _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecVPulseWidth

PacemakerWithoutHistoryVIIncLRL _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncLRL

PacemakerWithoutHistoryVDecLRL _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecLRL

PacemakerWithoutHistoryVIIncHRL _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVIncHRL

PacemakerWithoutHistoryVDecHRL _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVDecHRL

PacemakerWithoutHistoryVSetTargetCCI _____
PacemakerWithoutHistoryIntExtPacingWithHistoryIFr
PacingWithHistoryIVSetTargetCCI

PacemakerWithoutHistoryIntExtSensingFr == Δ *PacemakerWithoutHistoryGbl*

PacemakerWithoutHistoryVIIncASensitivity _____
PacemakerWithoutHistoryIntExtSensingFr
SensingVIIncASensitivity

PacemakerWithoutHistoryVDecASensitivity

PacemakerWithoutHistoryIntExtSensingFr

SensingVDecASensitivity

PacemakerWithoutHistoryVIncVSensitivity

PacemakerWithoutHistoryIntExtSensingFr

SensingVIncVSensitivity

PacemakerWithoutHistoryVDecVSensitivity

PacemakerWithoutHistoryIntExtSensingFr

SensingVDecVSensitivity

PacemakerWithoutHistoryVIncABlanking

PacemakerWithoutHistoryIntExtSensingFr

SensingVIncABlanking

PacemakerWithoutHistoryVDecABlanking

PacemakerWithoutHistoryIntExtSensingFr

SensingVDecABlanking

PacemakerWithoutHistoryVIncVBlanking

PacemakerWithoutHistoryIntExtSensingFr

SensingVIncVBlanking

PacemakerWithoutHistoryVDecVBlanking

PacemakerWithoutHistoryIntExtSensingFr

SensingVDecVBlanking

PacemakerWithoutHistoryIntExtResponseWithHistoryIFr == ΔPacemakerWithoutHistoryGbl

PacemakerWithoutHistoryVIIncURL _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVIncURL

PacemakerWithoutHistoryVDecURL _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVDecURL

PacemakerWithoutHistoryVIIncARP _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVIncARP

PacemakerWithoutHistoryVDecARP _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVDecARP

PacemakerWithoutHistoryVATROn _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVATROn

PacemakerWithoutHistoryVATROff _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVATROff

PacemakerWithoutHistoryVIIncATRDURATION _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVIncATRDURATION

PacemakerWithoutHistoryVDecATRDURATION _____
PacemakerWithoutHistoryIntExtResponseWithHistoryIFr
ResponseWithHistoryIVDecATRDURATION

PacemakerWithoutHistoryVIIncATRFbTm _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVIncATRFbTm] _____

PacemakerWithoutHistoryVDecATRFbTm _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVDecATRFbTm] _____

PacemakerWithoutHistoryVIIncATREEntryC _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVIncATREEntryC] _____

PacemakerWithoutHistoryVDecATREEntryC _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVDecATREEntryC] _____

PacemakerWithoutHistoryVIIncATRExitC _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVIncATRExitC] _____

PacemakerWithoutHistoryVDecATRExitC _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVDecATRExitC] _____

PacemakerWithoutHistoryVIIncVRP _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVIncVRP] _____

PacemakerWithoutHistoryVDecVRP _____
[*PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
ResponseWithHistoryIVDecVRP] _____

PacemakerWithoutHistoryVIIncPVARP _____
| *PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
| *ResponseWithHistoryIVIncPVARP*
| _____

PacemakerWithoutHistoryVDecPVARP _____
| *PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
| *ResponseWithHistoryIVDecPVARP*
| _____

PacemakerWithoutHistoryVIIncPVARPExt _____
| *PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
| *ResponseWithHistoryIVIncPVARPExt*
| _____

PacemakerWithoutHistoryVDecPVARPExt _____
| *PacemakerWithoutHistoryIntExtResponseWithHistoryIFr*
| *ResponseWithHistoryIVDecPVARPExt*
| _____

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr == Δ *PacemakerWithoutHistoryGbl*

PacemakerWithoutHistoryInit _____
| *PacemakerWithoutHistoryGbl*
| *PacingWithHistoryIInit*
| *SensingInit*
| *ResponseWithHistoryIInit*
| *BatteryInit*
| _____

PacemakerWithoutHistoryVSetBradycardiaMode _____
| *PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr*
| *PacingWithHistoryIVSetBradycardiaMode*
| *SensingVSetBradycardiaMode*
| *ResponseWithHistoryIVSetBradycardiaMode*
| *BatteryVSetBradycardiaMode*
| _____

PacemakerWithoutHistoryVUpdateBatteryStatus _____
| *PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr*
| *PacingWithHistoryIVUpdateBatteryStatus*
| *SensingVUpdateBatteryStatus*
| *ResponseWithHistoryIVUpdateBatteryStatus*
| *BatteryVUpdateBatteryStatus*
| _____

PacemakerWithoutHistoryVExecuteEvent

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVExecuteEvent
SensingVExecuteEvent
ResponseWithHistoryIVExecuteEvent

PacemakerWithoutHistoryVTriggerAPace

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVTriggerAPace
SensingVTriggerAPace
ResponseWithHistoryIVTriggerAPace

PacemakerWithoutHistoryVTriggerVPace

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVTriggerVPace
SensingVTriggerVPace
ResponseWithHistoryIVTriggerVPace

PacemakerWithoutHistoryVIInhibitAPace

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVInhibitAPace
SensingVIInhibitAPace
ResponseWithHistoryIVInhibitAPace

PacemakerWithoutHistoryVIInhibitVPace

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVInhibitVPace
SensingVIInhibitVPace
ResponseWithHistoryIVInhibitVPace

PacemakerWithoutHistoryVTrackASense

PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVTrackASense
SensingVTrackASense
ResponseWithHistoryIVTrackASense

PacemakerWithoutHistoryVTimedAPace _____
PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVTimedAPace

PacemakerWithoutHistoryVTimedVPace _____
PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
PacingWithHistoryIVTimedVPace

PacemakerWithoutHistoryVASenseOnly _____
PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
SensingVASenseOnly
ResponseWithHistoryIVASenseOnly

PacemakerWithoutHistoryVVSenseOnly _____
PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
SensingVVSenseOnly
ResponseWithHistoryIVVSenseOnly

PacemakerWithoutHistoryVASenseRefractory _____
PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
SensingVASenseRefractory
ResponseWithHistoryIVASenseRefractory

PacemakerWithoutHistoryVVSenseRefractory _____
PacemakerWithoutHistoryMergeExtPacingWithHistoryIFr
SensingVVSenseRefractory
ResponseWithHistoryIVVSenseRefractory

C.18 Package Pacemaker

section **Pacemaker** parents *ZOO_Toolkit*, *Model_Preamble*, *History*, *PacemakerWithoutHistory*

PacemakerGblSt _____
HistoryGbl
PacemakerWithoutHistoryGbl

PacemakerGbl _____
PacemakerGblSt

PacemakerIntExtPacemakerWithoutHistoryFr == Δ *PacemakerGbl* \wedge \exists *HistoryGbl*

PacemakerInit _____
PacemakerGbl
PacemakerWithoutHistoryInit

PacemakerVIncRSDown _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncRSDown

PacemakerVDecRSDown _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecRSDown

PacemakerVIncRSUp _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncRSUp

PacemakerVDecRSUp _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecRSUp

PacemakerVDynamicDelayOn _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDynamicDelayOn

PacemakerVDynamicDelayOff _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDynamicDelayOff

PacemakerVIncFixedAVI

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncFixedAVI

PacemakerVDecFixedAVI

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecFixedAVI

PacemakerVIncMinDynAVDelay

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncMinDynAVDelay

PacemakerVDecMinDynAVDelay

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecMinDynAVDelay

PacemakerVIncMaxDynAVDelay

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncMaxDynAVDelay

PacemakerVDecMaxDynAVDelay

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecMaxDynAVDelay

PacemakerVIncSAVDelayOffset

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncSAVDelayOffset

PacemakerVDecSAVDelayOffset

PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecSAVDelayOffset

PacemakerVIncAPulseAmplReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncAPulseAmplReg

PacemakerVDecAPulseAmplReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecAPulseAmplReg

PacemakerVIncVPulseAmplReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncVPulseAmplReg

PacemakerVDecVPulseAmplReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecVPulseAmplReg

PacemakerVIncAPulseAmplUReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncAPulseAmplUReg

PacemakerVDecAPulseAmplUReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecAPulseAmplUReg

PacemakerVIncVPulseAmplUReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVIncVPulseAmplUReg

PacemakerVDecVPulseAmplUReg _____
PacemakerIntExtPacemakerWithoutHistoryFr
PacemakerWithoutHistoryVDecVPulseAmplUReg

PacemakerVIncAPulseWidth _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncAPulseWidth

PacemakerVDecAPulseWidth _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecAPulseWidth

PacemakerVIncVPulseWidth _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncVPulseWidth

PacemakerVDecVPulseWidth _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecVPulseWidth

PacemakerVIncLRL _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncLRL

PacemakerVDecLRL _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecLRL

PacemakerVIncHRL _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncHRL

PacemakerVDecHRL _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecHRL

PacemakerVSetTargetCCI

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVSetTargetCCI

PacemakerVIncASensitivity

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncASensitivity

PacemakerVDecASensitivity

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecASensitivity

PacemakerVIncVSensitivity

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncVSensitivity

PacemakerVDecVSensitivity

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecVSensitivity

PacemakerVIncABlanking

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncABlanking

PacemakerVDecABlanking

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecABlanking

PacemakerVIncVBlanking

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncVBlanking

PacemakerVDecVBlanking _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecVBlanking

PacemakerVIncURL _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncURL

PacemakerVDecURL _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecURL

PacemakerVIncARP _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncARP

PacemakerVDecARP _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecARP

PacemakerVATROn _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVATROn

PacemakerVATROff _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVATROff

PacemakerVIncATRDuration _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncATRDuration

PacemakerVDecATRDuration _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecATRDuration

PacemakerVIncATRFbTm _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncATRFbTm

PacemakerVDecATRFbTm _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecATRFbTm

PacemakerVIncATREntryC _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncATREntryC

PacemakerVDecATREntryC _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecATREntryC

PacemakerVIncATRExitC _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncATRExitC

PacemakerVDecATRExitC _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecATRExitC

PacemakerVIncVRP _____

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncVRP

PacemakerVDecVRP

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecVRP

PacemakerVIncPVARP

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncPVARP

PacemakerVDecPVARP

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecPVARP

PacemakerVIncPVARPExt

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVIncPVARPExt

PacemakerVDecPVARPExt

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVDecPVARPExt

PacemakerVSetBradycardiaMode

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVSetBradycardiaMode

PacemakerVUpdateBatteryStatus

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVUpdateBatteryStatus

PacemakerVExecuteEvent

PacemakerIntExtPacemakerWithoutHistoryFr

PacemakerWithoutHistoryVExecuteEvent

PacemakerJoinExt1Fr == Δ *PacemakerGbl*

PacemakerVASenseRefractory _____

PacemakerJoinExt1Fr

PacemakerWithoutHistoryVASenseRefractory

HistoryCreateHistoryAEv[*am!*/*am?*, *augm!*/*augm?*, *mRe!*/*mRe?*, *mhy!*/*mhy?*, *msr!*/*msr?*, *mUp!*/*mUp?*, *mDo*

PacemakerVASenseOnly _____

PacemakerJoinExt1Fr

PacemakerWithoutHistoryVASenseOnly

HistoryCreateHistoryAEv[*am!*/*am?*, *augm!*/*augm?*, *mRe!*/*mRe?*, *mhy!*/*mhy?*, *msr!*/*msr?*, *mUp!*/*mUp?*, *mDo*

PacemakerVTriggerAPace _____

PacemakerJoinExt1Fr

PacemakerWithoutHistoryVTriggerAPace

HistoryCreateHistoryAEv[*am!*/*am?*, *augm!*/*augm?*, *mRe!*/*mRe?*, *mhy!*/*mhy?*, *msr!*/*msr?*, *mUp!*/*mUp?*, *mDo*

PacemakerVIInhibitAPace _____

PacemakerJoinExt1Fr

PacemakerWithoutHistoryVIInhibitAPace

HistoryCreateHistoryAEv[*am!*/*am?*, *augm!*/*augm?*, *mRe!*/*mRe?*, *mhy!*/*mhy?*, *msr!*/*msr?*, *mUp!*/*mUp?*, *mDo*

PacemakerVTrackASense _____

PacemakerJoinExt1Fr

PacemakerWithoutHistoryVTrackASense

HistoryCreateHistoryAEv[*am!*/*am?*, *augm!*/*augm?*, *mRe!*/*mRe?*, *mhy!*/*mhy?*, *msr!*/*msr?*, *mUp!*/*mUp?*, *mDo*

PacemakerVTimedAPace _____

PacemakerJoinExt1Fr

PacemakerWithoutHistoryVTimedAPace

HistoryCreateHistoryAEv[*am!*/*am?*, *augm!*/*augm?*, *mRe!*/*mRe?*, *mhy!*/*mhy?*, *msr!*/*msr?*, *mUp!*/*mUp?*, *mDo*

PacemakerJoinExt2Fr == Δ *PacemakerGbl*

PacemakerVVSenseRefractory

PacemakerJoinExt2Fr

PacemakerWithoutHistoryVVSenseRefractory

HistoryCreateHistoryVEv[vm!/vm?, augm!/augm?, mRe!/mRe?, mhy!/mhy?, msr!/msr?, mUp!/mUp?, mDo

PacemakerVVSenseOnly

PacemakerJoinExt2Fr

PacemakerWithoutHistoryVVSenseOnly

HistoryCreateHistoryVEv[vm!/vm?, augm!/augm?, mRe!/mRe?, mhy!/mhy?, msr!/msr?, mUp!/mUp?, mDo

PacemakerVTriggerVPace

PacemakerJoinExt2Fr

PacemakerWithoutHistoryVTriggerVPace

HistoryCreateHistoryVEv[vm!/vm?, augm!/augm?, mRe!/mRe?, mhy!/mhy?, msr!/msr?, mUp!/mUp?, mDo

PacemakerVTimedVPace

PacemakerJoinExt2Fr

PacemakerWithoutHistoryVTimedVPace

HistoryCreateHistoryVEv[vm!/vm?, augm!/augm?, mRe!/mRe?, mhy!/mhy?, msr!/msr?, mUp!/mUp?, mDo

PacemakerVIInhibitVPace

PacemakerJoinExt2Fr

PacemakerWithoutHistoryVIInhibitVPace

HistoryCreateHistoryVEv[vm!/vm?, augm!/augm?, mRe!/mRe?, mhy!/mhy?, msr!/msr?, mUp!/mUp?, mDo

Appendix D

ZOO Toolkit

section *ZOO_Toolkit* **parents** *standard_toolkit*

[*OBJ*]

relation(*opt_*)

$$\begin{array}{c}
 \models [X] = \boxed{\begin{array}{l} \text{opt_} : \mathbb{P}(\mathbb{P} X) \\ \text{the} : \mathbb{P} X \rightarrow X \\ \hline \forall S : \mathbb{P} X \bullet \text{opt } S \Leftrightarrow (\exists x : X \bullet S = \{x\}) \vee S = \{\} \\ \forall x : X \bullet \text{the } \{x\} = x \end{array}}
 \end{array}$$

$$\begin{array}{c}
 \models [L] = \boxed{\begin{array}{l} \Sigma : (L \nrightarrow \mathbb{Z}) \rightarrow \mathbb{Z} \\ \hline \Sigma \{\} = 0 \\ \forall l : L; n : \mathbb{Z} \bullet \Sigma \{(l, n)\} = n \\ \forall l : L; n : \mathbb{Z}; S : L \nrightarrow \mathbb{Z} \mid \neg l \in \text{dom } S \bullet \Sigma(\{(l, n)\} \cup S) = n + \Sigma S \end{array}}
 \end{array}$$

MultTy ::= *mm* | *mo* | *om* | *mzo* | *zom* | *mlo* | *lom* | *lolo* | *loo* | *olo* | *lozo* | *zolo*
 | *oo* | *zozo* | *zoo* | *ozo* | *ms* | *sm* | *ss* | *so* | *os* | *szo* | *zos*

relation(*mult_*)

$= [X, Y] =$
$mult_ : \mathbb{P}((X \leftrightarrow Y) \times \mathbb{P} X \times \mathbb{P} Y \times MultTy \times \mathbb{F}\mathbb{N} \times \mathbb{F}\mathbb{N})$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, mm, s_1, s_2)) \Leftrightarrow r \in sx \leftrightarrow sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, mo, s_1, s_2)) \Leftrightarrow r \in sx \rightarrow sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, om, s_1, s_2)) \Leftrightarrow r \sim \in sy \rightarrow sx$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, mzo, s_1, s_2)) \Leftrightarrow r \in sx \nrightarrow sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, zom, s_1, s_2)) \Leftrightarrow r \sim \in sy \nrightarrow sx$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, mlo, s_1, s_2)) \Leftrightarrow r \in sx \leftrightarrow sy \wedge \text{dom } r = sx$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, lom, s_1, s_2)) \Leftrightarrow r \in sx \leftrightarrow sy \wedge \text{ran } r = sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, lolo, s_1, s_2)) \Leftrightarrow r \in sx \leftrightarrow sy \wedge \text{dom } r = sx \wedge \text{ran } r = sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, loo, s_1, s_2)) \Leftrightarrow r \in sx \nrightarrow sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, olo, s_1, s_2)) \Leftrightarrow r \sim \in sy \nrightarrow sx$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, lozo, s_1, s_2)) \Leftrightarrow r \in sx \nrightarrow sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, zolo, s_1, s_2)) \Leftrightarrow r \sim \in sy \nrightarrow sx$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, oo, s_1, s_2)) \Leftrightarrow r \in sx \rightarrowtail sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, zozo, s_1, s_2)) \Leftrightarrow r \in sx \leftrightharpoons sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, zoo, s_1, s_2)) \Leftrightarrow r \in sx \rightarrowtail sy$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F}\mathbb{N} \bullet$
$(mult(r, sx, sy, ozo, s_1, s_2)) \Leftrightarrow r \sim \in sy \rightarrowtail sx$

$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, ms, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, mm, s_1, s_2))$ $\wedge (\forall x : \text{dom } r \bullet \#(\{x\} \triangleleft r) \in s_1)$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, sm, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, mm, s_1, s_2))$ $\wedge (\forall y : \text{ran } r \bullet \#(r \triangleright \{y\}) \in s_1)$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, ss, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, ms, s_1, \{\}))$ $\wedge (mult(r, sx, sy, sm, s_2, \{\}))$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, so, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, mo, s_1, s_2))$ $\wedge (mult(r, sx, sy, sm, s_1, s_2))$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, os, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, om, \{\}, \{\}))$ $\wedge (mult(r, sx, sy, ms, s_1, \{\}))$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, szo, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, mzo, \{\}, \{\}))$ $\wedge (mult(r, sx, sy, sm, s_1, \{\}))$
$\forall r : X \leftrightarrow Y; sx : \mathbb{P} X; sy : \mathbb{P} Y; s_1, s_2 : \mathbb{F} \mathbb{N} \bullet$ $(mult(r, sx, sy, zos, s_1, s_2)) \Leftrightarrow (mult(r, sx, sy, zom, \{\}, \{\}))$ $\wedge (mult(r, sx, sy, ms, s_1, \{\}))$
